

Aluminum Commode with Wheels Assembled



FIGURE A



FIGURE B

- Brass push button allows for easy, tool-free assembly and disassembly of back (Figure A)
- Non-skid, rust-resistant, swivel casters (2 with locks) (Figure B)
- Aluminum construction is lightweight and durable
- Angled back for support and comfort
- Full length plastic arms for added comfort
- Comes with a 7.5 qt commode bucket, splash shield and cover
- splash shield acts as pail holder

PRODUCTS

Description

11112-2	Commode Alum, 2/cs
---------	--------------------

SPECIFICATIONS

Description

Width	19.5" between arms 23.5" Outside
Seat Dimensions	13.75" (W) x 16.5" (D)
Height	19.25"
Weight	12 lbs.
Carton Shipping Weight	28 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

ACCESSORIES

Description

RTL12085	Biodegradable Sanitary Commode Liner, 12/bx, 6 bx/cs
----------	--