



## Serenity<sup>13</sup> COLLECTION

	Medium	Medium
	SR-358M <sup>4,16</sup>	SR-525M <sup>11</sup>
	3-Position Chaise Lounger 5'4"- 6'0"	Infinite-Position Zero-gravity Trendelenburg Chaise Lounger 5'4"- 6'0"
Weight Capacity	375 lbs.	375 lbs.
Seat-to-Floor	20"	20"
Seat Depth	20"	20"
Seat Width	20"	20"
Top of Back to Seat	24.5"	26.5"
Back Style	Blown Fiber Button Back	Blown Fiber Split Back
Fully Padded Chaise	Yes	Yes
Heat & Massage Option	No	No
Total Comfort Seating <sup>5</sup>	Yes	No
Distance From Wall	17"	29"
Overall Chair Width	35"	33"
Chair Weight	117 lbs.	126 lbs.
Hand Control	Easy Touch with USB Port	Easy Touch Infinite-Position with USB Port
Warranty	7-Year Prorated	7-Year Prorated



## ELEGANCE COLLECTION

	Medium	Petite	Small	Medium	Large
	LC-450C <sup>4,8,10,18</sup>	LC-570P <sup>4,8,18</sup>	LC-570S <sup>4,8,18</sup>	LC-570M <sup>4,8,18</sup>	LC-570L <sup>4,11,18</sup>
	3-Position Chaise Lounger 5'4"- 6'0"	3-Position Chaise Lounger 5'2"- & Below	3-Position Chaise Lounger 5'3" & Below	3-Position Chaise Lounger 5'4"- 6'0"	3-Position Chaise Lounger 5'9"- 6'2"
Weight Capacity	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.
Seat-to-Floor	20"	18"	19"	20"	21"
Seat Depth	20"	18"	19"	20"	22"
Seat Width	20"	18"	19"	20"	22"
Top of Back to Seat	27"	22.5"	24.5"	26.5"	27"
Back Style	Blown Fiber Split Back	Blown Fiber Pillow Back	Blown Fiber Pillow Back	Blown Fiber Waterfall Back	Blown Fiber Waterfall Back
Fully Padded Chaise	Yes	Yes	Yes	Yes	Yes
Heat & Massage Option	Yes	Yes	Yes	Yes	Yes
Total Comfort Seating <sup>5</sup>	Yes	Yes	Yes	Yes	Yes
Distance From Wall	17"	15.5"	16.5"	19"	24"
Overall Chair Width	38"	35"	35"	36"	38"
Chair Weight	119 lbs.	106 lbs.	110 lbs.	113 lbs.	114 lbs.
Hand Control	Easy Touch with USB Port	Easy Touch with USB Port	Easy Touch with USB Port	Easy Touch with USB Port	Easy Touch with USB Port
Warranty	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated



## Classic COLLECTION

	Medium	Medium
	LC-220 <sup>19</sup>	LC-250 <sup>19</sup>
	2-Position 5'4"- 6'0"	3-Position Chaise Lounger 5'4"- 6'0"
Weight Capacity	375 lbs.	375 lbs.
Seat-to-Floor	20"	20"
Seat Depth	20"	20"
Seat Width	20"	20"
Top of Back to Seat	24.5"	26"
Back Style	Blown Fiber Sewn Pillow Back	Blown Fiber Split Back
Fully Padded Chaise	No	Yes
Heat & Massage Option	No	No
Total Comfort Seating <sup>5</sup>	No	Yes
Distance From Wall	11"	18"
Overall Chair Width	32.5"	34"
Chair Weight	114 lbs.	118 lbs.
Hand Control	Easy Touch	Easy Touch
Warranty	2-Year Limited	2-Year Limited



## INFINITY COLLECTION

	Medium	Petite Wide	Small	Medium	Large	Medium	Large
	LC-108 <sup>10, 11</sup>	LC-525PW <sup>9,11</sup>	LC-525S <sup>9</sup>	LC-525M <sup>9,11</sup>	LC-525L <sup>9,11</sup>	LC-580M <sup>4,9,17</sup>	LC-580L <sup>4,9,17</sup>
	Infinite-Position Zero-gravity Trendelenburg Chaise Lounger 5'4"- 6'0"	Infinite-Position Zero-gravity Trendelenburg Chaise Lounger 5'3" & Below	Infinite-Position Zero-gravity Trendelenburg Chaise Lounger 5'3" & Below	Infinite-Position Zero-gravity Trendelenburg Chaise Lounger 5'4"- 6'0"	Infinite-Position Zero-gravity Trendelenburg Chaise Lounger 5'9"- 6'2"	Infinite-Position Zero-gravity Trendelenburg Chaise Lounger 5'4"- 6'0"	Infinite-Position Zero-gravity Trendelenburg Chaise Lounger 5'9"- 6'2"
Weight Capacity	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.
Seat-to-Floor	19"	19"	19"	20"	21"	20"	21"
Seat Depth	19"	19"	18"	20"	22"	20"	22"
Seat Width	20"	24"	18"	20"	22"	20"	22"
Top of Back to Seat	29"	25.5"	25.5"	26.5"	27.5"	27.5"	28"
Back Style	Sewn Pillow Back	Blown Fiber Split Back	Blown Fiber Split Back	Blown Fiber Split Back	Blown Fiber Split Back	Blown Fiber Biscuit Back	Blown Fiber Biscuit Back
Fully Padded Chaise	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Heat & Massage Option	No	Yes	Yes	Yes	Yes	Yes	Yes
Total Comfort Seating <sup>5</sup>	No	No	No	No	No	No	No
Distance From Wall	24"	28"	28"	29"	30"	29"	30"
Overall Chair Width	31.5"	38"	32"	33"	35"	37"	39"
Chair Weight	136 lbs.	130 lbs.	123 lbs.	126 lbs.	128 lbs.	130 lbs.	132 lbs.
Hand Control	Easy Touch Infinite-Position with USB Port	Easy Touch Infinite-Position with USB Port	Easy Touch Infinite-Position with USB Port	Easy Touch Infinite-Position with USB Port	Easy Touch Infinite-Position with USB Port	Easy Touch Infinite-Position with USB Port	Easy Touch Infinite-Position with USB Port
Warranty	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated




## HERITAGE COLLECTION

	Petite	Petite Wide	Small	Medium	Large	Heavy-duty	Very Heavy-duty
	LC-358P <sup>5,18</sup>	LC-358PW <sup>4,8,18</sup>	LC-358S <sup>9,18</sup>	LC-358M <sup>4,8,18</sup>	LC-358L <sup>4,8,18</sup>	LC-358XL <sup>1,6</sup>	LC-358XXL <sup>1,6</sup>
	3-Position Chaise Lounger 5'2"- & Below	3-Position Chaise Lounger 5'3" & Below	3-Position Chaise Lounger 5'3" & Below	3-Position Chaise Lounger 5'4"- 6'0"	3-Position Chaise Lounger 5'9"- 6'2"	3-Position Chaise Lounger 5'6"- 6'1"	2-Position Chaise Lounger 5'6"- 6'1"
Weight Capacity	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	500 lbs.	600 lbs.
Seat-to-Floor	18"	19"	19"	20"	21"	20"	20"
Seat Depth	18"	19"	19"	20"	22"	22"	22"
Seat Width	18"	24"	18"	20"	22"	26"	30"
Top of Back to Seat	24"	23.5"	24"	24.5"	25.5"	29"	29"
Back Style	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back
Fully Padded Chaise	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Heat & Massage Option	Yes	Yes	Yes	Yes	Yes	Yes <sup>15</sup>	Yes <sup>15</sup>
Total Comfort Seating <sup>5</sup>	Yes	Yes	Yes	Yes	Yes	No	No
Distance From Wall	17.5"	17"	17"	17"	17.5"	18"	10"
Overall Chair Width	34.5"	39"	35"	35"	37"	39.5"	43"
Chair Weight	112 lbs.	118 lbs.	118 lbs.	117 lbs.	120 lbs.	159 lbs.	168 lbs.
Hand Control	Easy Touch with USB Port	Easy Touch with USB Port	Easy Touch with USB Port	Easy Touch with USB Port	Easy Touch with USB Port	Easy Touch with USB Port	Easy Touch with USB Port
Warranty	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated



## essential COLLECTION

	Super Petite	Medium	Medium
	LC-102	LC-106	LC-107 <sup>11</sup>
	3-Position Chaise Lounger 4'8"- 5'4"	3-Position Chaise Lounger 5'2"- 5'10"	Infinite-Position Chaise Lounger 5'4"- 6'0"
Weight Capacity	250 lbs.	325 lbs.	325 lbs.
Seat-to-Floor	16"	18"	20"
Seat Depth	18"	19"	20"
Seat Width	16"	20"	20"
Top of Back to Seat	23.5"	27.5"	27.5"
Back Style	Dual Pillow Back	Sewn Pillow Back	Sewn Pillow Back
Fully Padded Chaise	Yes	Yes	Yes
Heat & Massage Option	No	No	No
Total Comfort Seating <sup>5</sup>	No	No	No
Distance From Wall	16"	18"	26"
Overall Chair Width	25.5"	32.75"	35"
Chair Weight	96 lbs.	96 lbs.	116 lbs.
Hand Control	Easy Touch	Easy Touch	Easy Touch Infinite-Position with USB Port
Warranty	1-Year Limited	1-Year Limited	1-Year Limited



## SPECIALTY COLLECTION

	Medium
	LC-805
	2-Position Chaise Lounger Wall Hugger 5'4"- 6'0"
Weight Capacity	300 lbs.
Seat-to-Floor	19"
Seat Depth	21"
Seat Width	20"
Top of Back to Seat	27"
Back Style	Blown Fiber Sewn Split Back
Fully Padded Chaise	Yes
Heat & Massage Option	Yes
Total Comfort Seating <sup>5</sup>	Yes
Distance From Wall	4"
Overall Chair Width	35"
Chair Weight	118 lbs.
Hand Control	Easy Touch with USB Port
Warranty	2-Year Limited



## HERITAGE COLLECTION

	Petite	Petite Wide	Small	Medium	Large	Heavy-duty	Very Heavy-duty
	LC-358P <sup>5,18</sup>	LC-358PW <sup>4,8,18</sup>	LC-358S <sup>9,18</sup>	LC-358M <sup>4,8,18</sup>	LC-358L <sup>4,8,18</sup>	LC-358XL <sup>1,6</sup>	LC-358XXL <sup>1,6</sup>
	3-Position Chaise Lounger 5'2"- & Below	3-Position Chaise Lounger 5'3" & Below	3-Position Chaise Lounger 5'3" & Below	3-Position Chaise Lounger 5'4"- 6'0"	3-Position Chaise Lounger 5'9"- 6'2"	3-Position Chaise Lounger 5'6"- 6'1"	2-Position Chaise Lounger 5'6"- 6'1"
Weight Capacity	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	500 lbs.	600 lbs.
Seat-to-Floor	18"	19"	19"	20"	21"	20"	20"
Seat Depth	18"	19"	19"	20"	22"	22"	22"
Seat Width	18"	24"	18"	20"	22"	26"	30"
Top of Back to Seat	24"	23.5"	24"	24.5"	25.5"	29"	29"
Back Style	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back
Fully Padded Chaise	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Heat & Massage Option	Yes	Yes	Yes	Yes	Yes	Yes <sup>15</sup>	Yes <sup>15</sup>
Total Comfort Seating <sup>5</sup>	Yes	Yes	Yes	Yes	Yes	No	No
Distance From Wall	17.5"	17"	17"	17"	17.5"	18"	10"
Overall Chair Width	34.5"	39"	35"	35"	37"	39.5"	43"
Chair Weight	112 lbs.	118 lbs.	118 lbs.	117 lbs.	120 lbs.	159 lbs.	168 lbs.
Hand Control	Easy Touch with USB Port	Easy Touch with USB Port	Easy Touch with USB Port	Easy Touch with USB Port	Easy Touch with USB Port	Easy Touch with USB Port	Easy Touch with USB Port
Warranty	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated



## HERITAGE COLLECTION

	Petite	Petite Wide	Small	Medium	Large	Heavy-duty	Very Heavy-duty
	LC-358P <sup>5,18</sup>	LC-358PW <sup>4,8,18</sup>	LC-358S <sup>9,18</sup>	LC-358M <sup>4,8,18</sup>	LC-358L <sup>4,8,18</sup>	LC-358XL <sup>1,6</sup>	LC-358XXL <sup>1,6</sup>
	3-Position Chaise Lounger 5'2"- & Below	3-Position Chaise Lounger 5'3" & Below	3-Position Chaise Lounger 5'3" & Below	3-Position Chaise Lounger 5'4"- 6'0"	3-Position Chaise Lounger 5'9"- 6'2"	3-Position Chaise Lounger 5'6"- 6'1"	2-Position Chaise Lounger 5'6"- 6'1"
Weight Capacity	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	500 lbs.	600 lbs.
Seat-to-Floor	18"	19"	19"	20"	21"	20"	20"
Seat Depth	18"	19"	19"	20"	22"	22"	22"
Seat Width	18"	24"	18"	20"	22"	26"	30"
Top of Back to Seat	24"	23.5"	24"	24.5"	25.5"	29"	29"
Back Style	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back
Fully Padded Chaise	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Heat & Massage Option	Yes	Yes	Yes	Yes	Yes	Yes <sup>15</sup>	Yes <sup>15</sup>
Total Comfort Seating <sup>5</sup>	Yes	Yes	Yes	Yes	Yes	No	No
Distance From Wall	17.5"	17"	17"	17"	17.5"	18"	10"
Overall Chair Width	34.5"	39"	35"	35"	37"	39.5"	43"
Chair Weight	112 lbs.	118 lbs.	118 lbs.	117 lbs.	120 lbs.	159 lbs.	168 lbs.
Hand Control	Easy Touch with USB Port	Easy Touch with USB Port	Easy Touch with USB Port	Easy Touch with USB Port	Easy Touch with USB Port	Easy Touch with USB Port	Easy Touch with USB Port
Warranty	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated











## HERITAGE COLLECTION

	Petite	Petite Wide	Small	Medium	Large	Heavy-duty	Very Heavy-duty
	LC-358P <sup>5,18</sup>	LC-358PW <sup>4,8,18</sup>	LC-358S <sup>9,18</sup>	LC-358M <sup>4,8,18</sup>	LC-358L <sup>4,8,18</sup>	LC-358XL <sup>1,6</sup>	LC-358XXL <sup>1,6</sup>
	3-Position Chaise Lounger 5'2"- & Below	3-Position Chaise Lounger 5'3" & Below	3-Position Chaise Lounger 5'3" & Below	3-Position Chaise Lounger 5'4"- 6'0"	3-Position Chaise Lounger 5'9"- 6'2"	3-Position Chaise Lounger 5'6"- 6'1"	2-Position Chaise Lounger 5'6"- 6'1"
Weight Capacity	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	500 lbs.	600 lbs.
Seat-to-Floor	18"	19"	19"	20"	21"	20"	20"
Seat Depth	18"	19"	19"	20"	22"	22"	22"
Seat Width	18"	24"	18"	20"	22"	26"	30"
Top of Back to Seat	24"	23.5"	24"	24.5"	25.5"	29"	29"
Back Style	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back
Fully Padded Chaise	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Heat & Massage Option	Yes	Yes	Yes	Yes	Yes	Yes <sup>15</sup>	Yes <sup>15</sup>
Total Comfort Seating <sup>5</sup>	Yes	Yes	Yes	Yes	Yes	No	No
Distance From Wall	17.5"	17"	17"	17"	17.5"	18"	10"
Overall Chair Width	34.5"	39"	35"	35"	37"	39.5"	43"
Chair Weight	112 lbs.	118 lbs.	118 lbs.	117 lbs.	120 lbs.	159 lbs.	168 lbs.
Hand Control	Easy Touch with USB Port	Easy Touch with USB Port	Easy Touch with USB Port	Easy Touch with USB Port	Easy Touch with USB Port	Easy Touch with USB Port	Easy Touch with USB Port
Warranty	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated



## HERITAGE COLLECTION

								
	Lifted	Seated	Reading	Napping	Relaxing	Sleeping	Trendelenburg	Zero-Gravity
2-Position	The 2-position TV recline is achievable on all Pride® Lift Chairs. This is a comfortable position for performing many relaxing activities including reading, watching TV, working on your tablet or for conversation with family and friends.							
3-Position	The 3-position full recline is available on most Pride® Lift Chair models and puts you in a position ideal for napping or relaxing.							
Infinite Position	The infinite-position sleep recline is perfect for taking a nap - whether quick or extended - or even for a full night's sleep. Infinite-position Pride® Lift Chairs enable the backrest and footrest to move independently from one another so that a variety of positions can be reached.							
Trendelenburg	The Trendelenburg position, elevates one's feet higher than his or her head for increased comfort and relaxation.							
Zero gravity	The zero gravity position, which can be reached on infinite-position Pride® Lift Chair models, positions one's thighs and lower legs higher than his or her torso.							