



2-in-1 Locking, Raised Toilet Seat with Tool-free Removable Arms



FIGURE A



FIGURE B



RETAIL BOX

- Designed for individuals who have difficulty sitting down or getting up from the toilet
- Adds 5" to height of toilet seat
- Tool-free arms can be removed or added (Figure A - B)
- Arm sleeves and receivers are made of metal
- Heavy-duty molded plastic construction
- Standard locking mechanism ensures safety
- Lightweight and portable
- Easily attaches to toilet bowl
- No tools required for installation
- Easy-to-clean
- Fits most toilets

PRODUCTS

	Description
12008KDR	Knocked Down, 1/cs

SPECIFICATIONS

	Description
Seat Dimensions	15" (W) x 16.5" (D) x 5" (H)
Width Between Arms	18.5"
Weight	5 lbs.
Carton Shipping Weight	8 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime