



## Aluminum Shower Chair and Commode with Casters Knocked Down



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Easy-to-assemble frame
- Padded, open front, vinyl toilet seat (Figure A)
- Aluminum construction is lightweight and durable
- 3" non-skid, rust resistant, swivel casters (2 rear locking) (Figure B)
- Brass push button allows for easy assembly (Figure C)
- Removable, tool-free back
- Plastic arms for added comfort
- Comes with removable 12 qt commode bucket and cover (Figure D)

### PRODUCTS

#### Description

11114KD-1	Knock Down, 1/cs
-----------	------------------

### SPECIFICATIONS

#### Description

Width	18" between arms   21.5" Outside
Seat Dimensions	16" (W) x 16" (D)
Height	21"
Weight	14 lbs.
Carton Shipping Weight	18 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

### ACCESSORIES

#### Description

RTL12085	Biodegradable Sanitary Commode Liner, 12/bx, 6 bx/cs
----------	--