

FUZE T50



pelg
MOBILITY

FUZE T50



Features

50° of manual centre-of-gravity tilt	Narrow overall width
Adjustable axle depth and seat height	Short wheelbase for maneuverability
Adjustable seat width and seat depth	Pediatric version available (Fuze JR)
Constant center of gravity tilt	Power Tilt available
Durable Aluminum frame	Transit tie-down

With its effortless center-of-gravity pivot, the Fuze T50 manual tilt-in-space wheelchair offers up to 50° of tilt. It is designed for users who need pressure relief or routine repositioning for postural or access reasons. Offering a wide range of adjustability and options, the T50 is one of the most adaptable wheelchairs available.

A Adjustable Seat Height

13" to 20" of front seat height to optimize occupant tilt performance.



B Adjustable Axle and Seat Position

Adjust the center of gravity location to optimize mobility requirements.



C Adjustable Seat Depth and Back Rest

A large selection of seat depth options and up to 36 degrees of back angle adjustment for general positioning and comfort.



D Adjustable Chair Width

On-board adjustability of the chair needs to accommodate occupants growing need.



ADJ SEAT WIDTH
12"–22"
FIXED SEAT WIDTH
15"–24"



SEAT DEPTH
15"–22"
ALSO AVAILABLE
23" & 24"



FRONT SEAT HEIGHT
13"–20"



BACK ANGLE
84°–120°



DYNAMIC TILT
Up to 50°



MAX CAPACITY ADJUSTABLE FRAME
250 lbs/113 kg
MAX CAPACITY FIXED FRAME
350 lbs/158 kg



CHAIR WEIGHT
46 lbs/21 kg



OVERALL WIDTH
7.5" (RW 12/16)
8.25" (20+)