



BRODA[®]

Comfort is Essential

ELITE
TILT RECLINER

INDUSTRY LEADING POSITIONING VERSATILITY AND FUNCTION

Enhanced positioning and innovative features provide exceptional comfort. Superior durability accommodates use in multiple care settings.



- 1 INFINITELY ADJUSTABLE SEAT TILT
- 2 INFINITELY ADJUSTABLE BACK RECLINE (TRENDLENBURG CAPABLE)
- 3 HEIGHT ADJUSTABLE ARMRESTS
- 4 SWING AWAY REMOVABLE ARMS
- 5 REMOVABLE WINGS FOR UPPER LATERAL SUPPORT
- 6 SPLIT FLIP DOWN FOOTREST
- 7 MULTI-DIRECTIONAL LOCKING REAR CASTERS
- 8 INDEPENDENT ELEVATING, ARTICULATING, LENGTH ADJUSTABLE LEGREST

PATIENT BENEFITS



Posterior tilt effectively opens the diaphragm improving digestion, oxygenation, and organ function.



Tilt in Space feature facilitates positioning for comfort, pain management and pain relief.



Lay flat feature facilitates supine transfers, reducing the risk associated with falls and injuries.



Provides ease of maneuverability to accommodate various activities (Middle Mag Wheel Option Shown).

CAREGIVER AND FACILITY BENEFITS

- Facilitates positioning for comfort, pain management and improved respiratory function
- Optimum safety and functionality withstands severe agitation and involuntary movements
- Tilt and recline provide ease of use frequent re-positioning, helping improve patient outcomes
- Proper seating alignment decreases patient fatigue and postural deviations reducing the use of restraints
- Infinitely adjustable positioning cradles and soothes patient to reduce agitation and provide contentment

COMFORT TENSION SEATING® PROVIDES

- Long term seating comfort
- Conforming support for proper pelvic positioning
- Improved air flow, decreasing heat and moisture build up
- Pressure redistribution that aids in maintaining skin integrity
- Accommodation for conditions such as kyphosis and lordosis
- Decreased pressure, improved blood flow and oxygenation of the tissue

SPECIFICATIONS

Frame	16 Gauge tubular steel frame
Seat & Backrest	BRODA Comfort Tension Seating®
Degrees of Tilt	40°
Degrees of Recline	90°(Trendelenburg capable)
Mag Wheel Option	20" or 22" Front or rear mount (50cm, 55cm) 14" rear or mid mount (35cm)
Standard Seat Width	16", 18", 20" (40cm, 45cm, 50cm)
Custom Seat Width	22", 24", 26", 28" (55cm, 60cm, 65cm, 70cm)
Standard Seat Depth	17", 18.5" (43cm, 46cm)
Custom Seat Depth	15.5", 20" (39cm, 50cm)
Standard Seat Height	20.5" (51cm)
Custom Seat Height	19", 22", 23.5" (48cm, 55cm, 59cm)
Front Casters	5" Front total lock (12.5cm)
Rear Casters	6" Multi-directional locking (15cm)
Weight Capacity	350lbs. / 159kg
Standard Padding	Headrest cushion, shoulder bolsters, lower side pads, upper side pads, seat and back pad with cover, calf pad, sole pads
ADP Device Code*	WMA5B-020, WMA5B-010

* Specific configurations apply. Consult your rep for details

COMMON ACCESSORIES

- Mag wheel package
- Aluminum front casters
- Twin wheel front casters
- Huntington's padding package (HSP)
- Padded armrests
- Full wrap vinyl seat and back covers
- BRODA app (additional positioning package)
- IV Pole
- O2 Holder
- ABS Legrest
- Solid or clear tray
- Half tray
- Padded tray cover
- Thigh belt



CTS Strapping flexes for individualized fit & superior comfort.

UPHOLSTERY OPTIONS



RENAISSANCE
BURGUNDY



RENAISSANCE
CHARCOAL



RENAISSANCE
METEOR



SILVERTEX
SAPPHIRE BLUE

FRAME OPTIONS



GREY
HAMMER



BLACK



HEALTHCARE DIVISION

info@brodaseating.com

p. 1-844-552-7632

Brodaseating.com

©Broda USA, INC.

8/16