TRILLIUM **SERIES**

The Trillium Series is the most comfortable walker on the market today.

With the new padded seat cushion and wide back strap pad, users can sit in comfort and walk in style. The Trillium comes with an all-new braking system, allowing for effortless handle height adjustment and easy maintenance.

FFATURES.

Seat 1.5" thick padded seat cushion for extra comfort Brake System Concealed cable / Easy to adjust and maintain Collapsible Basket Foldable & removable Wheels Lightweight 8" EVA tires / Offset rear forks for wider gait Wide Back Strap Pad Extra wide for added comfort *Colors Chinatown Blue, Carbon Black, Pink



SPECIFICATIONS.

Regular

14 lbs 18" Seat Width 21" Seat Height 26" Width x 25.5" Length 31.5" - 34.5" Handle Height 300 lbs Weight Capacity

Tall

14 lbs 18" Seat Width 24" Seat Height 26" Width x 25.5" Length 34.5" - 37.5" Handle Height 350 lbs Weight Capacity

Order Codes

WIR63B Trillium Regular-Carbon Black WIR64B Trillium Regular-Chinatown Blue WIR65B Trillium Regular-Pink

OPTIONAL ACCESSORIES.

- B35XP Cup Holder
- **B13IR Curb Climber**
- B14IB Slow Down Brake
 - **B21IR Cane Holder**



WIT63B Trillium Tall-Carbon Black WIT64B Trillium Tall-Chinatown Blue WIT65B Trillium Tall-Pink

- B38XRPH4 Phone Holder (4.7")
- B38XRPH5 Phone Holder (5.5")
- **B36XRV Flower Vase**



EVOLUTION TECHNOLOGIES INC. 2530 Davies Avenue. Port Coguitlam, BC V3C2J9

tf +1 800.556.2558 info@evolution.ca

EVOLUTION TECHNOLOGIES USA INC. 2657 Delta Ring Road Ferndale WA 98248

tf +1 800 839 2558 info@evolutionwalker.com

www.evolutionwalker.com Rev3.6W_20170109_SL