


Gravity 7 Long Term Care Pressure Redistribution Mattress



HCPCS: E0184



15876

- High quality therapeutic foam mattress provides pressure redistribution and shear/friction reduction
- Drive's most popular therapeutic foam mattress maximizes value and comfort with 3 full layers of latex free foam
- Full width top layer uses unique die-cutting to optimize zoned pressure redistribution
- Bottom layer provides full length horizontally scored articulation cuts extending the durability and life of the mattress
- High density foam perimeters provide solid edge support and facilitate safe patient transfers
- Available with or without raised side rails
- Visco foam heel section offloads pressure from vulnerable heels
- Removable multi-stretch cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- Meets flammability standard CFR 1632 and CFR 1633

PRODUCTS

PRODUCTS	Description
15876	Multi-Layered/Multi-zoned Foam Mattress, E0184, 1/ea
15770	Multi-Layered/Multi-zoned Foam Mattress, E0184, 1/ea
15784	Multi-Layered/Multi-zoned Foam Mattress, E0184, 1/ea
15877	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out, E0184, 1/ea
15777	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out, E0184, 1/ea
15785	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out, E0184, 1/ea

SPECIFICATIONS

SPECIFICATIONS	Description
Cover Material	Stretch, Fluid-Resistant Nylon Top Waterproof Vinyl Bottom
Mattress Layers of Foam	3
Side Rails	Yes
Dedicated Heel Section	Yes
Weight Capacity	350 lbs.
Warranty	7 Year Limited Foam 1 Year Cover

ACCESSORIES

ACCESSORIES	Description
15770C	Replacement Cover, 1/ea
15777C	Replacement Cover, 1/ea
15784C	Replacement Cover, 1/ea
15785C	Replacement Cover, 1/ea
15876C	Replacement Cover, 1/ea
15877C	Replacement Cover, 1/ea