



#### GO BEYOND!

Motion Composites believes in empowering wheelchair users. Why? Because clunky and inefficient equipment is unfortunately too common. Which means you're less mobile and less efficient than you could be — and should be.

So we did something about it. We use the world's most advanced technologies to enhance people's mobility and improve their lives.

Case in point — the ultra lightweight Helio C2. It offers unrivaled strength and durability through the use of cutting edge carbon fiber technology. With its superior fatigue resistance, extreme temperature tolerances and exceptional durability, carbon fiber is one of the strongest yet lightest materials on earth. That's why it's the ultimate wheelchair manufacturing material.

With its innovative materials and advanced engineering, the Helio C2 is one of the lightest, most technologically advanced wheelchair in its category.





## **SHELIOC**2

### What a folding wheelchair should be.

Look closely. From frame to materials to every detail of the Helio C2 is optimized to help you go farther. We've ratcheted up performance and brought down weight so you actually expend less energy to propel the Helio C2.

Result? You gain in mobility while reducing the risk of chronic shoulder and joint injuries.

The Helio C2's innovative engineering begins with its symmetrical crossbrace and one-piece side frame. This reduces weight and distributes forces equally throughout the frame for superior energy transfer and propulsion efficiency. And, with a vertical axle plate to shave even





## Significantly lighter than the competition.

The Helio C2's outstanding performance and unrivaled lightness has positively changed how people think about wheelchairs. With a transport weight of only 9.5 lbs. (4.3 kg), the Helio C2 is significantly lighter than any other wheelchair in its category.

### The lowest seat-to-floor height.

The innovative frame design lets you achieve a super low  $13 \frac{1}{2}$ " (34.3 cm) seat-to-floor height with a 4" caster on a standard frame.

#### Better design for better mobility.

Even an ultra-lightweight carbon fiber frame can get weighed down without the right components. That's why every single accessory, bolt and bracket of the Helio C2 is designed with weight reduction in mind. From the Flip-Back Convertible Armrests to the forged aluminum axle plate, lightness is a top priority.

# Get more mileage out of your energy.

Everybody benefits from an easier to propel wheelchair. The Symmetrical Molded Crossbrace, oval shaped tubes and oversized pivot axles work in concert, while our unique Ultrarigid Folding System maximizes frame rigidity and energy conservation.







High-Modulus Carbon T700 — One of the lightest and most rigid material available, also renowned for its vibration dampening properties.



**Rigid Unibody Frame** — A unibody frame is much stronger and needs less maintenance than a standard two-part frame. It also reduces weight while maximizing propulsion efficiency.



Symmetrical Molded Crossbrace 3D — Entirely symmetrical carbon fiber crossbrace for reduced torsion and better energy distribution throughout the frame.



**Ultrarigid Folding System** — High-precision tolerances and oversized pivot axles for best-in-class propulsion efficiency.



Forged Vertical Axleplate — Offers the industry's most precise rear wheel adjustability. The vertical mounting maximizes rigidity and responsiveness.



**Evolve Caster Housing** — Integrated into the frame for rock-solid durability, the Evolve Caster Housing offers easy and precise infinite adjustments.



Flip-back Convertible Armrest — Easy to adjust and made of composite materials for strength and reduced weight, it is the only armrest that converts from flipback to single post without the use of tools.



Anti-flutter system - Reduces flutter through the simple turn of a screw. .



#### Technical Specifications

#### Structure

Frame	Folding Unibody Frame
Material	T700 High-Modulus Carbon Fiber
Transport weight	9.5 lbs. (4.3 kg) (16*16, w/o rear wheels & footrest)
Weight capacity	265 lbs. (120 kg) HD Kit: 350 lbs. (159 kg)
Dimensions	
Width	14" (35.6 cm) to 20" (50.8 cm) 18" (45.7 cm) to 22" (55.9 cm) (Heavy Duty Kit)
Depth	14" (35.6 cm) to 20" (50.8 cm)
Front seat-to-floor height	(w/4" caster) (w/8" caster) 13 ½" (34.3 cm) to 21½" (54.6 cm)
Rear seat-to-floor height	(w/20" wheel) (w/26" wheel) 12 ½" (31.8 cm) to 21¼" (54 cm)

#### **Motion Composites**

Quebec, Canada

T (866) 650-6555 F (888) 966-6555 info@motioncomposites.com www.motioncomposites.com





Newton Accessories — Parts and accessories designed to be lighter, with improved functionality and style.



Seat depth 14" to 20" 35.6 cm to 50.8 cm



Seat width 14" to 22" 35.6 to 55.9 cm



Front seat to floor 13" to 211/4" 33.0 cm to 54.0 cm



Rear seat to floor 12 1/2" to 20" 31.8 cm to 50.8 cm



Back height 9" to 21" 22.9 cm to 53.3 cm



Footrest angle 60°, 70°, 80°,90°



Back angle 70° to 110°

Weight limit

265 lbs. (120 kg)

(HD Kit: 350 lbs.



Wheel camber 0°, 3°, 6°





Overall width from +7 1/2" from +19.1 cm



Armrest height 8" to 14" 20.3 cm to 35.6 cm





C2 brochure pages