



FABRIC









ONYX BLACK



VOLCANO RED





PLATINUM SILVER



PURPLE RAIN



PHANTOM GREY

21", 24" - 29"



REALTREE AP CAMO



POP STAR PINK

460, 485, 510 mm

460, 485, 510 mm

335, 410, 460 mm

470, 520, 570, 620 mm 470 - 620 mm

530, 600 - 730 mm

-	DE		\sim	АΤ.	\sim	NS
	ΡЕ			<u>^</u>		1

Maximum User Weight	450 lbs	205 kg	Seat to Floor Height	(Elevator)	18", 19", 20"	460, 48
Maximum Speed	5 MPH	8 km/h		(Seat Tube)	18", 19", 20"	460, 48
Driving Range ¹	17.8 miles	28.6 km	Tilt Options	(Power)	0°	- 45°
Base Width	26"	660 mm		(Manual)	No ma	nual option
Base Length	38"	965 mm	Recline Options	(Power)	85°	' - 150°
Minimum Turning Radius	22"	555 mm		(Manual)	85°	' - 120°
Ground Clearance	3"	75 mm	Legrest Elevation Options	(Power)	85°	' - 175°
Optional Seat Elevator	8″	200mm		(Manual)	85°	' - 175°
Battery Type	Group 24, 0	Group 34	Armrest Pad Lengths		13", 16", 18"	335, 410
Drive Electronics 120A R-n		R-net	Seat Widths		19", 21", 23", 25"	470, 520, 5
Maximum Slope ²	7.5° Up / 7.	.5° Down	Seat Depths		19" - 25"	470 -

¹ Actual driving range will vary based on driving and battery conditions.
² Refer to Owner Manual for safe driving practices

www.permobil.com

M300HD090116

Backrest Heights



INTRODUCING THE NEW M300 CORPUS® HD, a heavy duty power wheelchair designed specifically for users up to 450 lbs. The HD has extremely sturdy arm and legrests as well as wider footplates for additional room, and it features the world-renowned Corpus seating system for the ultimate in comfort with an unmatched power tilt range of 0-45° and an unprecedented recline of 85-150°. The M300 base is highly maneuverable and has been both extended and widened for the HD model, while the suspension has been strengthened for an overall more stable platform. The ultra-strong actuators allow users to safely elevate their feet above their hearts, while the durable mesh-style upholstery helps keep moisture away from the body. With the legendary power and toughness of Permobil, the M300 Corpus HD is **BUILT TO WITHSTAND THE RIGORS OF ANY HEAVY DUTY LIFESTYLE.**

www.permobil.com







EVERYTHING YOU'VE COME TO EXPECT FROM THE CORPUS® SEATING SYSTEM. ONLY STRONGER.

0-45° Power Tilt 85-150° Power Recline Power Elevating Legrests




Adjustable legrests Extremely durable, sturdy legrests



Adjustable footplates Oversized, independently adjustable footplates

Strong actuators Incredibly strong actuators for seat functions

Power tilt, recline and elevating legrests allow users' feet to be positioned above the heart.



Suspension Independently adjustable suspension for a smooth ride under any condition