




STANDER

For a different perspective !



prodije

Since 1997

It Grows

Because children grow so fast, the ProdiJe Stander is adjustable both in width and height to provide comfort as they grow.



It Reclines

The locking gas cylinder mechanism allows for easy angle adjustment from 0 to 90 degrees without any danger for the child.




It's Adaptable



SUPINE

The unique aluminium frame concept allows for either a prone or a supine positioning depending on the needs of the child.

PRONE



The Difference

Easy to disassemble

In a split second, the stander separates into 2 sections thereby reducing the weight and facilitating transportation without losing or modifying the initial adjustments.

The Prodiye stander is a high quality product, resistant and appealing. A refurbishment is possible at an affordable price. Please contact an authorized distributor for more information.

Junior Stander

Minimum Floor/head Height : 115 cm (45")

Maximum Floor/head Height : 170 cm (67")

Maximum Weight : 69 kg (150 lbs)



prodiye

1242, De Lanaudière, Joliette
(Québec) Canada J6E 3P1
Phone : 450 752-8818
Fax : 450 752-1644

www.prodiye.ca

PATENT US-2004-0267175-A1 / CDN PATENT 2,433,816

Copyright, any unauthorized use or copying of the contents or part thereof is prohibited / PRINTED IN CANADA