



Stand-Assist Lift



FIGURE A



FIGURE B



FIGURE C



FIGURE D




FIGURE E

- Provides stable assistance in standing, transferring and toileting
- Four sling attachment points enable the lift to accommodate a variety of stand-up sling designs
- Extra-wide, non-skid foot plate provides a stable base for foot positioning (Figure A)
- Soft, padded adjustable leg support includes a support strap to secure legs to the support (Figure A)
- Easy-roll 4" casters mean stability and easy maneuverability over a variety of surfaces (Figure B)
- 4.5" under-bed clearance provides access under most long term care beds
- 24V DC motor provides power and reliability
- Motor has an emergency stop button and manual lowering feature
- Fast charging, rechargeable battery is removable to eliminate lift downtime
- Battery has a charged and low-charge indicator Audible tone when charge is complete
- Battery charge meter indicates the battery power status while in use
- Audible tone when battery power is below 50%
- Includes removable, rechargeable battery item #13240BATT (Figure C), charging station #13255 (Figure D) and wall-mount bracket item #13260 (Figure E)

PRODUCTS	Description
13246	Stand-Assist-Lift, 1/cs
SPECIFICATIONS	Description
Boom Height	26" - 81"
Base Width	40" (Open) 24.5" (Closed)
Base Clearance	4.5"
Casters	3" (Front) 4" (Rear)
Weight	113 lbs.
Carton Shipping Weight	128 lbs.
Weight Limit	400 lbs.
PARTS	Description
13255	External Charger, 1/ea
13240BATT	Rechargeable Battery, 1/ea
13260	Wall-Mounted Bracket, 1/ea
13240A	Actuator, 1/ea
13240BC	Battery Cell * *Need to order 2 each, 1/ea
13240CAL	Caster with Lock, 1/ea
13240-CC-SL	Cradle Clips Spring-Loaded
13240HC	Hand Control, 1/ea
13240SBH	Spreader Bar Handle, 1/ea
1324601	Sling Hook Screw
STDS2A438-L	Arm Assembly, 1/ea

PARTS	Description
1324603	Sling Bar Foam
1324604	Push Handle Foam, 1/ea
1324605	Knee Pad Assembly
1324606	Knee Pad with Screws
1324607	Knee Pad Strap
1324608	Bolt Caps
1324609	Knee Pad Locking Lever
1324610	Foot Plate
1324611	Knee Pad Bushing
1324612	Spreader Bar
1324613	UpperActuatorMountngPinComplete13246
1324614	Control Box
1324615	KneePadAngleAdjLockingLever
1324616	KneePadAngleAdjustmentRod
1324617	Mast Locking Lever
1324618	Mast Locking Bolt
1324619	Non-Locking Caster
1324620	Base Actuating Rod
1324621	Receiver Plate Assembly
1324622	Base Assembly
1324623	Upper Mast
1324624	Lower Mast
1324625	LowerActuatorMtgPin&SnapRing
1324626	Charger Stand
1324627	Knee Pad Assembly Receiver
13255PCH	Power Cord Hanger
CPCFUSE	Fuse
13255PC	Power Cord for 13255 external charger, 1/ea
1324602	Sling Knob