



Advanced Wound Care Supplies Catalogue

2018-10-26

Category Table of Contents

Absorbent Dressings..... 3	Dressings - Suture / Staple Re- movers 34	Pediatric Catheters.....74
Acute Wounds 7	Dressings - Tapes35	Pleural Drainage.....74
Antimicrobials Absorbent 8	Equipment - Respiratory Aids .35	Protease Inhibitors.....76
Antimicrobials Calcium Alginate 12	Foam Adhesive Dressings38	Respiratory Supplies - Tracheotomy 77
Antimicrobials Charcoal..... 13	Foam Non Adhesive Dressings43	Safety Needles & Syringes78
Antimicrobials Gauze..... 14	Gauze & Non-Woven Dressings 45	Solutions.....90
Antimicrobials Non-Absorbent16	Gloves, Gowns, Masks 53	Swabs and Wipes95
Barrier Wipes/Creams/ Moisturizers 20	Hydrocolloid Dressings57	Tapes99
Calcium Alginate Dressings..... 21	Hydrofibre Dressings..... 60	Total Contact Casting (TCC)...104
Catheters 25	Hydrogels 61	Transparent Films107
Charcoal Dressings..... 27	Incontinence Products 63	Trays & Kits.....111
Compression Bandages 28	IV / Central Lines - Misc 63	Tubular Bandages.....112
Conforming Bandages..... 32	Miscellaneous 63	Venous Support - Pastes118
	Non Adherent Dressing..... 72	Wound Miscellaneous.....118

DUP89545

Absorbent Dressings

Xtrasorb Classic 10cm x 12cm

Manufacturer: Derma Science **Manufacturer Code:** DUP89545

Features: XTRASORB® is a super-absorbent wound dressing. The sachet-style dressing has a non-adherent contact layer to allow for atraumatic removal during dressing changes, and a strike-through protection layer on the outer facing side. Inside this sachet is a super-absorbent polymer fiber core. XTRASORB®'s super-absorbent polymers bind fluid and convert it to a gel. Once fluid enters the dressing and converts to a gel, very little can be re-released, even under compression.

Usage: XTRASORB® "Classic," is for highly exuding wounds when even a foam dressing is not capable of handling the level of fluid. This makes the dressing ideal for use on leg ulcers, where compression bandages can be left in place for up to seven days, as well as on other wounds with copious wound fluid where minimizing dressing changes is a goal.

Precautions:

Size: 10x12cm

Standard Box Quantity: Each



DUP89569

Absorbent Dressings

Xtrasorb Classic 15cm x 22cm

Manufacturer: Derma Science **Manufacturer Code:** DUP89569

Features: XTRASORB® is a super-absorbent wound dressing. The sachet-style dressing has a non-adherent contact layer to allow for atraumatic removal during dressing changes, and a strike-through protection layer on the outer facing side. Inside this sachet is a super-absorbent polymer fiber core. XTRASORB®'s super-absorbent polymers bind fluid and convert it to a gel. Once fluid enters the dressing and converts to a gel, very little can be re-released, even under compression.

Usage: XTRASORB® "Classic," is for highly exuding wounds when even a foam dressing is not capable of handling the level of fluid. This makes the dressing ideal for use on leg ulcers, where compression bandages can be left in place for up to seven days, as well as on other wounds with copious wound fluid where minimizing dressing changes is a goal.

Precautions:

Size: 15x22cm

Standard Box Quantity: Each



MOL285780

Absorbent Dressings

Mesalt 7.5cm x 7.5cm

Manufacturer: Molnlyke **Manufacturer Code:** 285780

Features: Mesalt is made of an absorbent, viscose/polyester nonwoven dressing impregnated with sodium chloride.

Usage: The exudate from the wound releases the sodium chloride from the dressing. Mesalt stimulates the cleansing of heavily draining and discharging wounds in the inflammatory phase by absorbing exudate. Change rate is B.I.D to daily if necessary.

Precautions: Mesalt is not to be used on low exuding or dry wounds.

Size: 7.5x7.5cm

Standard Box Quantity: Each



MOL677001

Absorbent Dressings

Mesorb 10 x 10cm

Manufacturer: Molnlyke **Manufacturer Code:** 677001

Features: Mesorb is an absorbent, low-adherent dressing. It contains an absorbent core of cellulose pulp, sandwiched between a non-woven polypropylene wound contact layer and a textured fluid repellent non-woven backing layer.

Usage: If clinically indicated, the wound should be cleaned in accordance with current nursing practice, and the surrounding skin thoroughly dried. A suitable primary dressing should be placed on the surface of the wound if required and a suitably sized Mesorb applied, ensuring an overlap of at least one centimetre on to the surrounding skin. Changes will normally be determined by the amount of exudate produced.

Precautions: Mesorb should not be used on individuals with a known sensitivity to the dressing or one of its components. The dressing must not be cut, or shaped to fit around drains etc., as this will lead to the loss of cellulose pulp.

Size: 10x10cm

Standard Box Quantity: 1



MOL677201

Absorbent Dressings

Mesorb 10 x 15cm

Manufacturer: Molnlyke **Manufacturer Code:** 677201

Features: Mesorb is an absorbent, low-adherent dressing. It contains an absorbent core of cellulose pulp, sandwiched between a non-woven polypropylene wound contact layer and a textured fluid repellent non-woven backing layer.

Usage: If clinically indicated, the wound should be cleaned in accordance with current nursing practice, and the surrounding skin thoroughly dried. A suitable primary dressing should be placed on the surface of the wound if required and a suitably sized Mesorb applied, ensuring an overlap of at least one centimetre on to the surrounding skin. Changes will normally be determined by the amount of exudate produced.

Precautions: Mesorb should not be used on individuals with a known sensitivity to the dressing or one of its components. The dressing must not be cut, or shaped to fit around drains etc., as this will lead to the loss of cellulose pulp.

Size: 10x15cm

Standard Box Quantity: 1



MOL677401

Absorbent Dressings

Mesorb 10 x 20cm

Manufacturer: Molnlyke **Manufacturer Code:** 677401

Features: Mesorb is an absorbent, low-adherent dressing. It contains an absorbent core of cellulose pulp, sandwiched between a non-woven polypropylene wound contact layer and a textured fluid repellent non-woven backing layer.

Usage: If clinically indicated, the wound should be cleaned in accordance with current nursing practice, and the surrounding skin thoroughly dried. A suitable primary dressing should be placed on the surface of the wound if required and a suitably sized Mesorb applied, ensuring an overlap of at least one centimetre on to the surrounding skin. Changes will normally be determined by the amount of exudate produced.

Precautions: Mesorb should not be used on individuals with a known sensitivity to the dressing or one of its components. The dressing must not be cut, or shaped to fit around drains etc., as this will lead to the loss of cellulose pulp.

Size: 10x20cm

Standard Box Quantity: 1



MOL677601

Absorbent Dressings

Mesorb 15 x 20cm

Manufacturer: Molnlyke **Manufacturer Code:** 677601

Features: Mesorb is an absorbent, low-adherent dressing. It contains an absorbent core of cellulose pulp, sandwiched between a non-woven polypropylene wound contact layer and a textured fluid repellent non-woven backing layer.

Usage: If clinically indicated, the wound should be cleaned in accordance with current nursing practice, and the surrounding skin thoroughly dried. A suitable primary dressing should be placed on the surface of the wound if required and a suitably sized Mesorb applied, ensuring an overlap of at least one centimetre on to the surrounding skin. Changes will normally be determined by the amount of exudate produced.

Precautions: Mesorb should not be used on individuals with a known sensitivity to the dressing or one of its components. The dressing must not be cut, or shaped to fit around drains etc., as this will lead to the loss of cellulose pulp.

Size: 15x20cm

Standard Box Quantity: 1



MOL677701

Absorbent Dressings

Mesorb 20 x 25cm

Manufacturer: Molnlyke **Manufacturer Code:** 677701

Features: Mesorb is an absorbent, low-adherent dressing. It contains an absorbent core of cellulose pulp, sandwiched between a non-woven polypropylene wound contact layer and a textured fluid repellent non-woven backing layer.

Usage: If clinically indicated, the wound should be cleaned in accordance with current nursing practice, and the surrounding skin thoroughly dried. A suitable primary dressing should be placed on the surface of the wound if required and a suitably sized Mesorb applied, ensuring an overlap of at least one centimetre on to the surrounding skin. Changes will normally be determined by the amount of exudate produced.

Precautions: Mesorb should not be used on individuals with a known sensitivity to the dressing or one of its components. The dressing must not be cut, or shaped to fit around drains etc., as this will lead to the loss of cellulose pulp.

Size: 20x25cm

Standard Box Quantity: 1



MOL677801

Absorbent Dressings

Mesorb 20 x 30cm

Manufacturer: Molnlyke **Manufacturer Code:** 677801

Features: Mesorb is an absorbent, low-adherent dressing. It contains an absorbent core of cellulose pulp, sandwiched between a non-woven polypropylene wound contact layer and a textured fluid repellent non-woven backing layer.

Usage: If clinically indicated, the wound should be cleaned in accordance with current nursing practice, and the surrounding skin thoroughly dried. A suitable primary dressing should be placed on the surface of the wound if required and a suitably sized Mesorb applied, ensuring an overlap of at least one centimetre on to the surrounding skin. Changes will normally be determined by the amount of exudate produced.

Precautions: Mesorb should not be used on individuals with a known sensitivity to the dressing or one of its components. The dressing must not be cut, or shaped to fit around drains etc., as this will lead to the loss of cellulose pulp.

Size: 20x30cm

Standard Box Quantity: 1



SN5999004120

Absorbent Dressings

Exu-Dry 10 x 15cm

Manufacturer: Smith & Nephew **Manufacturer Code:** 5999004120

Features: Exudry is a unique one-piece wound dressing designed to cover the wound with reduced friction.

Usage: Indicated for exudate absorption and the management of partial to full-thickness wounds.

Precautions: If the wound continues to get larger after the first few dressing changes, consult an appropriate healthcare professional.

Size: 10 x 15cm

Standard Box Quantity: 1



SN5999006

Absorbent Dressings

Exudry 15 x 23cm

Manufacturer: Smith & Nephew **Manufacturer Code:** 5999006

Features: Exudry is a unique one-piece wound dressing designed to cover the wound with reduced friction.

Usage: Indicated for exudate absorption and the management of partial to full-thickness wounds.

Precautions: If the wound continues to get larger after the first few dressing changes, consult an appropriate healthcare professional.

Size: 15x23cm

Standard Box Quantity: 1



SN5999009

Absorbent Dressings

Exudry 23 x 38cm

Manufacturer: Smith & Nephew **Manufacturer Code:** 5999009

Features: Exudry is a unique one-piece wound dressing designed to cover the wound with reduced friction.

Usage: Indicated for exudate absorption and the management of partial to full-thickness wounds.

Precautions: If the wound continues to get larger after the first few dressing changes, consult an appropriate healthcare professional.

Size: 23x38cm

Standard Box Quantity: 1



MOL670820

Acute Wounds

Mepore Pro 6 x 7cm

Manufacturer: Molnlyke **Manufacturer Code:** 670820

Features: Mepore Pro is a self-adhesive, absorbent and semi-permeable dressing with an outertransparent film layer to protect the wound from water and contamination. Its skin-friendly, water-based, solvent-free polyacrylate adhesive provides secure yet gentle fixation. It is bacterial proof, viral proof and waterproof to protect the patient and the clinician.

Usage: Mepore Pro is appropriate for wounds with low to moderate exudate levels. These include surgical incisions, sutured wounds, Laproscopy and Arthroscopy puncture sites, abrasions, traumatic wounds and cannula sites.

Precautions: Care should be taken to ensure that Mepore Pro is not applied under tension, to prevent shearing forces causing damage to the skin. This is particularly important when the dressing is applied over joints. It should not be applied to patients who are sensitive to acrylic adhesives or to individuals who have very fragile or easily damaged skin.

Size: 6 x 7cm

Standard Box Quantity: 1



MOL670920

Acute Wounds

Mepore Pro 9 X 10cm

Manufacturer: Molnlyke **Manufacturer Code:** 670920

Features: Mepore Pro is a self-adhesive, absorbent and semi-permeable dressing with an outertransparent film layer to protect the wound from water and contamination. Its skin-friendly, water-based, solvent-free polyacrylate adhesive provides secure yet gentle fixation. It is bacterial proof, viral proof and waterproof to protect the patient and the clinician.

Usage: Mepore Pro is appropriate for wounds with low to moderate exudate levels. These include surgical incisions, sutured wounds, Laproscopy and Arthroscopy puncture sites, abrasions, traumatic wounds and cannula sites.

Precautions: Care should be taken to ensure that Mepore Pro is not applied under tension, to prevent shearing forces causing damage to the skin. This is particularly important when the dressing is applied over joints. It should not be applied to patients who are sensitive to acrylic adhesives or to individuals who have very fragile or easily damaged skin.

Size: 9 X 10cm

Standard Box Quantity: 1



MOL671220

Acute Wounds

Mepore Pro 9 x 25cm

Manufacturer: Molnlyke **Manufacturer Code:** 671220

Features: Mepore Pro is a self-adhesive, absorbent and semi-permeable dressing with an outertransparent film layer to protect the wound from water and contamination. Its skin-friendly, water-based, solvent-free polyacrylate adhesive provides secure yet gentle fixation. It is bacterial proof, viral proof and waterproof to protect the patient and the clinician.

Usage: Mepore Pro is appropriate for wounds with low to moderate exudate levels. These include surgical incisions, sutured wounds, Laproscopy and Arthroscopy puncture sites, abrasions, traumatic wounds and cannula sites.

Precautions: Care should be taken to ensure that Mepore Pro is not applied under tension, to prevent shearing forces causing damage to the skin. This is particularly important when the dressing is applied over joints. It should not be applied to patients who are sensitive to acrylic adhesives or to individuals who have very fragile or easily damaged skin.

Size: 9 x 25cm

Standard Box Quantity: 1



COL9625

Antimicrobials Absorbent

Biatain Foam AG 15 X 15cm Non Adhesive**Manufacturer:** Coloplast **Manufacturer Code:** 9625

Features: Non-Adhesive Biatain Ag dressings are prepared with hydro-activated silver, which is released into the wound during wear. This small innovation has shown excellent results: infection in and around the wound is reduced by 50 per cent within four weeks, and the odour fades or completely disappears in just one week. Biatain Ag gives effective protection against bacteria and provides the best wound bed conditions for healing acute and chronic wounds.

Usage: Biatain Ag is suitable for the treatment of leg ulcers, pressure ulcers, diabetic foot ulcers, smaller second-degree burns, donor sites, postoperative wounds, skin abrasions. Biatain Ag releases silver ions according to the amount of exudate present in the wound. The more exudate the more silver is released. The release persists for up to 7 days.

Precautions: A healthcare professional should frequently inspect and manage diabetic wounds and wounds, which are solely or partially caused by arterial insufficiency, in accordance with local standards.

Size: 15x15cm**Standard Box Quantity:** 1

CON403708


Antimicrobials Absorbent

Aquacel AG Hydrofiber 10 x 10cm**Manufacturer:** ConvaTec **Manufacturer Code:** 403708

Features: AQUACEL® Ag dressing incorporating unique Hydrofiber® Technology with 1.2% (w/w) silver combines the favorable gelling characteristics of Hydrofiber® Technology with the broad-spectrum antimicrobial properties of ionic silver (Ag+).

Usage: It is a versatile primary dressing indicated for use on moderately and highly exuding chronic and acute wounds where there is an infection* or an increased risk of infection.

Precautions: Should not be used on individuals that have a known hypersensitivity to silver.

Size: 10x10cm**Standard Box Quantity:** 1

CON403710

Antimicrobials Absorbent

Aquacel AG Hydrofiber 14.5 x 14.5cm**Manufacturer:** ConvaTec **Manufacturer Code:** 403710

Features: AQUACEL® Ag dressing incorporating unique Hydrofiber® Technology with 1.2% (w/w) silver combines the favorable gelling characteristics of Hydrofiber® Technology with the broad-spectrum antimicrobial properties of ionic silver (Ag+).

Usage: It is a versatile primary dressing indicated for use on moderately and highly exuding chronic and acute wounds where there is an infection* or an increased risk of infection.

Precautions: Should not be used on individuals that have a known hypersensitivity to silver.

Size: 14.5x14.5cm**Standard Box Quantity:** 1

CON403771

Antimicrobials Absorbent

Aquacel AG Hydrofiber Ribbon/Rope 2 x 45cm

Manufacturer: ConvaTec **Manufacturer Code:** 403771

Features: AQUACEL® Ag dressing incorporating unique Hydrofiber® Technology with 1.2% (w/w) silver combines the favorable gelling characteristics of Hydrofiber® Technology with the broad-spectrum antimicrobial properties of ionic silver (Ag+).

Usage: It is a versatile primary dressing indicated for use on moderately and highly exuding chronic and acute wounds where there is an infection* or an increased risk of infection.

Precautions: Should not be used on individuals that have a known hypersensitivity to silver.

Size: 2x45cm

Standard Box Quantity: 1



CON420642

Antimicrobials Absorbent

Aquacel AG Foam Dressing (Non-Adhesive) 10 X 10cm

Manufacturer: ConvaTec **Manufacturer Code:** 420642

Features: AQUACEL® Ag Foam dressing offers the healing benefits of a Hydrofiber® interface layer, the comfort of foam and the bacteria-killing power of ionic silver.

Usage: Management of both chronic and acute wounds, such as: leg ulcers, pressure ulcers (stage II-IV) and diabetic ulcers; surgical wounds (post-operative, donor sites, dermatological); partial-thickness (second-degree) burns; traumatic or surgical wounds left to heal by secondary intention such as dehiscent surgical incisions; surgical wounds that heal by primary intent, such as dermatological and surgical incisions (e.g., orthopedic and vascular); local management of wounds that are prone to bleeding, such as wounds that have been mechanically or surgically debrided and donor sites; painful wounds; abrasions; lacerations; minor cuts; minor scalds and burns.

Precautions: AQUACEL® Ag Foam dressings should not be used on individuals who are sensitive to or who have had an allergic reaction to the dressing or its components.

Size: 10 X 10cm

Standard Box Quantity: 1



CON420681

Antimicrobials Absorbent

Aquacel AG Foam Dressing ADHESIVE 10 X 10cm

Manufacturer: ConvaTec **Manufacturer Code:** 420681

Features: AQUACEL® Ag Foam dressing offers the healing benefits of a Hydrofiber® interface layer, the comfort of foam and the bacteria-killing power of ionic silver.

Usage: Management of both chronic and acute wounds, such as: leg ulcers, pressure ulcers (stage II-IV) and diabetic ulcers; surgical wounds (post-operative, donor sites, dermatological); partial-thickness (second-degree) burns; traumatic or surgical wounds left to heal by secondary intention such as dehiscent surgical incisions; surgical wounds that heal by primary intent, such as dermatological and surgical incisions (e.g., orthopedic and vascular); local management of wounds that are prone to bleeding, such as wounds that have been mechanically or surgically debrided and donor sites; painful wounds; abrasions; lacerations; minor cuts; minor scalds and burns.

Precautions: AQUACEL® Ag Foam dressings should not be used on individuals who are sensitive to or who have had an allergic reaction to the dressing or its components.

Size: 10 X 10cm

Standard Box Quantity: 1



HOLBT0906

Antimicrobials Absorbent

Hydrofera Blue Classic Tunnel Dressing 9mm x 14cm

Manufacturer: Hollister **Manufacturer Code:** HBT0906

Features: Provides a protective antibacterial tunneling dressing that inhibits the growth of microorganisms. Physically binds harmful endotoxins which can aid in patient comfort. Effective against MRSA, VRE and Candida. Highly Absorptive - pulls bacteria-laden exudate up and away from the wound which may facilitate healing and aid in patient comfort. Similar absorbency to other highly absorptive dressings.

Usage: Intended as an external dressing for pressure ulcers, venous stasis ulcers, arterial ulcers, donor sites, abrasions, lacerations, superficial burns, postsurgical incisions, diabetic ulcers, and other wounds caused by trauma.

Precautions: Compatible for use with enzymatic debridement agents and growth factors. Not indicated for third-degree burns.

Size: 9mm x 14cm (2-pk)

Standard Box Quantity: Pk/2



HOLHB2214

Antimicrobials Absorbent

Hydrofera Blue Classic Dressing 5 x 5cm

Manufacturer: Hollister **Manufacturer Code:** HB2214

Features: Provides a protective antibacterial cover that inhibits the growth of microorganisms. Physically binds harmful endotoxins which can aid in patient comfort. Effective against MRSA, VRE and Candida. Highly Absorptive - pulls bacteria-laden exudate up and away from the wound which may facilitate healing and aid in patient comfort. Similar absorbency to other highly absorptive dressings.

Usage: Intended as an external dressing for pressure ulcers, venous stasis ulcers, arterial ulcers, donor sites, abrasions, lacerations, superficial burns, postsurgical incisions, diabetic ulcers, and other wounds caused by trauma.

Precautions: Compatible for use with enzymatic debridement agents and growth factors. Not indicated for third-degree burns.

Size: 5cm x 5cm

Standard Box Quantity: 1



HOLHB4414

Antimicrobials Absorbent

Hydrofera Blue Classic Dressing 10 x 10cm

Manufacturer: Hollister **Manufacturer Code:** HB4414CA

Features: Provides a protective antibacterial cover that inhibits the growth of microorganisms. Physically binds harmful endotoxins which can aid in patient comfort. Effective against MRSA, VRE and Candida. Highly Absorptive - pulls bacteria-laden exudate up and away from the wound which may facilitate healing and aid in patient comfort. Similar absorbency to other highly absorptive dressings.

Usage: Intended as an external dressing for pressure ulcers, venous stasis ulcers, arterial ulcers, donor sites, abrasions, lacerations, superficial burns, postsurgical incisions, diabetic ulcers, and other wounds caused by trauma.

Precautions: Compatible for use with enzymatic debridement agents and growth factors. Not indicated for third-degree burns.

Size: 10cm x 10cm

Standard Box Quantity: 1



Hydrofera Blue Antimicrobial Dressing 15x15cm**Manufacturer:** Hollister **Manufacturer Code:** HB6614

Features: Provides a protective antibacterial cover that inhibits the growth of microorganisms. Physically binds harmful endotoxins which can aid in patient comfort. Effective against MRSA, VRE and Candida. Highly Absorptive - pulls bacteria-laden exudate up and away from the wound which may facilitate healing and aid in patient comfort. Similar absorbency to other highly absorptive dressings.

Usage: Intended as an external dressing for pressure ulcers, venous stasis ulcers, arterial ulcers, donor sites, abrasions, lacerations, superficial burns, postsurgical incisions, diabetic ulcers, and other wounds caused by trauma.

Precautions: Compatible for use with enzymatic debridement agents and growth factors. Not indicated for third-degree burns.

Size: 15cm x 15cm**Standard Box Quantity:** Each**Iodosorb Ointment 10gm (tube)****Manufacturer:** Smith & Nephew **Manufacturer Code:** 66060630

Features: Iodosorb ointment contains 50% cadexomer iodine and 50% ointment base. Iodosorb changes colour from brown to yellow when the dressing requires changing.

Usage: Iodosorb* is used for the treatment of chronic leg ulcers. It helps promote wound healing and patient comfort.

Precautions: Should not be used in patients with known or suspected iodine sensitivity, Hashimoto's thyroiditis or prior history of Graves' Disease, or in the case of non-toxic nodular goitre. Should not be used pregnant or lactating women.

Size: 10gm**Standard Box Quantity:** 1**Inadine 5 x 5cm****Manufacturer:** Acelyt **Manufacturer Code:** P01481

Features: INADINE* Dressing is a topical wound dressing impregnated with an ointment containing 10% Povidone iodine (PVP-1). The dressing also contains polyethylene glycol and purified water. Povidone Iodine controls bacteria by reducing bacterial load and inhibits infection.

Usage: INADINE Dressing is designed to protect the wound, even if infected. INADINE dressing is indicated for the management of ulcerative wounds and may also be used for the prevention of infection in minor burns and traumatic skin loss injuries.

Precautions: Safe for use on adults and children and in combination with systemic antibiotics.

Size: 5x5cm**Standard Box Quantity:** 1

Inadine 9.5 x 9.5cm**Manufacturer:** Acelity **Manufacturer Code:** P01512

Features: INADINE® Dressing is a topical wound dressing impregnated with an ointment containing 10% Povidone iodine (PVP-1). The dressing also contains polyethylene glycol and purified water. Povidone Iodine controls bacteria by reducing bacterial load and inhibits infection.

Usage: INADINE Dressing is designed to protect the wound, even if infected. INADINE dressing is indicated for the management of ulcerative wounds and may also be used for the prevention of infection in minor burns and traumatic skin loss injuries.

Precautions: Safe for use on adults and children and in combination with systemic antibiotics.

Size: 9.5x9.5cm**Standard Box Quantity:** 1

COL3760


Antimicrobials Calcium Alginate

Biatain Alginate AG Dressing 10X10cm**Manufacturer:** Coloplast **Manufacturer Code:** 3760

Features: A highly absorbent antimicrobial pad which releases silver ions in the presence of wound exudate as an effective broad spectrum antimicrobial agent and aid in odour reduction. As wound exudate is absorbed, alginate forms a gel, which assists in maintaining a moist environment for optimal wound healing, and allows intact removal.

Usage: Indicated for moderate to heavily exuding wounds such as pressure ulcers, leg ulcers, diabetic ulcers, second-degree burns, grafts, donor sites or traumatic wounds. May be used on infected wounds under the discretion of a healthcare professional. Can be used for cavity wounds where it should be gently layered into the wound, as the gel will conform to the wound shape.

Precautions: Do not use for clients with known sensitivity or allergy to silver or alginates.

Size: 10X10cm**Standard Box Quantity:** 1

COL3780

Antimicrobials Calcium Alginate

Biatain Alginate AG Ribbon Rope 3x44cm**Manufacturer:** Coloplast **Manufacturer Code:** 3780

Features: A highly absorbent antimicrobial pad which releases silver ions in the presence of wound exudate as an effective broad spectrum antimicrobial agent and aid in odour reduction. As wound exudate is absorbed, alginate forms a gel, which assists in maintaining a moist environment for optimal wound healing, and allows intact removal.

Usage: For highly exuding, deep cavity wounds. May be used on infected wounds under the discretion of a healthcare professional. Can be used for cavity wounds where it should be gently layered into the wound, as the gel will conform to the wound shape.

Precautions: Do not use for clients with known sensitivity or allergy to silver or alginates.

Size: 3x44cm**Standard Box Quantity:** 1

SYSMAS065

Antimicrobials Charcoal

Actisorb Silver 220 Activated Charcoal Dressing 9.5cm x 6.5 cm

Manufacturer: Acelity **Manufacturer Code:** MAS065

Features: Consists principally of activated carbon impregnated with metallic silver produced by a heating a specially treated fine viscose fabric under carefully controlled conditions.

Usage: Actisorb Silver 220 may be used in the treatment of most types of chronic wounds but is particularly recommended for the management of malodorous, infected wounds, including fungating lesions, faecal fistulae, infected pressure sores and heavily exuding leg ulcers.

Precautions: Actisorb Silver 220 should not be used on patients who are sensitive to nylon, and should be used with care as a primary dressing on wounds that have a tendency to dry out. In these situations adherence of the dressing to the surface of the wound may become a problem.

Size: 9.5x6.5cm

Standard Box Quantity: 1



SYSMAS105

Antimicrobials Charcoal

Actisorb Silver 220 Activated Charcoal Dressing 10.5cm x 10.5 cm

Manufacturer: Acelity **Manufacturer Code:** MAS105

Features: Consists principally of activated carbon impregnated with metallic silver produced by a heating a specially treated fine viscose fabric under carefully controlled conditions.

Usage: Actisorb Silver 220 may be used in the treatment of most types of chronic wounds but is particularly recommended for the management of malodorous, infected wounds, including fungating lesions, faecal fistulae, infected pressure sores and heavily exuding leg ulcers.

Precautions: Actisorb Silver 220 should not be used on patients who are sensitive to nylon, and should be used with care as a primary dressing on wounds that have a tendency to dry out. In these situations adherence of the dressing to the surface of the wound may become a problem.

Size: 10.5x10.5 cm

Standard Box Quantity: 1



SYSMAS190

Antimicrobials Charcoal

Actisorb Silver 220 Activated Charcoal Dressing 19.5cm x 10.5cm

Manufacturer: Acelity **Manufacturer Code:** MAS190

Features: Consists principally of activated carbon impregnated with metallic silver produced by a heating a specially treated fine viscose fabric under carefully controlled conditions.

Usage: Actisorb Silver 220 may be used in the treatment of most types of chronic wounds but is particularly recommended for the management of malodorous, infected wounds, including fungating lesions, faecal fistulae, infected pressure sores and heavily exuding leg ulcers.

Precautions: Actisorb Silver 220 should not be used on patients who are sensitive to nylon, and should be used with care as a primary dressing on wounds that have a tendency to dry out. In these situations adherence of the dressing to the surface of the wound may become a problem.

Size: 19x10.5cm

Standard Box Quantity: 1



COV2506

Antimicrobials Gauze

AMD Gauze 5 x 5cm

Manufacturer: Covidien **Manufacturer Code:** 2506

Features: Kerlix™ AMD antimicrobial dressings super sponges are impregnated with PHMB (Polyhexamethylene Biguanide), a powerful yet safe antiseptic that has a broad range effectiveness against gram positive and gram negative microorganisms including some multi-drug resistant strains such as MRSA.

Usage: Can be used as a primary dressing to minimize bacterial log counts in the wound bed. Can be used as a secondary dressing to eliminate growth within the gauze and protect against bacterial proliferation through the gauze.

Precautions: Dakin's Solution™* may inactivate PHMB. This will impact dressing efficacy but should not be dangerous to patient health.

Size: 5x5cm

Standard Box Quantity: PKG/2



COV2533

Antimicrobials Gauze

AMD Gauze 10 x 10cm

Manufacturer: Covidien **Manufacturer Code:** 2533

Features: Kerlix™ AMD antimicrobial dressings super sponges are impregnated with PHMB (Polyhexamethylene Biguanide), a powerful yet safe antiseptic that has a broad range effectiveness against gram positive and gram negative microorganisms including some multi-drug resistant strains such as MRSA.

Usage: Can be used as a primary dressing to minimize bacterial log counts in the wound bed. Can be used as a secondary dressing to eliminate growth within the gauze and protect against bacterial proliferation through the gauze.

Precautions: Dakin's Solution™* may inactivate PHMB. This will impact dressing efficacy but should not be dangerous to patient health.

Size: 10x10cm

Standard Box Quantity: PKG/2



COV3332

Antimicrobials Gauze

AMD Gauze 11.4cm x 3.7m

Manufacturer: Covidien **Manufacturer Code:** 3332

Features: Kerlix™ AMD rolls are impregnated with PHMB (Polyhexamethylene Biguanide), a powerful yet safe antiseptic that has a broad range effectiveness against gram positive and gram negative microorganisms including some multi-drug resistant strains such as MRSA. In addition to its antimicrobial capabilities, Kerlix™ AMD rolls provide fast-wicking action, superior aeration, and excellent absorbency.

Usage: Can be used as a primary dressing to minimize bacterial log counts in the wound bed. Can be used as a secondary dressing to eliminate growth within the gauze and protect against bacterial proliferation through the gauze.

Precautions: Dakin's Solution™* may inactivate PHMB. This will impact dressing efficacy but should not be dangerous to patient health.

Size: 11.4cm x 3.7m roll

Standard Box Quantity: 1



COV7831AMD

Antimicrobials Gauze

AMD Ribbon Gauze.6cm (1/4 inchx1yrd)

Manufacturer: Covidien **Manufacturer Code:** 7831AMD

Features: AMD packing strips are impregnated with PHMB (Polyhexamethylene Biguanide), a powerful yet safe antiseptic that has a broad range effectiveness against gram positive and gram negative microorganisms including some multi-drug resistant strains such as MRSA.

Usage: AMD antimicrobial dressings resist bacterial colonization, reduce bacterial penetration, limit cross-contamination, and maintain balance of native skin flora.

Precautions: Dakin's Solution™ may inactivate PHMB. This will impact dressing efficacy but should not be dangerous to patient health.

Size: 1/4 inch x 1yd

Standard Box Quantity: 1



COV7832AMD

Antimicrobials Gauze

AMD Ribbon Gauze 1.3cm (1/2 inchx1yrd)

Manufacturer: Covidien **Manufacturer Code:** 7832AMD

Features: AMD packing strips are impregnated with PHMB (Polyhexamethylene Biguanide), a powerful yet safe antiseptic that has a broad range effectiveness against gram positive and gram negative microorganisms including some multi-drug resistant strains such as MRSA.

Usage: AMD antimicrobial dressings resist bacterial colonization, reduce bacterial penetration, limit cross-contamination, and maintain balance of native skin flora.

Precautions: Dakin's Solution™ may inactivate PHMB. This will impact dressing efficacy but should not be dangerous to patient health.

Size: 1/2 inch x 1yd

Standard Box Quantity: 1



COV7833AMD

Antimicrobials Gauze

AMD Ribbon Gauze 2.5cm (1 inchx1yrd)

Manufacturer: Covidien **Manufacturer Code:** 7833AMD

Features: AMD packing strips are impregnated with PHMB (Polyhexamethylene Biguanide), a powerful yet safe antiseptic that has a broad range effectiveness against gram positive and gram negative microorganisms including some multi-drug resistant strains such as MRSA.

Usage: AMD antimicrobial dressings resist bacterial colonization, reduce bacterial penetration, limit cross-contamination, and maintain balance of native skin flora.

Precautions: Dakin's Solution™ may inactivate PHMB. This will impact dressing efficacy but should not be dangerous to patient health.

Size: 1 inch x 1yd

Standard Box Quantity: 1



SM405

Antimicrobials Non-Absorbent

Medihoney Gel 50g

Manufacturer: Derma Science **Manufacturer Code:** DUP405

Features: Medihoney dressings have been shown to have powerful and long-lasting antimicrobial activity – similar to that of ionic silver – but without the associated toxicity of ionic silver.

Usage: For use in wounds with light to moderate exudate and hard to dress areas including tunneled or undermined wounds, or sinus tracts.

Precautions: Due to the high osmotic activity of MEDIHONEY® dressings, when appropriate it is recommended to protect the skin with a skin protectant, and to initially cover the dressing with a highly absorbent secondary dressing.

Size: 50 gm tube

Standard Box Quantity: Each



SM793

Antimicrobials Non-Absorbent

Medihoney Calcium Alginate Rope 1.9 x 30.5cm

Manufacturer: Derma Science **Manufacturer Code:** DUP793

Features: Medihoney Calcium Alginate is a honey-based dressing for moderately to heavily exuding wounds. The dressing contains 95% Active Leptospermum honey and can be used to cover a shallow wound or fill a deep wound, while providing a moist environment conducive to wound healing. Leptospermum Honey assists in promoting an optimal healing environment. It has a low pH which is beneficial to chronic wounds.

Usage: For use in superficial to deep wounds with moderate to heavy exudate.

Precautions: Latex free.

Size: 1.9cm x 30cm

Standard Box Quantity: 1



SM794

Antimicrobials Non-Absorbent

Medihoney Calcium Alginate 5 x 5cm

Manufacturer: Derma Science **Manufacturer Code:** DUP794

Features: Medihoney Calcium Alginate is a honey-based dressing for moderately to heavily exuding wounds. The dressing contains 95% Active Leptospermum honey and can be used to cover a shallow wound or fill a deep wound, while providing a moist environment conducive to wound healing. Leptospermum Honey assists in promoting an optimal healing environment. It has a low pH which is beneficial to chronic wounds.

Usage: For use in superficial to deep wounds with moderate to heavy exudate.

Precautions: Latex free.

Size: 5cm x 5cm

Standard Box Quantity: 1



Medihoney Calcium Alginate 10 x 10cm**Manufacturer:** Derma Science **Manufacturer Code:** DUP795

Features: Medihoney Calcium Alginate is a honey-based dressing for moderately to heavily exuding wounds. The dressing contains 95% Active Leptospermum honey and can be used to cover a shallow wound or fill a deep wound, while providing a moist environment conducive to wound healing. Leptospermum Honey assists in promoting an optimal healing environment. It has a low pH which is beneficial to chronic wounds.

Usage: For use in superficial to deep wounds with moderate to heavy exudate.

Precautions: Latex free.

Size: 10cm x 10cm

Standard Box Quantity: 1

**Bactigras 10 x 10cm****Manufacturer:** Smith & Nephew **Manufacturer Code:** 66003650

Features: Bactigras* is a medicated paraffin gauze dressing and features soft paraffin base, sterile leno weave presentation and contains Chlorhexidine acetate 0.5%, an antiseptic with a broad spectrum.

Usage: Bactigras is indicated for adjunctive treatment and prevention of infection in skin loss lesions, including wounds, burns and ulcers.

Precautions: Due to the cationic properties of chlorhexidine, Bactigras is incompatible with anionic surfactants and other anionic compounds. One layer of tulle gras only should be applied directly to the surface of lesions being treated, and this should be done with forceps under sterile conditions. Avoid contact with eyes, the middle ear, meninges and brain. Bactigras should not be used if the seal or wrapper is broken.

Size: 10x10cm

Standard Box Quantity: 1

**Acticoat Flex 7 (Up to 7 days) 5 x 5cm****Manufacturer:** Smith & Nephew **Manufacturer Code:** 66800395

Features: The ACTICOAT™ FLEX 7 dressing consists of a flexible, low adherent polyester layer coated with nanocrystalline silver.

Usage: Acticoat FLEX 7 is indicated for use on partial and full thickness wounds. This includes First and second degree burns, Covering of grafts, Surgical sites, Venous ulcers, Pressure ulcers and Diabetic ulcers.

Precautions: For external use only, not compatible with oil-based products such as petrolatum, Avoid contact with electrodes or conductive gels during electronic measurements, not intended for use on 3rd degree burns and should only be used on premature infants (less than 37 weeks gestation) when the clinical benefit outweighs potential risks.

Size: 5 x 5cm

Standard Box Quantity: 1



SN66800397

Antimicrobials Non-Absorbent

Acticoat Flex 7 (Up to 7 days) 10 x 10cm

Manufacturer: Smith & Nephew **Manufacturer Code:** 66800397

Features: The ACTICOAT™ FLEX 7 dressing consists of a flexible, low adherent polyester layer coated with nanocrystalline silver.

Usage: Acticoat FLEX 7 is indicated for use on partial and full thickness wounds. This includes First and second degree burns, Covering of grafts, Surgical sites, Venous ulcers, Pressure ulcers and Diabetic ulcers.

Precautions: For external use only, not compatible with oil-based products such as petrolatum, Avoid contact with electrodes or conductive gels during electronic measurements, not intended for use on 3rd degree burns and should only be used on premature infants (less than 37 weeks gestation) when the clinical benefit outweighs potential risks.

Size: 10 x 10cm

Standard Box Quantity: 1



SN66800399

Antimicrobials Non-Absorbent

Acticoat Flex 3 10 x 10 cm

Manufacturer: Smith & Nephew **Manufacturer Code:** 66800399

Features: The ACTICOAT™ FLEX 3 dressing consists of a flexible, low adherent polyester layer coated with nanocrystalline silver.

Usage: Acticoat Flex 3 is indicated for use on partial and full thickness wounds. This includes First and second degree burns, Covering of grafts, Surgical sites, Venous ulcers, Pressure ulcers and Diabetic ulcers.

Precautions: For external use only, not compatible with oil-based products such as petrolatum, Avoid contact with electrodes or conductive gels during electronic measurements, not intended for use on 3rd degree burns and should only be used on premature infants (less than 37 weeks gestation) when the clinical benefit outweighs potential risks.

Size: 10x10cm

Standard Box Quantity: Each



SN66800420

Antimicrobials Non-Absorbent

Acticoat Flex 7 (up to 7 days) 10 x 20cm

Manufacturer: Smith & Nephew **Manufacturer Code:** 66800420

Features: The ACTICOAT™ FLEX 7 dressing consists of a flexible, low adherent polyester layer coated with nanocrystalline silver.

Usage: Acticoat FLEX 7 is indicated for use on partial and full thickness wounds. This includes First and second degree burns, Covering of grafts, Surgical sites, Venous ulcers, Pressure ulcers and Diabetic ulcers.

Precautions: For external use only, not compatible with oil-based products such as petrolatum, Avoid contact with electrodes or conductive gels during electronic measurements, not intended for use on 3rd degree burns and should only be used on premature infants (less than 37 weeks gestation) when the clinical benefit outweighs potential risks.

Size: 10 x 20cm

Standard Box Quantity: 1



SN7456

Antimicrobials Non-Absorbent

Bactigras 5 x 5cm

Manufacturer: Smith & Nephew **Manufacturer Code:** 7456

Features: Bactigras* is a medicated paraffin gauze dressing and features soft paraffin base, sterile leno weave presentation and contains Chlorhexidine acetate 0.5%, an antiseptic with a broad spectrum.

Usage: Bactigras is indicated for adjunctive treatment and prevention of infection in skin loss lesions, including wounds, burns and ulcers.

Precautions: Due to the cationic properties of chlorhexidine, Bactigras is incompatible with anionic surfactants and other anionic compounds. One layer of tulle gras only should be applied directly to the surface of lesions being treated, and this should be done with forceps under sterile conditions. Avoid contact with eyes, the middle ear, meninges and brain. Bactigras should not be used if the seal or wrapper is broken.

Size: 5x5cm

Standard Box Quantity: 1



SN7461

Antimicrobials Non-Absorbent

Bactigras 15 x 20cm

Manufacturer: Smith & Nephew **Manufacturer Code:** 7461

Features: Bactigras* is a medicated paraffin gauze dressing and features soft paraffin base, sterile leno weave presentation and contains Chlorhexidine acetate 0.5%, an antiseptic with a broad spectrum.

Usage: Bactigras is indicated for adjunctive treatment and prevention of infection in skin loss lesions, including wounds, burns and ulcers.

Precautions: Due to the cationic properties of chlorhexidine, Bactigras is incompatible with anionic surfactants and other anionic compounds. One layer of tulle gras only should be applied directly to the surface of lesions being treated, and this should be done with forceps under sterile conditions. Avoid contact with eyes, the middle ear, meninges and brain. Bactigras should not be used if the seal or wrapper is broken.

Size: 15x20cm

Standard Box Quantity: 1



SYSCAD7011

Antimicrobials Non-Absorbent

Silvercel Non Adherent 11cm X 11cm

Manufacturer: Acelyt **Manufacturer Code:** CAD7011

Features: SILVERCEL™ Hydro-Alginate Dressing with Silver is a sterile, non-woven pad composed of a high G (guluronic acid) alginate, carboxy methyl cellulose (CMC) and silver coated nylon fibres. This unique composition manages exudate in infected or heavily colonised wounds. The silver fibres combat wound microorganisms, which reduces the bacterial colonization of the wound and inhibits infection - providing a powerful antimicrobial solution.

Usage: SILVERCEL™ Hydro-Alginate Antimicrobial Dressing with Silver is intended for use in the management of all moderate to heavily exuding partial and full thickness chronic wounds including: decubitus (pressure) ulcers, venous ulcers, diabetic ulcers, donor sites, traumatic and surgical wounds.

Precautions: SILVERCEL™ Dressing is not indicated for use on the following: Third-degree burns, Patients with a known sensitivity to alginates or silver, To control heavy bleeding.

Size: 11x11cm

Standard Box Quantity: Each



Cavilon No Sting Barrier Spray 28ml btl**Manufacturer:** 3M **Manufacturer Code:** 3346E

Features: A fast-drying, non-sticky, alcohol-free liquid barrier film that forms a breathable, transparent coating on the skin. Protects damaged or intact skin from body fluids, adhesive trauma, friction, and incontinence. The film is hypoallergenic, non-cytotoxic and will not sting even when applied to damaged or denuded skin.

Usage: A hypoallergenic, non-cytotoxic liquid film which will not sting even when applied to damaged or denuded skin.

Precautions: If Cavilon No Sting Barrier Film is applied to an area with skin folds or other skin-to-skin contact, ensure that these areas of skin are separated to allow coating to dry thoroughly before allowing skin to return to normal position.

Size: 28ml**Standard Box Quantity:** 1**Sween Cream 190ml****Manufacturer:** Coloplast **Manufacturer Code:** 7068

Features: The non-occlusive moisturizers and natural vitamins A&D, immediately help restore dry, fragile skin that's at risk for cracking, especially around contact points like heels and elbows. Non-occlusive formula absorbs quickly into the skin and will not interfere with the adhesion of tape or ostomy appliances.

Usage: Use Sween Cream® for: Cracked and flaky dry skin; Red skin due to minor burns; Radiation/chemotherapy skin reactions; Peri-wound skin; Perineal Erythema; Surgical incisions - reduces itching and pulling discomfort; Relieves itching under tape; Callouses; Any reddened skin.

Precautions: For external use only. If irritation occurs, consult a physician.

Size: 190ml**Standard Box Quantity:** 1**Critic Aid Cream Clear 71gm****Manufacturer:** Coloplast **Manufacturer Code:** 7573

Features: A clear petrolatum based moisture barrier ointment containing dimethicone. Absorbs skin moisture allowing it to adhere to denuded skin.

Usage: For light to heavy skin irritation due to incontinence.

Precautions: Do not use this medication if you are allergic to zinc, dimethicone, lanolin, cod liver oil, petroleum jelly, parabens, mineral oil, or wax. Zinc oxide topical will not treat a bacterial or fungal infection. Avoid getting this medication in your mouth or eyes.

Size: 71g**Standard Box Quantity:** 1


SN59431279

Barrier Wipes/Creams/Moisturizers

Secura Cream (Zinc) 78gm(tube)**Manufacturer:** Smith & Nephew **Manufacturer Code:** 59431279**Features:** Restores moisture to very dry skin, Contains petrolatum and Vitamin E, Conditions and soothes dry skin.**Usage:** Moisturizes and protects dry skin.**Precautions:** For external use only. If irritation occurs, consult a physician.**Size:** 78g**Standard Box Quantity:** 1


COL3705

Calcium Alginate Dressings

Biatain Alginate 5cm x 5cm**Manufacturer:** Coloplast **Manufacturer Code:** 3705**Features:** Biatain Alginate is a highly absorbent alginate dressing for moderately to heavily exuding wounds of any size, shape and depth. The design and structural integrity of the dressing allow you to cut the dressing to fit to any size, shape or depth of a wound. The dressing will conform to the wound bed which is important for the wound healing process.**Usage:** Biatain Alginate can be used on all moderate to heavily exuding partial to full thickness wounds, including leg ulcers, pressure ulcers, diabetic foot ulcers, donor sites.**Precautions:** Biatain Alginate should not be used on dry or lightly exuding wounds or to control heavy bleeding.**Size:** 5x5cm**Standard Box Quantity:** Each

COL3710

Calcium Alginate Dressings

Biatain Alginate 10cm x 10cm**Manufacturer:** Coloplast **Manufacturer Code:** 3710**Features:** Biatain Alginate is a highly absorbent alginate dressing for moderately to heavily exuding wounds of any size, shape and depth. The design and structural integrity of the dressing allow you to cut the dressing to fit to any size, shape or depth of a wound. The dressing will conform to the wound bed which is important for the wound healing process.**Usage:** Biatain Alginate can be used on all moderate to heavily exuding partial to full thickness wounds, including leg ulcers, pressure ulcers, diabetic foot ulcers, donor sites.**Precautions:** Biatain Alginate should not be used on dry or lightly exuding wounds or to control heavy bleeding.**Size:** 10x10cm**Standard Box Quantity:** Each

Biatain Alginate Rope 16"**Manufacturer:** Coloplast **Manufacturer Code:** 3740**Features:** Biatain Alginate dressing and filler are highly absorbent alginate dressings with vertical absorption to reduce maceration. On contact with wound exudate, Biatain Alginate forms a translucent, cohesive, soft gel for managing moderate to heavily exuding wounds.**Usage:** Can be used for any wound shape and cavity wounds, as the gel will conform to the shape of the wound. Should be used with a secondary cover dressing.**Precautions:****Size:** 2.5cm X 44cm**Standard Box Quantity:** Each**Kaltostat 2 gm Rope****Manufacturer:** ConvaTec **Manufacturer Code:** 168117**Features:** KALTOSTAT® dressing is an absorbent, calcium sodium alginate (80% Ca and 20% Na) with a high Guluronic acid content. When in contact with wound exudate, calcium ions in the dressing exchange with sodium ions in the exudate causing the dressing to transform from a dry, fibrous state to a firm, moist gel. The dressing provides a moist environment helping to create optimal wound healing conditions and maintains its integrity allowing for easy, atraumatic removal.**Usage:** KALTOSTAT® dressing is indicated as a primary dressing for moderate to highly exuding wounds and for wounds with minor bleeding. Chronic wounds: leg ulcers, pressure ulcers, diabetic ulcers and fungating lesions. Acute wounds: donor sites, abrasions, lacerations and post-surgical wounds.**Precautions:** This product does not contain latex.**Size:** 2 gm**Standard Box Quantity:** 1**Kaltostat 5 x 5cm****Manufacturer:** ConvaTec **Manufacturer Code:** 168210**Features:** KALTOSTAT® dressing is an absorbent, calcium sodium alginate (80% Ca and 20% Na) with a high Guluronic acid content. When in contact with wound exudate, calcium ions in the dressing exchange with sodium ions in the exudate causing the dressing to transform from a dry, fibrous state to a firm, moist gel. The dressing provides a moist environment helping to create optimal wound healing conditions and maintains its integrity allowing for easy, atraumatic removal.**Usage:** KALTOSTAT® dressing is indicated as a primary dressing for moderate to highly exuding wounds and for wounds with minor bleeding. Chronic wounds: leg ulcers, pressure ulcers, diabetic ulcers and fungating lesions. Acute wounds: donor sites, abrasions, lacerations and post-surgical wounds.**Precautions:** This product does not contain latex.**Size:** 5x5cm**Standard Box Quantity:** 1

CON168212

Calcium Alginate Dressings

Kaltostat 7.5 x 12cm

Manufacturer: ConvaTec **Manufacturer Code:** 168212

Features: KALTOSTAT® dressing is an absorbent, calcium sodium alginate (80% Ca and 20% Na) with a high Guluronic acid content. When in contact with wound exudate, calcium ions in the dressing exchange with sodium ions in the exudate causing the dressing to transform from a dry, fibrous state to a firm, moist gel. The dressing provides a moist environment helping to create optimal wound healing conditions and maintains its integrity allowing for easy, atraumatic removal.

Usage: KALTOSTAT® dressing is indicated as a primary dressing for moderate to highly exuding wounds and for wounds with minor bleeding. Chronic wounds: leg ulcers, pressure ulcers, diabetic ulcers and fungating lesions. Acute wounds: donor sites, abrasions, lacerations and post-surgical wounds.

Precautions: This product does not contain latex.

Size: 7.5x12cm

Standard Box Quantity: 1



CON168214

Calcium Alginate Dressings

Kaltostat 10 x 19cm

Manufacturer: ConvaTec **Manufacturer Code:** 168214

Features: KALTOSTAT® dressing is an absorbent, calcium sodium alginate (80% Ca and 20% Na) with a high Guluronic acid content. When in contact with wound exudate, calcium ions in the dressing exchange with sodium ions in the exudate causing the dressing to transform from a dry, fibrous state to a firm, moist gel. The dressing provides a moist environment helping to create optimal wound healing conditions and maintains its integrity allowing for easy, atraumatic removal.

Usage: KALTOSTAT® dressing is indicated as a primary dressing for moderate to highly exuding wounds and for wounds with minor bleeding. Chronic wounds: leg ulcers, pressure ulcers, diabetic ulcers and fungating lesions. Acute wounds: donor sites, abrasions, lacerations and post-surgical wounds.

Precautions: This product does not contain latex.

Size: 10x19cm

Standard Box Quantity: 1



MOL285280

Calcium Alginate Dressings

Mesalt Ribbon 2 x 100cm

Manufacturer: Molnlycke **Manufacturer Code:** 285280

Features: Mesalt is made of an absorbent, viscose/polyester nonwoven dressing impregnated with sodium chloride.

Usage: The exudate from the wound releases the sodium chloride from the dressing. Mesalt stimulates the cleansing of heavily draining and discharging wounds in the inflammatory phase by absorbing exudate. Change rate is B.I.D to daily if necessary.

Precautions: Mesalt is not to be used on low exuding or dry wounds.

Size: 2cm x 100cm

Standard Box Quantity: 1



Mesalt 10 x 10 cm**Manufacturer:** Molnlyke **Manufacturer Code:** 286080**Features:** Mesalt is made of an absorbent, viscose/polyester nonwoven dressing impregnated with sodium chloride.**Usage:** The exudate from the wound releases the sodium chloride from the dressing. Mesalt stimulates the cleansing of heavily draining and discharging wounds in the inflammatory phase by absorbing exudate. Change rate is B.I.D to daily if necessary.**Precautions:** Mesalt is not to be used on low exuding or dry wounds.**Size:** 10cm x 10cm**Standard Box Quantity:** 1**Algisite M 5 x 5cm****Manufacturer:** Smith & Nephew **Manufacturer Code:** 66000519**Features:** Algisite M is a calcium-alginate dressing which forms a soft, integral gel when it comes into contact with wound exudate. Algisite M can absorb 13 times its own weight in exudate.**Usage:** Algisite M may be used for the management of minor conditions such as Lacerations, Abrasions, Skin tears and Minor burns. Under the care of a healthcare professional Algisite M may be used for the management of Leg Ulcers, Pressure ulcers, Diabetic foot ulcers and Surgical wounds.**Precautions:** Ensure the packing in cavities, sinuses etc. has been removed. The dressing may adhere if used on very lightly exuding wounds. If the dressing is not easily removed, moisten it with saline, then remove. Do not use on patients with a known allergy to alginates.**Size:** 5x5cm**Standard Box Quantity:** 1**Algisite M 2 x 30cm****Manufacturer:** Smith & Nephew **Manufacturer Code:** 66000522**Features:** Algisite M is a calcium-alginate dressing which forms a soft, integral gel when it comes into contact with wound exudate. Algisite M can absorb 13 times its own weight in exudate.**Usage:** Algisite M may be used for the management of minor conditions such as Lacerations, Abrasions, Skin tears and Minor burns. Under the care of a healthcare professional Algisite M may be used for the management of Leg Ulcers, Pressure ulcers, Diabetic foot ulcers and Surgical wounds.**Precautions:** Ensure the packing in cavities, sinuses etc. has been removed. The dressing may adhere if used on very lightly exuding wounds. If the dressing is not easily removed, moisten it with saline, then remove. Do not use on patients with a known allergy to alginates.**Size:** 2x30cm**Standard Box Quantity:** 1

BARVUPD1012

Catheters

Securement Statlock for Catheter 10-12FR

Manufacturer: Bard **Manufacturer Code:** VUPD1012

Features: Advanced technologies now combine with the StatLock® Stabilization Device design sophistication to deliver skin-safe, suture-free, releasable stabilization of even the most lubricious catheters. The StatLock® Universal Plus Stabilization Device comes in four sizes, to fit virtually all percutaneous drainage catheters.

Usage: For suture-free stabilization of biliary and nephrostomy percutaneous drainage, and peritoneal lavage catheters.

Precautions: More effective than traditional forms of tape and suture stabilization to achieve positive outcomes.

Size: 10-12FR

Standard Box Quantity: 1



BARVUPD1214

Catheters

Securement Statlock for Catheter 12-14FR

Manufacturer: Bard **Manufacturer Code:** VUPD1214

Features: Advanced technologies now combine with the StatLock® Stabilization Device design sophistication to deliver skin-safe, suture-free, releasable stabilization of even the most lubricious catheters. The StatLock® Universal Plus Stabilization Device comes in four sizes, to fit virtually all percutaneous drainage catheters.

Usage: For suture-free stabilization of biliary and nephrostomy percutaneous drainage, and peritoneal lavage catheters.

Precautions: More effective than traditional forms of tape and suture stabilization to achieve positive outcomes.

Size: 12-14FR

Standard Box Quantity: 1



BARVUPD1416

Catheters

Securement Statlock for Catheter 14-16FR

Manufacturer: Bard **Manufacturer Code:** VUPD1416

Features: Advanced technologies now combine with the StatLock® Stabilization Device design sophistication to deliver skin-safe, suture-free, releasable stabilization of even the most lubricious catheters. The StatLock® Universal Plus Stabilization Device comes in four sizes, to fit virtually all percutaneous drainage catheters.

Usage: For suture-free stabilization of biliary and nephrostomy percutaneous drainage, and peritoneal lavage catheters.

Precautions: More effective than traditional forms of tape and suture stabilization to achieve positive outcomes.

Size: 14-16FR

Standard Box Quantity: 1



Securement Statlock for Catheter 6-8FR**Manufacturer:** Bard **Manufacturer Code:** VUPD68

Features: Advanced technologies now combine with the StatLock® Stabilization Device design sophistication to deliver skin-safe, suture-free, releasable stabilization of even the most lubricious catheters. The StatLock® Universal Plus Stabilization Device comes in four sizes, to fit virtually all percutaneous drainage catheters.

Usage: For suture-free stabilization of biliary and nephrostomy percutaneous drainage, and peritoneal lavage catheters.

Precautions: More effective than traditional forms of tape and suture stabilization to achieve positive outcomes.

Size: 6-8FR**Standard Box Quantity:** 1**Flexitrack Anchoring Device****Manufacturer:** ConvaTec **Manufacturer Code:** 37449

Features: FLEXI-TRAK® Anchoring Device is a skin-friendly, modern acrylic adhesive that securely holds tubes, drains and catheters in place on the skin.

Usage: ConvaTec FLEXI TRAK Large Anchoring Device For anchoring tubes and catheters to skin, such as urinary catheters, gastrostomy tubes, chest tubes. Adhesive.

Precautions: This product does not contain latex.

Size: 38mm X 100mm**Standard Box Quantity:** 1**Suction Catheter Coude 12FR w Port****Manufacturer:** Medline **Manufacturer Code:** DYND41801

Features: Valve to help prevent contact with fluids and helping to prevent the release of microbial contamination during suctioning process.

Usage: Suction catheters used to remove secretions from the mouth, oropharynx, trachea and bronchi.

Precautions: Sterile technique required.

Size: 12FR**Standard Box Quantity:** 1

SL501014

Catheters

Suction Straight 14FR (Disposable)

Manufacturer: Med-Rx **Manufacturer Code:** 50-1014

Features: Valve to help prevent contact with fluids and helping to prevent the release of microbial contamination during suctioning process.

Usage: Suction catheters used to remove secretions from the mouth, oropharynx, trachea and bronchi.

Precautions: Sterile technique required.

Size: 14FR

Standard Box Quantity: 1



URO7000

Catheters

Catheter Clamp

Manufacturer: Urocare **Manufacturer Code:** 7000-00

Features: Urocare Adjustable Thumb Clamp. Two position.

Usage: May be used with clear-vinyl, latex or white rubber drainage or extension tubing as a flow control device.

Precautions:

Size: Universal

Standard Box Quantity: 1



SN7064

Charcoal Dressings

Carbonet 10 x 10cm

Manufacturer: Smith & Nephew **Manufacturer Code:** 7064

Features: Carbonet dressings have wound contact layer and absorbent padding layer. Absorbs bacteria and offensive odours. They utilise the important, well-proven benefits of an activated charcoal layer, with the additional features of absorbency and low-adherence.

Usage: Discharging infected and malodorous wounds such as Venous leg ulcers, Fungating carcinomas (eg breast), Faecal fistulae and Infected surgical wounds.

Precautions: Latex Free.

Size: 10x10cm

Standard Box Quantity: 1



Carbonet 10 x 20cm**Manufacturer:** Smith & Nephew **Manufacturer Code:** 7065**Features:** Carbonet dressings have wound contact layer and absorbent padding layer. Absorbs bacteria and offensive odours. They utilise the important, well-proven benefits of an activated charcoal layer, with the additional features of absorbency and low-adherence.**Usage:** Discharging infected and malodorous wounds such as Venous leg ulcers, Fungating carcinomas (eg breast), Faecal fistulae and Infected surgical wounds.**Precautions:** Latex Free.**Size:** 10x20cm**Standard Box Quantity:** 1**6 inch Coban 2 Layer Compression Kit****Manufacturer:** 3M **Manufacturer Code:** 20096**Features:** 3M™ Coban™ 2 Layer Compression System is designed to consistently deliver comfortable, therapeutic compression for venous insufficiency patients with an ABPI of 0.8 – 1.0 mmHg. This patented two-layer compression system is clinically proven to: Be safe, comfortable and effective, Provide sustained compression for up to 7 days, easy to apply, enable patients to wear their ordinary footwear and clothing.**Usage:** For the treatment of venous leg ulcers, lymphedema and other conditions where compression therapy is appropriate. The secure sleeve provides the appropriate amount of sustained compression to optimize patients' muscle movement to move fluid improving venous return and reducing edema.**Precautions:** Make sure you do not wrap the limb too tightly. If a limb is swelling, wrapping too tightly could cut off circulation.**Size:** 6 inches**Standard Box Quantity:** 1**Coban 2 Lite 1"****Manufacturer:** 3M **Manufacturer Code:** 20721**Features:** Coban 2 - Lite Compression Layer only, 1 in x 3.8 yd (2.5 cm x 3.5 m)**Usage:** Used over 3M™ Coban™ 2 Lite Compression System - comfort foam layer (ordered separately)**Precautions:** Does not contain latex.**Size:** 1 inch x 3.5M**Standard Box Quantity:** Each

3M2094

Compression Bandages

4 inch Coban 2 Layer Compression Kit

Manufacturer: 3M **Manufacturer Code:** 2094

Features: 3M™ Coban™ 2 Layer Compression System is designed to consistently deliver comfortable, therapeutic compression for venous insufficiency patients with an ABPI of 0.8 – 1.0 mmHg. This patented two-layer compression system is clinically proven to: Be safe, comfortable and effective, Provide sustained compression for up to 7 days, easy to apply, enable patients to wear their ordinary footwear and clothing.

Usage: For the treatment of venous leg ulcers, lymphedema and other conditions where compression therapy is appropriate. The secure sleeve provides the appropriate amount of sustained compression to optimize patients' muscle movement to move fluid improving venous return and reducing edema.

Precautions: Make sure you do not wrap the limb too tightly. If a limb is swelling, wrapping too tightly could cut off circulation.

Size: 4 inches

Standard Box Quantity: 1



3M2794

Compression Bandages

Coban 2 Layer Lite Compression Kit

Manufacturer: 3M **Manufacturer Code:** 2794-HC

Features: 3M™ Coban™ 2 Layer Lite Compression System reduces the risk of tissue damage on mixed venous and arterial disease patients with an ABPI of 0.5 – 0.8 mmHg. This patented two-layer compression system is clinically proven to be safe, comfortable and effective, provide sustained compression for up to 7 days, be preferred by patients for comfort, be easy to apply, enable patients to wear their ordinary footwear and clothing.

Usage: For the treatment of venous leg ulcers, lymphedema and other conditions where compression therapy is appropriate. The secure sleeve provides the appropriate amount of sustained compression to optimize patients' muscle movement to move fluid, improving venous return and reducing edema.

Precautions: Make sure you do not wrap the limb too tightly. If a limb is swelling, wrapping too tightly could cut off circulation.

Size: 4 inches

Standard Box Quantity: 1



BSN0102600

Compression Bandages

Comprilan Stretch Bandage 6cm x 5m

Manufacturer: BSN **Manufacturer Code:** 102600

Features: Provides clinically relevant levels of compression including graduated compression when applied with a spiral technique. Two bandages applied one over another will give a high level of compression useable for the treatment of venous leg ulcers. 100% cotton for easy application and is well tolerated. Can be washed and reapplied to the same patient.

Usage: Indicated to assist in the management of venous and mixed aetiology leg ulcers. Effective in helping to control oedma associated with venous disorders.

Precautions: Should not be used on patients with ABPI of less than 0.5 or with advanced small vessel disease. Failure to detect significantly reduced arterial flow can result in pressure necrosis, amputation or even death. Apply as early in the day as possible and remove at night.

Size: 6cm x 5m

Standard Box Quantity: 1



Comprilan Stretch Bandage 8cm x 5m**Manufacturer:** BSN **Manufacturer Code:** 102700

Features: Provides clinically relevant levels of compression including graduated compression when applied with a spiral technique. Two bandages applied one over another will give a high level of compression useable for the treatment of venous leg ulcers. 100% cotton for easy application and is well tolerated. Can be washed and reapplied to the same patient.

Usage: Indicated to assist in the management of venous and mixed aetiology leg ulcers. Effective in helping to control oedma associated with venous disorders.

Precautions: Should not be used on patients with ABPI of less than 0.5 or with advanced small vessel disease. Failure to detect significantly reduced arterial flow can result in pressure necrosis, amputation or even death. Apply as early in the day as possible and remove at night.

Size: 8cm x 5m**Standard Box Quantity:** 1**Comprilan Stretch Bandage 10cm x 5m****Manufacturer:** BSN **Manufacturer Code:** 102800

Features: Provides clinically relevant levels of compression including graduated compression when applied with a spiral technique. Two bandages applied one over another will give a high level of compression useable for the treatment of venous leg ulcers. 100% cotton for easy application and is well tolerated. Can be washed and reapplied to the same patient.

Usage: Indicated to assist in the management of venous and mixed aetiology leg ulcers. Effective in helping to control oedma associated with venous disorders.

Precautions: Should not be used on patients with ABPI of less than 0.5 or with advanced small vessel disease. Failure to detect significantly reduced arterial flow can result in pressure necrosis, amputation or even death. Apply as early in the day as possible and remove at night.

Size: 10cm x 5m**Standard Box Quantity:** 1**Comprilan 12cm x 5m****Manufacturer:** BSN **Manufacturer Code:** 1029

Features: Provides clinically relevant levels of compression including graduated compression when applied with a spiral technique. Two bandages applied one over another will give a high level of compression useable for the treatment of venous leg ulcers. 100% cotton for easy application and is well tolerated. Can be washed and reapplied to the same patient.

Usage: Indicated to assist in the management of venous and mixed aetiology leg ulcers. Effective in helping to control oedma associated with venous disorders.

Precautions: Should not be used on patients with ABPI of less than 0.5 or with advanced small vessel disease. Failure to detect significantly reduced arterial flow can result in pressure necrosis, amputation or even death. Apply as early in the day as possible and remove at night.

Size: 12cm x 5m**Standard Box Quantity:** Each

CON650948

Compression Bandages

Surepress Absorbent Padding 10cm x 3m

Manufacturer: Convatec **Manufacturer Code:** 650948

Features: A soft non-woven and absorbent bandage composed of a blend of viscose / bi-component and super-absorbent fibres. This bandage protects the skin, provides additional absorption and even distribution of pressure under the SurePress® High Compression Bandage.

Usage: For use under SurePress® High Compression Bandages. The SurePress® Absorbent Padding helps to distribute graduated compression and protect skin from exudates leaking from the primary dressing.

Precautions: Latex Free.

Size: 10 cm x 3M

Standard Box Quantity: Each



MED618

Compression Bandages

Tensor With Clips 10cm

Manufacturer: Medicom **Manufacturer Code:** 618

Features: An elasticized bandage used to provide support to a sprained or strained limb. This type of bandage may also be known as an ace bandage or ace wrap.

Usage: This bandage design consists of an elasticized knit. The bandage can be wrapped around an ankle, wrist, or knee and it will conform, holding tension and providing support. Tensor bandages are held in place with small, toothed metal clips.

Precautions: Make sure you do not wrap the limb too tightly. If a limb is swelling, wrapping too tightly could cut off circulation.

Size: 10cm

Standard Box Quantity: 1



SM72414

Compression Bandages

Dufore 4-Layer Compression Bandage Kit

Manufacturer: Derma Science **Manufacturer Code:** DUP72414

Features: This system is designed to sustain effective graduated compression for up to 7 days. The system is comprised of four layers. The first layer, a padding absorption layer, provides padding for bony prominences, distributes compression evenly over the leg and absorbs exudate. The second layer provides short-stretch compression by resisting expansion of the calf muscle during limb movement. It works with the third and fourth layer elastic compression bandages to effectively minimize edema and increase venous return. A wound contact dressing is included with the system to help maintain a moist wound healing environment and to protect the wound from lint or outside contaminants.

Usage: Sustained compression has been recognized as the most important factor in managing venous leg ulcers.

Precautions: Latex Free.

Size: 4 inches

Standard Box Quantity: 1



SN4948

Compression Bandages

Viscopaste Bandage 7.5cm

Manufacturer: Smith & Nephew **Manufacturer Code:** 4948

Features: Viscopaste* is a general purpose first-line treatment for wet and dry leg ulcers. Viscopaste bandages soothe and protect the skin and are easily removed. Impregnated with zinc paste.

Usage: The softness of the Viscopaste bandages means they are simple to apply and easy to remove (soft and pliable after application). Bandages will not dry out during wear time. Bandages provide a moist wound healing environment over the ulcer.

Precautions: Where venous insufficiency exists, Viscopaste should be used under graduated compression bandaging, after first assessing the patient to exclude arterial disease. The use of Doppler ultrasound is recommended for this purpose. Compression bandaging should not be used on patients who have arterial disease (ABPI<0.8) and is not recommended for use on diabetic patients with advanced small vessel disease.

Size: 7.5cm

Standard Box Quantity: 1



SN66000415

Compression Bandages

Profore Lite Compression Kit

Manufacturer: Smith & Nephew **Manufacturer Code:** 66000415

Features: A 3-Layers sustained compression bandage system. Each system pack contains the following: One Profore wound contact layer. One Profore #1 natural padding bandage One Profore #2 light conformable dressing. One Profore #4 flexible cohesive bandage."

Usage: PROFORE is indicated for the management of venous leg ulcers and associated conditions. PROFORE Lite is specifically designed for conditions in which there is both venous and arterial involvement.

Precautions: The standard kit contains latex, however there are latex free formulation kits available for all ankle circumferences and compression variants which are specially formulated with no latex containing ingredients. Do not use on patients with an ankle: brachial pressure index (ABPI) of less than 0.6 (Profore Lite) or on diabetic patients with advanced small vessel disease.

Size:

Standard Box Quantity: 1



3M1583

Conforming Bandages

Coban 3 inch Self Adherent Roll

Manufacturer: 3M **Manufacturer Code:** 1583K

Features: Coban consists of a water-vapour permeable, non-woven polyester fabric containing longitudinal strands of elastane. The bandage is able to stick to itself but not to skin or clothing.

Usage: Coban is applied using a standard bandaging technique. Working from the inner aspect of the leg, the bandage is applied in the form of a spiral ensuring a 50% overlap between turns. Coban does not require pins or tapes to hold it in place. Coban may be left undisturbed for extended periods, (up to two weeks has been reported).

Precautions: Care should be taken to ensure that it is not applied with excessive tension. If the bandage is too tight or causes the patient any discomfort, it should be removed and reapplied.

Size: 7.5cm X 4.6m

Standard Box Quantity: 1



3M1584

Conforming Bandages

Coban 4 inch Self Adherent Roll

Manufacturer: 3M **Manufacturer Code:** 1584K

Features: Coban consists of a water-vapour permeable, non-woven polyester fabric containing longitudinal strands of elastane. The bandage is able to stick to itself but not to skin or clothing.

Usage: Coban is applied using a standard bandaging technique. Working from the inner aspect of the leg, the bandage is applied in the form of a spiral ensuring a 50% overlap between turns. Coban does not require pins or tapes to hold it in place. Coban may be left undisturbed for extended periods, (up to two weeks has been reported).

Precautions: Care should be taken to ensure that it is not applied with excessive tension. If the bandage is too tight or causes the patient any discomfort, it should be removed and reapplied.

Size: 10cm X 4.6m

Standard Box Quantity: 1



MED205

Conforming Bandages

Conform Bandage 5cm (roll)

Manufacturer: MDA Group **Manufacturer Code:** 205

Features: A conforming bandage is a bandage designed to mold itself to the shape of the area where it is applied. This helps the bandage stay in place, providing extra security.

Usage: For securing wound dressings.

Precautions: Should the patient develop pain or pale, cool or numb extremities distal to the dressing, the bandages should be promptly removed.

Size: 5cm

Standard Box Quantity: 1



MED305

Conforming Bandages

Conform Bandage 15cm (roll)

Manufacturer: MDA Group **Manufacturer Code:** MED305

Features: A conforming bandage is a bandage designed to mold itself to the shape of the area where it is applied. This helps the bandage stay in place, providing extra security.

Usage: For securing wound dressings.

Precautions: Should the patient develop pain or pale, cool or numb extremities distal to the dressing, the bandages should be promptly removed.

Size: 7.5cm

Standard Box Quantity: Each



MED405

Conforming Bandages

Conform Bandage 10cm (roll)

Manufacturer: MDA Group **Manufacturer Code:** 405

Features: A conforming bandage is a bandage designed to mold itself to the shape of the area where it is applied. This helps the bandage stay in place, providing extra security.

Usage: For securing wound dressings.

Precautions: Should the patient develop pain or pale, cool or numb extremities distal to the dressing, the bandages should be promptly removed.

Size: 10cm

Standard Box Quantity: 1



MED605

Conforming Bandages

Conform Bandage 15cm (roll)

Manufacturer: MDA Group **Manufacturer Code:** 605

Features: A conforming bandage is a bandage designed to mold itself to the shape of the area where it is applied. This helps the bandage stay in place, providing extra security.

Usage: For securing wound dressings.

Precautions: Should the patient develop pain or pale, cool or numb extremities distal to the dressing, the bandages should be promptly removed.

Size: 15cm

Standard Box Quantity: 1



DYN4521

Dressings - Suture / Staple Removers

Suture Removal Kit (VAC Only, 1 per client)

Manufacturer: Medline **Manufacturer Code:** 4521

Features: Components are packaged in sequence of use, allowing the clinician to progress efficiently through procedures. Suture Removal Tray With Metal Littauer Scissors Contains: Scissors (Wire Littauer), Forceps (Metal), Gauze (3" x 3").

Usage: Used to remove sutures.

Precautions: Sterile

Size: Kit

Standard Box Quantity: Each



SM81121

Dressings - Tapes

Steri-Strips, 6/Pkg, 1/2"

Manufacturer: Meditron **Manufacturer Code:** DKC81121

Features: Wound closure strips are made of 100% polyamide material coated with a hypoallergenic adhesive which allows almost pain free closure of the wound.

Usage: Does not traumatise the edges of the wound allowing almost pain free wound closure. Reliably and securely adapts to the edges of the wound.

Precautions: Do not use if individual package is not intact. Single use. **STORAGE** Keep away from light in a dry place below 25 ° C.

Size: 12.7mm X 102mm

Standard Box Quantity: Each



COV8DIC

Equipment - Respiratory Aids

Inner Cannulas, Disposable

Manufacturer: Covidien **Manufacturer Code:** 8DIC

Features: Shiley™ inner cannula disposable (DIC) features integral 15mm snap-lock connector for quick, safe disposal after single use

Usage:

Precautions:

Size: Size #8

Standard Box Quantity: Each



GLE1075KIT

Equipment - Respiratory Aids

Aerosol Kit, Trach (Nebulizer Cup, Mask & 7ft tubing)

Manufacturer: Glenwood **Manufacturer Code:** GK1808

Features: Includes tracheostomy mask and tubing.

Usage: For the administration of medication in aerosol form for tracheostomy patients.

Precautions:

Size: 7 feet tubing

Standard Box Quantity: Each



Aerosol Kit, Adult (Nebulizer Cup, Mask & 7ft tubing)**Manufacturer:** Glenwood **Manufacturer Code:** GK1705**Features:** Includes adult aerosol mask and tubing.**Usage:** For the administration of medication in aerosol form.**Precautions:****Size:** Adult, 7 feet tubing**Standard Box Quantity:** Each**Suction Tip, Yankauer****Manufacturer:** Med-Rx **Manufacturer Code:** 30-3064**Features:** Solid one-piece molded construction, shatter proof transparent plastic, with a universal connector, smooth atraumatic bulb with smooth eyes, and easy-to-grasp ribbed handle at proximal end.**Usage:** For suctioning procedures.**Precautions:****Size:** 10in**Standard Box Quantity:** Each**Suction, Catheters, 10 French****Manufacturer:** Med-Rx **Manufacturer Code:** 50-1010**Features:** Valve to help prevent contact with fluids and helping to prevent the release of microbial contamination during suctioning process.**Usage:** Suction catheters used to remove secretions from the mouth, oropharynx, trachea and bronchii.**Precautions:** Sterile technique required.**Size:** 10FR**Standard Box Quantity:** Each

SL501012

Equipment - Respiratory Aids

Suction, Catheters, 12 French

Manufacturer: Med-Rx **Manufacturer Code:** 50-1012

Features: Valve to help prevent contact with fluids and helping to prevent the release of microbial contamination during suctioning process.

Usage: Suction catheters used to remove secretions from the mouth, oropharynx, trachea and bronchii.

Precautions: Sterile technique required.

Size: 12FR

Standard Box Quantity: Each



SL501016

Equipment - Respiratory Aids

Suction, Catheters, 16 French

Manufacturer: Med-Rx **Manufacturer Code:** 50-1016

Features: Valve to help prevent contact with fluids and helping to prevent the release of microbial contamination during suctioning process.

Usage: Suction catheters used to remove secretions from the mouth, oropharynx, trachea and bronchii.

Precautions: Sterile technique required.

Size: 16FR

Standard Box Quantity: Each



SL522014

Equipment - Respiratory Aids

Suction, Catheter, Coude Tip w/valve, 14 F

Manufacturer: Med-Rx **Manufacturer Code:** 52-2014

Features: Valve to help prevent contact with fluids and helping to prevent the release of microbial contamination during suctioning process.

Usage: Suction catheters used to remove secretions from the mouth, oropharynx, trachea and bronchii.

Precautions: Sterile technique required.

Size: 14FR

Standard Box Quantity: Each



Suction, 6' Connecting Tubing**Manufacturer:** Med-Rx **Manufacturer Code:** 70-7072**Features:** extension/connection tubing for suction tubing, 7/32in x 72in**Usage:****Precautions:** Use aseptic technique. Follow physicians instructions in opening package and use of contents.**Size:** 7/32 X 72 inches**Standard Box Quantity:** Each

COL33434

Foam Adhesive Dressings

Biatain Silicone Foam Dressing w Border 7.5 x 7.5cm**Manufacturer:** Coloplast **Manufacturer Code:** 33434**Features:** Soft perforated silicone adhesive wound contact layer ensures a secure fit to keep the dressing in place. The soft and flexible materials allow Biatain Silicone to follow body movements enhancing the patient's ability to go about with their daily lives. The foam pad at skin level makes the dressing comfortable to wear - even under compression. The silicone adhesive wound contact layer ensures minimal pain upon dressing removal.**Usage:** Biatain Silicone is the versatile solution for moist wounds such as leg ulcers, non-infected diabetic foot ulcers, second degree burns, donor sites, post operative wounds, skin abrasions or pressure ulcers.**Precautions:** The dressing should not be left in place for more than 7 days without clinical inspection by a physician.**Size:** 7.5x7.5cm**Standard Box Quantity:** 1

COL33435

Foam Adhesive Dressings

Biatain Adhesive With Border 10cm X 10cm**Manufacturer:** Coloplast **Manufacturer Code:** 33435**Features:** Soft perforated silicone adhesive wound contact layer ensures a secure fit to keep the dressing in place. The soft and flexible material follows body movements enhancing the patient's ability to go about with their daily lives. The foam pad at skin level makes the dressing comfortable to wear even under compression. The silicone adhesive wound contact layer ensures atraumatic removal.**Usage:** For moist wounds such as leg ulcers, non-infected diabetic foot ulcers, second degree burns, donor sites, post operative wounds, skin abrasions or pressure ulcers.**Precautions:** The dressing should not be left in place for more than 7 days without clinical inspection by a physician.**Size:** 10 x 10cm**Standard Box Quantity:** Each

COL33436

Foam Adhesive Dressings

Biatain Silicone Foam Dressing w Border 12.5 x 12.5cm**Manufacturer:** Coloplast **Manufacturer Code:** 33436

Features: Soft perforated silicone adhesive wound contact layer ensures a secure fit to keep the dressing in place. The soft and flexible materials allow Biatain Silicone to follow body movements enhancing the patient's ability to go about with their daily lives. The foam pad at skin level makes the dressing comfortable to wear - even under compression. The silicone adhesive wound contact layer ensures minimal pain upon dressing removal.

Usage: Biatain Silicone is the versatile solution for moist wounds such as leg ulcers, non-infected diabetic foot ulcers, second degree burns, donor sites, post operative wounds, skin abrasions or pressure ulcers.

Precautions: The dressing should not be left in place for more than 7 days without clinical inspection by a physician.

Size: 12.5x12.5cm**Standard Box Quantity:** 1

COL33438

Foam Adhesive Dressings

Biatain Silicone Foam Dressing w Border 17.5 x 17.5cm**Manufacturer:** Coloplast **Manufacturer Code:** 33438

Features: Soft perforated silicone adhesive wound contact layer ensures a secure fit to keep the dressing in place. The soft and flexible materials allow Biatain Silicone to follow body movements enhancing the patient's ability to go about with their daily lives. The foam pad at skin level makes the dressing comfortable to wear - even under compression. The silicone adhesive wound contact layer ensures minimal pain upon dressing removal.

Usage: Biatain Silicone is the versatile solution for moist wounds such as leg ulcers, non-infected diabetic foot ulcers, second degree burns, donor sites, post operative wounds, skin abrasions or pressure ulcers.

Precautions: The dressing should not be left in place for more than 7 days without clinical inspection by a physician.

Size: 17.5x17.5cm**Standard Box Quantity:** 1

COL33444

Foam Adhesive Dressings

Biatain® Silicone Lite Dressing 7.5 x 7.5cm**Manufacturer:** Coloplast **Manufacturer Code:** 33444

Features: Biatain® Silicone Lite is a thin and ultra flexible foam dressing with gentle silicone adhesive ensuring minimal pain upon removal.

Usage: Biatain Silicone is the versatile solution for moist wounds such as leg ulcers, non-infected diabetic foot ulcers, second degree burns, donor sites, post operative wounds, skin abrasions or pressure ulcers.

Precautions: The Biatain Silicone dressing can not be cut, as this would mean that the SAP particles will be released into the wound.

Size: 7.5 x 7.5cm**Standard Box Quantity:** 1

COL33445

Foam Adhesive Dressings

Biatain Silicone Lite Foam Dressing w Border 10 x 10cm**Manufacturer:** Coloplast **Manufacturer Code:** 33445**Features:** Biatain® Silicone Lite is a thin and ultra flexible foam dressing with gentle silicone adhesive ensuring minimal pain upon removal.**Usage:** Biatain Silicone is the versatile solution for moist wounds such as leg ulcers, non-infected diabetic foot ulcers, second degree burns, donor sites, post operative wounds, skin abrasions or pressure ulcers.**Precautions:** The Biatain Silicone dressing can not be cut, as this would mean that the SAP particles will be released into the wound.**Size:** 10x10cm**Standard Box Quantity:** 1

COL33446

Foam Adhesive Dressings

Biatain Silicone Lite Foam Dressing w Border 12.5 x 12.5cm**Manufacturer:** Coloplast **Manufacturer Code:** 33446**Features:** Biatain Foam dressings provide superior exudate management for faster wound healing. It also offers high patient comfort due to a soft and flexible foam with bevelled edges. The superior absorption and retention properties of Biatain, minimises risk of leakage and maceration. Biatain offers a cost effective treatment due to long wear time. The semi-permeable top film is shower proof and provides a barrier function to bacteria.**Usage:** Biatain Silicone is the versatile solution for moist wounds such as leg ulcers, non-infected diabetic foot ulcers, second degree burns, donor sites, post operative wounds, skin abrasions or pressure ulcers.**Precautions:** The Biatain Silicone dressing can not be cut, as this would mean that the SAP particles will be released into the wound.**Size:** 12.5x12.5cm**Standard Box Quantity:** 1

COL3483

Foam Adhesive Dressings

Biatain Adhesive Foam Small Sacral Dressing w Border 17 x 17cm**Manufacturer:** Coloplast **Manufacturer Code:** 3483**Features:** Biatain Foam dressings provide superior exudate management for faster wound healing. It also offers high patient comfort due to a soft and flexible foam with bevelled edges. The superior absorption and retention properties of Biatain, minimises risk of leakage and maceration. Biatain offers a cost effective treatment due to long wear time. The semi-permeable top film is shower proof and provides a barrier function to bacteria.**Usage:** Biatain Adhesive Foam Dressing / Shapes Is suitable for wounds where an anatomically shaped dressing is needed. including pressure ulcers, second degree burns, donor sites, postoperative wounds and skin abrasions. Biatain foam dressings may be used on patients who are in treatment for a local or systemic infection under the discretion of a physician. Suitable for use in combination with compression therapy.**Precautions:****Size:** 17x17cm**Standard Box Quantity:** 1

Biatain Adhesive Foam Large Sacral Dressing w Border 23 x 23cm**Manufacturer:** Coloplast **Manufacturer Code:** 3485

Features: Biatain Foam dressings provide superior exudate management for faster wound healing. It also offers high patient comfort due to a soft and flexible foam with bevelled edges. The superior absorption and retention properties of Biatain, minimises risk of leakage and maceration. Biatain offers a cost effective treatment due to long wear time. The semi-permeable top film is shower proof and provides a barrier function to bacteria.

Usage: Biatain Adhesive Foam Dressing / Shapes Is suitable for wounds where an anatomically shaped dressing is needed. including pressure ulcers, second degree burns, donor sites, postoperative wounds and skin abrasions. Biatain foam dressings may be used on patients who are in treatment for a local or systemic infection under the discretion of a physician. Suitable for use in combination with compression therapy.

Precautions:**Size:** 23x23cm**Standard Box Quantity:** 1**Aquacel Foam Adhesive w Border 12.5 x 12.5cm****Manufacturer:** ConvaTec **Manufacturer Code:** 420619

Features: Aquacel adhesive dressing consists of a soft absorbent foam pad with a water-proof top layer which protects the wound and helps to manage the moisture vapour from the exudate, and a Hydrofiber wound interface which creates a soft gel that contours to the wound bed, maintains a moist environment and aids in autolytic debridement. A gentle silicone adhesive border provides skin-friendly adhesion for non-traumatic dressing removal.

Usage: Primary or secondary dressing for both chronic and acute wounds, such as leg ulcers, pressure ulcers (Stages II-IV) and diabetic ulcers; surgical wounds; partial thickness burns; traumatic or surgical wounds left to heal by secondary intention such as dehiscent surgical incisions; dermatological and surgical incisions.

Precautions: The dressing can be cut to shape for convenience. If the dressing is cut, additional tape may be required to secure it.

Size: 12.5x12.5cm**Standard Box Quantity:** 1**Aquacel Foam Adhesive w Border 17.5 x 17.5cm****Manufacturer:** ConvaTec **Manufacturer Code:** 420621

Features: Aquacel adhesive dressing consists of a soft absorbent foam pad with a water-proof top layer which protects the wound and helps to manage the moisture vapour from the exudate, and a Hydrofiber wound interface which creates a soft gel that contours to the wound bed, maintains a moist environment and aids in autolytic debridement. A gentle silicone adhesive border provides skin-friendly adhesion for non-traumatic dressing removal.

Usage: Primary or secondary dressing for both chronic and acute wounds, such as leg ulcers, pressure ulcers (Stages II-IV) and diabetic ulcers; surgical wounds; partial thickness burns; traumatic or surgical wounds left to heal by secondary intention such as dehiscent surgical incisions; dermatological and surgical incisions.

Precautions: The dressing can be cut to shape for convenience. If the dressing is cut, additional tape may be required to secure it.

Size: 17.5x17.5cm**Standard Box Quantity:** 1

CON420623

Foam Adhesive Dressings

Aquacel Foam Adhesive w Border 21 x 21cm

Manufacturer: ConvaTec **Manufacturer Code:** 420623

Features: Aquacel adhesive dressing consists of a soft absorbent foam pad with a water-proof top layer which protects the wound and helps to manage the moisture vapour from the exudate, and a Hydrofiber wound interface which creates a soft gel that contours to the wound bed, maintains a moist environment and aids in autolytic debridement. A gentle silicone adhesive border provides skin-friendly adhesion for non-traumatic dressing removal.

Usage: Primary or secondary dressing for both chronic and acute wounds, such as leg ulcers, pressure ulcers (Stages II-IV) and diabetic ulcers; surgical wounds; partial thickness burns; traumatic or surgical wounds left to heal by secondary intention such as dehiscent surgical incisions; dermatological and surgical incisions.

Precautions: The dressing can be cut to shape for convenience. If the dressing is cut, additional tape may be required to secure it.

Size: 21x21cm

Standard Box Quantity: 1



CON420680

Foam Adhesive Dressings

Aquacel Foam Adhesive w Border 10 x 10cm

Manufacturer: ConvaTec **Manufacturer Code:** 420680

Features: AQUACEL® Foam dressings consists of a soft absorbent foam pad with a water-proof top layer which protects the wound and helps to manage the moisture vapour from the exudate, and a Hydrofiber wound interface which creates a soft gel that contours to the wound bed, maintains a moist environment and aids in autolytic debridement without damaging newly formed tissue.

Usage: Primary or secondary dressing for both chronic and acute wounds, such as leg ulcers, pressure ulcers (Stages II-IV) and diabetic ulcers; surgical wounds; partial thickness burns; traumatic or surgical wounds left to heal by secondary intention such as dehiscent surgical incisions; dermatological and surgical incisions.

Precautions: The dressing can be cut to shape for convenience. If the dressing is cut, additional tape may be required to secure it.

Size: 10x10cm

Standard Box Quantity: 1



CON420804

Foam Adhesive Dressings

Aquacel Foam Adhesive w Border 8 x 8cm

Manufacturer: ConvaTec **Manufacturer Code:** 420804

Features: Aquacel adhesive dressing consists of a soft absorbent foam pad with a water-proof top layer which protects the wound and helps to manage the moisture vapour from the exudate, and a Hydrofiber wound interface which creates a soft gel that contours to the wound bed, maintains a moist environment and aids in autolytic debridement. A gentle silicone adhesive border provides skin-friendly adhesion for non-traumatic dressing removal.

Usage: Primary or secondary dressing for both chronic and acute wounds, such as leg ulcers, pressure ulcers (Stages II-IV) and diabetic ulcers; surgical wounds; partial thickness burns; traumatic or surgical wounds left to heal by secondary intention such as dehiscent surgical incisions; dermatological and surgical incisions.

Precautions: The dressing can be cut to shape for convenience. If the dressing is cut, additional tape may be required to secure it.

Size: 8x8cm

Standard Box Quantity: 1



Mepilex Transfer 15 x 20cm**Manufacturer:** Molnlycke **Manufacturer Code:** 294800

Features: Mepilex Transfer is a soft silicone exudate transfer dressing. It consists of a flexible and hydrophilic polyurethane foam coated with a soft silicone layer (Safetac Technology). Mepilex Transfer moves exudate away from the wound, minimizes the risk of maceration and enables management of difficult-to-dress wounds.

Usage: Mepilex Transfer is designed for a wide range of exuding and difficult-to-dress wounds. Mepilex Transfer can also be used as a protective layer on non-exuding wounds and/or large areas of fragile skin. Mepilex Transfer can be used under compression.

Precautions: Mepilex Transfer may be left in place for several days depending on the condition of the wound and the surrounding skin, or as indicated by accepted clinical practice. A change in dressing regimen can result in an initial change in the level of wound exudate. On infected wounds the dressing should be changed at least daily until infection is controlled.

Size: 15x20cm**Standard Box Quantity:** 1**Biatain 10cm x 10cm****Manufacturer:** Coloplast **Manufacturer Code:** 3410

Features: The unique Biatain foam is softer and more flexible due to an optimised method of processing. This helps the dressings conform to individual body contours, and thus generates a more pleasant experience for the patient when wearing the dressings.

Usage: Select a Biatain pad where the foam pad overlaps the wound by 1-2 cm and apply. A healthcare professional should constantly monitor infected wounds and the healing environment.

Precautions: Do not use Biatain in combination with oxidising solutions. Ensure that all evaporating solution is completely dried off prior to application of dressing.

Size: 10cm x 10cm**Standard Box Quantity:** Each**Biatain Soft-Hold Foam Dressing Borderless 10x10cm****Manufacturer:** Coloplast **Manufacturer Code:** 3470

Features: Biatain Foam dressings provide superior exudate management for faster wound healing. It also offers high patient comfort due to a soft and flexible foam with bevelled edges. The superior absorption and retention properties of Biatain, minimises risk of leakage and maceration. Biatain offers a cost effective treatment due to long wear time. The semi-permeable top film is shower proof and provides a barrier function to bacteria.

Usage: Biatain Soft-Hold Foam Dressing is suitable for difficult to dress wounds including leg ulcers, pressure ulcers, non-infected diabetic foot ulcers, second degree burns, donor sites, postoperative wounds and skin abrasions. Biatain foam dressings may be used on patients who are in treatment for a local or systemic infection under the discretion of a physician. Suitable for use in combination with compression therapy.

Precautions:**Size:** 10x10cm**Standard Box Quantity:** 1

COL3472

Foam Non Adhesive Dressings

Biatain Soft-Hold Foam Dressing Borderless 10x20cm**Manufacturer:** Coloplast **Manufacturer Code:** 3472

Features: Biatain Foam dressings provide superior exudate management for faster wound healing. It also offers high patient comfort due to a soft and flexible foam with bevelled edges. The superior absorption and retention properties of Biatain, minimises risk of leakage and maceration. Biatain offers a cost effective treatment due to long wear time. The semi-permeable top film is shower proof and provides a barrier function to bacteria.

Usage: Biatain Soft-Hold Foam Dressing is suitable for difficult to dress wounds including leg ulcers, pressure ulcers, non-infected diabetic foot ulcers, second degree burns, donor sites, postoperative wounds and skin abrasions. Biatain foam dressings may be used on patients who are in treatment for a local or systemic infection under the discretion of a physician. Suitable for use in combination with compression therapy.

Precautions:**Size:** 10x20cm**Standard Box Quantity:** 1

COL4110

Foam Non Adhesive Dressings

Biatain IBU 10 x 10cm**Manufacturer:** Coloplast **Manufacturer Code:** 4110

Features: Superior exudate management and continuous release of ibuprofen. Provides excellent moist wound healing and may reduce wound pain caused by tissue damage. High patient comfort due to a soft and flexible foam dressing.

Usage: For exuding wounds including leg, pressure and non-infected foot ulcers. For second-degree burns, donor sites, postoperative wounds and abrasions. May be used throughout healing in the presence of exudate to provide protection for the indicated wound types. Suitable for use in combination with compression therapy.

Precautions: Precautions are used with the attendant liver and kidney diseases, chronic heart failure, with dyspeptic symptoms before treatment, immediately after surgery, with indications of a history of gastrointestinal bleeding, allergic reactions associated with NSAIDs.

Size: 10cm X 10cm**Standard Box Quantity:** Each

CON420633

Foam Non Adhesive Dressings

Aquacel Foam 10 x 10cm**Manufacturer:** ConvaTec **Manufacturer Code:** 420633

Features: Aquacel foam dressing consists of a soft absorbent foam pad with a waterproof top layer which protects the wound and helps to manage the moisture vapour from the exudate, and a Hydrofiber wound interface which creates a soft gel that contours to the wound bed, maintains a moist environment and aids in autolytic debridement without damaging newly formed tissue.

Usage: Primary or secondary dressing for both chronic and acute wounds, such as leg ulcers, pressure ulcers (Stages II-IV) and diabetic ulcers; surgical wounds; partial thickness burns; traumatic or surgical wounds left to heal by secondary intention such as dehiscent surgical incisions; dermatological and surgical incisions.

Precautions: The dressing can be cut to shape for convenience. Additional tape, or other means of fixation, will be required to secure the dressing.

Size: 10 x 10 cm**Standard Box Quantity:** 1

CON420635

Foam Non Adhesive Dressings

Aquacel Foam 15 x 15cm

Manufacturer: ConvaTec **Manufacturer Code:** 420635

Features: Aquacel foam dressing consists of a soft absorbent foam pad with a waterproof top layer which protects the wound and helps to manage the moisture vapour from the exudate, and a Hydrofiber wound interface which creates a soft gel that contours to the wound bed, maintains a moist environment and aids in autolytic debridement without damaging newly formed tissue.

Usage: Primary or secondary dressing for both chronic and acute wounds, such as leg ulcers, pressure ulcers (Stages II-IV) and diabetic ulcers; surgical wounds; partial thickness burns; traumatic or surgical wounds left to heal by secondary intention such as dehiscent surgical incisions; dermatological and surgical incisions.

Precautions: The dressing can be cut to shape for convenience. Additional tape, or other means of fixation, will be required to secure the dressing.

Size: 15 x 15 cm

Standard Box Quantity: 1



CON420637

Foam Non Adhesive Dressings

Aquacel Foam 15 x 20cm

Manufacturer: ConvaTec **Manufacturer Code:** 420637

Features: Aquacel foam dressing consists of a soft absorbent foam pad with a waterproof top layer which protects the wound and helps to manage the moisture vapour from the exudate, and a Hydrofiber wound interface which creates a soft gel that contours to the wound bed, maintains a moist environment and aids in autolytic debridement without damaging newly formed tissue.

Usage: Primary or secondary dressing for both chronic and acute wounds, such as leg ulcers, pressure ulcers (Stages II-IV) and diabetic ulcers; surgical wounds; partial thickness burns; traumatic or surgical wounds left to heal by secondary intention such as dehiscent surgical incisions; dermatological and surgical incisions.

Precautions: The dressing can be cut to shape for convenience. Additional tape, or other means of fixation, will be required to secure the dressing.

Size: 15 x 20 cm

Standard Box Quantity: 1



AS2100CH

Gauze & Non-Woven Dressings

Gauze 10 x 10cm - Unsterile (200/pkg)

Manufacturer: Cardinal **Manufacturer Code:** AS2100CH

Features: Gauze 10 x 10cm - Unsterile.

Usage: For general purpose use.

Precautions: Non-sterile.

Size: 10 x 10cm

Standard Box Quantity: PKG/200



CHC4847

Gauze & Non-Woven Dressings

Gauze, Packing, Sterile 1x36"

Manufacturer: Derma Science **Manufacturer Code:** 4847

Features: 100% Cotton packing gauze. Prewashed, Fluffed Gauze.

Usage:

Precautions: Sterile.

Size: 1 x 36 inch

Standard Box Quantity: Each



CHC4853

Gauze & Non-Woven Dressings

Gauze, Packing, Sterile 1/2x55

Manufacturer: Derma Science **Manufacturer Code:** 4853

Features: 100% Cotton packing gauze. Xray Detectable. Prewashed, Fluffed Gauze.

Usage:

Precautions: Sterile.

Size: 1/2 x 54 inch

Standard Box Quantity: Each



CHC4869

Gauze & Non-Woven Dressings

Vaginal Packing 5 x 182cm (roll)

Manufacturer: Derma Science **Manufacturer Code:** DUP4869

Features: Radiopaque gauze packing strip.

Usage: To pack fistulae, absorb drainage or aid in hemostasis. Often soaked with saline, antibiotic or antifungal cream before insertion.

Precautions: Place with care to avoid distention or compression of vasculature.

Size: 5x182cm

Standard Box Quantity: 1



CHC82408

Gauze & Non-Woven Dressings

Gauze, Sterile (2/Pkg) , Cotton 12ply 4x4"

Manufacturer: Derma Science **Manufacturer Code:** 82408

Features: These woven sponges are made from 100% cotton.

Usage: For general use or for debriding, covering and packing wounds.

Precautions: Sterile. Latex-free.

Size: 10 x 10 cm

Standard Box Quantity: Pack/2



CHCMDS21408

Gauze & Non-Woven Dressings

Gauze, Non-Sterile, 200/4x4" Pkg

Manufacturer: Dumex **Manufacturer Code:** MDS21408

Features: Gauze 10 x 10cm - Unsterile.

Usage:

Precautions: Non-sterile.

Size: 10x10cm

Standard Box Quantity: Pack/200



COV1050

Gauze & Non-Woven Dressings

Telfa Sponge 7.6 x 10cm - Sterile

Manufacturer: Covidien **Manufacturer Code:** 1050

Features: Telfa Sponge 7.6 x 10cm - Sterile. Superior "ouchless" Telfa dressing won't disrupt healing tissue by sticking to wound. Each sterile dressing individually wrapped in peel-open envelope.

Usage: Ideal as primary dressing for lightly draining wounds. Absorbent cotton pad. Bonded on both sides with perforated non-adherent film; can be cut to fit any shape without separating.

Precautions: Absorbent core is made out of 100% USP grade absorbent cotton.

Size: 7.6x10cm

Standard Box Quantity: 1



COV1238

Gauze & Non-Woven Dressings

Telfa Sponge 7.6 x 20cm - Sterile

Manufacturer: Covidien **Manufacturer Code:** 1238

Features: Telfa Sponge 7.6 x 20cm - Sterile. Superior "ouchless" Telfa dressing won't disrupt healing tissue by sticking to wound. Each sterile dressing individually wrapped in peel-open envelope.

Usage: Ideal as primary dressing for lightly draining wounds. Absorbent cotton pad. Bonded on both sides with perforated non-adherent film; can be cut to fit any shape without separating.

Precautions: Absorbent core is made out of 100% USP grade absorbent cotton.

Size: 7.6x20cm

Standard Box Quantity: 1



COV28959

Gauze & Non-Woven Dressings

Telfa 20 x 25cm

Manufacturer: Covidien **Manufacturer Code:** 28959

Features: Superior "ouchless" Telfa dressing won't disrupt healing tissue by sticking to wound. Each sterile dressing individually wrapped in peel-open envelope.

Usage: Ideal as primary dressing for lightly draining wounds. Absorbent cotton pad. Bonded on both sides with perforated non-adherent film; can be cut to fit any shape without separating.

Precautions: Absorbent core is made out of 100% USP grade absorbent cotton.

Size: 20x25cm

Standard Box Quantity: 1



DUM97820

Gauze & Non-Woven Dressings

Abdominal Pad 20 x 182cm (6ft section)

Manufacturer: Dumex **Manufacturer Code:** DUP97820

Features: This absorbent rayon/polyester pad will not interfere with granulating tissue by sticking to the wound. The soft perforated film allows air to circulate, and the absorbent pad allows the fluid to collect. It can easily be cut to fit any shape without separating.

Usage: Excellent for controlling bleeding and preventing contamination of large wounds.

Precautions: Non sterile. Latex-free.

Size: 20 x 182cm

Standard Box Quantity: 1



DUP82812

Gauze & Non-Woven Dressings

Gauze 10 x 20cm 12-PLY (2/pkg)

Manufacturer: Derma Science **Manufacturer Code:** DUP82812

Features: Ducare woven sponges are made from 12-ply 100% cotton. 2 per sterile package.

Usage: They are for general use or for debriding, covering and packing wounds.

Precautions: Latex-free. Sterile.

Size: 20cm x 10cm

Standard Box Quantity: Pk/2



DUP84122

Gauze & Non-Woven Dressings

Gauze Sponge Sterile 5 x 5cm (2/pkg)

Manufacturer: Derma Science **Manufacturer Code:** DUP84122

Features: Gauze sponge - Sterile 5 x 5cm 4-ply. 2-pack

Usage: These non-woven sponges are ideal for wound cleansing and also can be used as a secondary dressing for added absorbency

Precautions: Sterile. Latex free.

Size: 5x5cm

Standard Box Quantity: PKG/2



DUP90208PK50

Gauze & Non-Woven Dressings

Gauze, Non-Sterile, 50/2x2"

Manufacturer: Derma Science **Manufacturer Code:** 90208

Features: These woven sponges are made from 100% cotton.

Usage: For general use or for debriding, covering and packing wounds.

Precautions: Sterile. Latex-free.

Size: 5 x 5 cm

Standard Box Quantity: Pack/50



MC2102

Gauze & Non-Woven Dressings

Gauze 5 x 5cm - Unsterile (200/pkg)

Manufacturer: Disposamed **Manufacturer Code:** MC2102

Features: Gauze 5 x 5cm - Unsterile. 200 pack.

Usage: For general purpose use.

Precautions: Non-sterile.

Size: 5x5cm, 200-pack

Standard Box Quantity: PKG/200



MED2800

Gauze & Non-Woven Dressings

Gauze, Sterile Cotton (2/Pkg) 8ply 2x2"

Manufacturer: Medicom **Manufacturer Code:** 2800

Features: Provides optimal protection, comfort with highly absorbent soft sponges. Easy-to-use convenient individual peel-open package design. Worry-free with sterile heat-sealed envelope to insure safety.

Usage: To cover wounds, burns etc. Gauze can be used in various required situations.

Precautions:

Size: 5cm x 5cm

Standard Box Quantity: Pack/2



ML21444

Gauze & Non-Woven Dressings

Gauze Sponge Sterile 10 x 10cm (2/pkg)

Manufacturer: Medline **Manufacturer Code:** NON21444

Features: Strong, absorbent, rayon and polyester formed fabric sponges produce less lint & are virtually non-adherent for less trauma to fragile wound tissue. Faster wicking ability.

Usage: To cover wounds, burns etc. Gauze can be used in various required situations.

Precautions: Sterile. Latex free.

Size: 10x10cm

Standard Box Quantity: PKG/2



ML21450

Gauze & Non-Woven Dressings

Abdominal Pad 12.5 X 22.5cm

Manufacturer: Medline **Manufacturer Code:** NON21450

Features: This sterile, absorbent rayon/polyester pad will not interfere with granulating tissue by sticking to the wound. The soft perforated film allows air to circulate, and the absorbent pad allows the fluid to collect. It can easily be cut to fit any shape without separating. Latex-free.

Usage: Excellent for controlling bleeding and preventing contamination of large wounds.

Precautions: Contents sterile unless package is opened or damaged.

Size: 12.5 x 22.5cm

Standard Box Quantity: 1



ML21454

Gauze & Non-Woven Dressings

Abdominal Pad 20 x 25cm

Manufacturer: Medline **Manufacturer Code:** NON21454

Features: This sterile, absorbent rayon/polyester pad will not interfere with granulating tissue by sticking to the wound. The soft perforated film allows air to circulate, and the absorbent pad allows the fluid to collect. It can easily be cut to fit any shape without separating. Latex-free.

Usage: Excellent for controlling bleeding and preventing contamination of large wounds.

Precautions: Contents sterile unless package is opened or damaged.

Size: 20 x 25cm

Standard Box Quantity: 1



ML255015

Gauze & Non-Woven Dressings

Gauze Packing 2.5 x 457cm (btl) 1"

Manufacturer: Medline **Manufacturer Code:** NON255015

Features: Sterile packing strips are ideal for use as nasal packing, drainage wicks, or for packing wounds.

Usage: All are packaged as one strip per bottle.

Precautions: Keep container closed.

Size: 2.5x457cm

Standard Box Quantity: 1



ML255125

Gauze & Non-Woven Dressings

Gauze Packing 1.2 x 457cm (btl) ½"

Manufacturer: Medline **Manufacturer Code:** NON255125

Features: Sterile packing strips are ideal for use as nasal packing, drainage wicks, or for packing wounds.

Usage: All are packaged as one strip per bottle.

Precautions: Keep container closed.

Size: 1.2x457cm

Standard Box Quantity: 1



ML255145

Gauze & Non-Woven Dressings

Gauze Packing .6 x 457cm (btl) ¼"

Manufacturer: Medline **Manufacturer Code:** NON255145

Features: Sterile packing strips are ideal for use as nasal packing, drainage wicks, or for packing wounds.

Usage: All are packaged as one strip per bottle.

Precautions: Keep container closed.

Size: .6x457cm

Standard Box Quantity: 1



ML2560

Gauze & Non-Woven Dressings

Sterile Fenestrated Drain Sponge 10 x 10cm 6ply (2/Pkg)

Manufacturer: Medline **Manufacturer Code:** NON256000


Features: Drain and I.V. Sponges. Pre-cut notch conforms well around drains, catheters, chest tubes, I.V.'s and tracheostomies. Sterile.


Usage: Conforms well around drains, catheters, chest tubes, I.V.'s and tracheostomies.

Precautions: These sponges are virtually lint free which minimizes the risk of wound contamination.

Size: 10 x 10cm

Standard Box Quantity: 1



Trache Sponge 7.6 x 7.6cm (1/Pkg)**Manufacturer:** Acelity **Manufacturer Code:** M12507**Features:** Topper Dressing Sponges provide ideal protection and absorbency for light-to-moderately draining wounds.**Usage:** Multi-layered material for extra absorbency. Easy release and is softer for greater patient comfort.**Precautions:** Single use.**Size:** 7.6x7.6cm**Standard Box Quantity:** 1

3M1860S


Gloves, Gowns, Masks

High Filtration Mask For TB Small**Manufacturer:** 3M **Manufacturer Code:** 1860S**Features:** This small-sized health care particulate respirator and surgical mask helps provide respiratory protection against certain airborne biological particles. It is disposable and fluid resistant to splash and spatter of blood and other infectious material.**Usage:** This healthcare respirator is designed to help provide respiratory protection for the wearer. It meets CDC guidelines for Mycobacterium tuberculosis exposure control. As a disposable particulate respirator, it is intended to reduce wearer exposure to certain airborne particles including those generated by electrocautery, laser surgery, and other powered medical instruments.**Precautions:****Size:** Small**Standard Box Quantity:** 1

AMG018244

Gloves, Gowns, Masks

Masks - Standard Regular**Manufacturer:** AMG **Manufacturer Code:** 018-244**Features:** Pleat design prevents fluid pooling. Made of fluid resistant materials to help prevent contamination. Ideal for sensitive skin – latex-free & fiberglass-free. Smooth texture reduces the chance of skin irritation. Soft inner layer. High breathability. Odourless. Aluminum nosepiece for adjustable fit.**Usage:** Isolation mask.**Precautions:** Latex-free & fiberglass-free.**Size:** Regular**Standard Box Quantity:** 1

Gowns Non Sterile Long sleeved Med/Large**Manufacturer:** Medline **Manufacturer Code:** NON24258**Features:** For general protection in patient settings.**Usage:** One Size fits all**Precautions:****Size:** Universal**Standard Box Quantity:** 1**Sterile Small (pr)****Manufacturer:** Misc **Manufacturer Code:** 484405**Features:** Examination Gloves, Sterile, Latex Free Vinyl Small.**Usage:** For examination procedures that require latex free gloves.**Precautions:** Vinyl gloves should not be used when handling chemotherapy drugs.**Size:** Small**Standard Box Quantity:** PAIR**Sterile Medium (pr)****Manufacturer:** Misc **Manufacturer Code:** 484406**Features:** Examination Gloves, Sterile, Latex Free Vinyl Medium.**Usage:** For examination procedures that require latex free gloves.**Precautions:** Vinyl gloves should not be used when handling chemotherapy drugs.**Size:** Medium**Standard Box Quantity:** PAIR

MAX484407

Gloves, Gowns, Masks

Sterile Large (pr)**Manufacturer:** Misc **Manufacturer Code:** 484407**Features:** Examination Gloves, Sterile, Latex Free Vinyl Large.**Usage:** For examination procedures that require latex free gloves.**Precautions:** Vinyl gloves should not be used when handling chemotherapy drugs.**Size:** Large**Standard Box Quantity:** PAIR

MAX484802

Gloves, Gowns, Masks

Non Sterile Powder Free Chemo - Med (pr)**Manufacturer:** Misc **Manufacturer Code:** 484802**Features:** Non-Sterile Powder-Free Nitrile Exam Gloves have durable thin mil nitrile material that provides maximum protection against chemicals, harmful abrasions, Type I allergic reactions to NRL and viral exposure. Comfortable and secure even during the most difficult procedures.**Usage:** Glove is approved for chemotherapy treatments.**Precautions:** Latex free, non-sterile, single use.**Size:** Medium**Standard Box Quantity:** PAIR

MAX484803

Gloves, Gowns, Masks

Non Sterile Powder Free Chemo - Lg (pr)**Manufacturer:** Misc **Manufacturer Code:** 484803**Features:** Non-Sterile Powder-Free Nitrile Exam Gloves have durable thin mil nitrile material that provides maximum protection against chemicals, harmful abrasions, Type I allergic reactions to NRL and viral exposure. Comfortable and secure even during the most difficult procedures.**Usage:** Glove is approved for chemotherapy treatments.**Precautions:** Latex free, non-sterile, single use.**Size:** Large**Standard Box Quantity:** PAIR

MED1155B

Gloves, Gowns, Masks

Non Sterile Small (pr)

Manufacturer: Misc **Manufacturer Code:** 1155B

Features: Examination Gloves, Non-Sterile, Latex Free Vinyl Small.

Usage: For examination procedures that require latex free gloves.

Precautions: Vinyl gloves should not be used when handling chemotherapy drugs.

Size: Small

Standard Box Quantity: PAIR



MED1195CP

Gloves, Gowns, Masks

Non Sterile Medium (pr)

Manufacturer: Misc **Manufacturer Code:** 1195C

Features: Examination Gloves, Non-Sterile, Latex Free Vinyl Medium.

Usage: For examination procedures that require latex free gloves.

Precautions: Vinyl gloves should not be used when handling chemotherapy drugs.

Size: Medium

Standard Box Quantity: PAIR



MED1195DP

Gloves, Gowns, Masks

Non Sterile Large (pr)

Manufacturer: Misc **Manufacturer Code:** 1195D

Features: Examination Gloves, Non-Sterile, Latex Free Vinyl Large.

Usage: For examination procedures that require latex free gloves.

Precautions: Vinyl gloves should not be used when handling chemotherapy drugs.

Size: Large

Standard Box Quantity: PAIR



3M90001

Hydrocolloid Dressings

Tegaderm Oval Regular 7 x 9cm

Manufacturer: 3M **Manufacturer Code:** 90001

Features: This advanced formulation hydrocolloid creates ideal conditions for optimum wound healing. The polysaccharide particles in the dressing swell, absorbing exudate and producing a soft gel. A layer of 3M™ Tegaderm™ Film Dressing covers the hydrocolloid, extending to a border on oval dressings. This gives additional security against dressing edge roll or leakage.

Usage: For moderate to heavily exuding wounds. Gives up to a complete week's wear time.

Precautions: This product does not contain latex.

Size: 7x9cm

Standard Box Quantity: 1



3M90003

Hydrocolloid Dressings

Tegaderm Hydrocolloid 10cm x 12cm

Manufacturer: 3M **Manufacturer Code:** 90003

Features: Tegaderm Hydrocolloid consists of an oval polyurethane film, coated with acrylic adhesive. The film is impermeable to exudate and micro-organisms and extends past the border of an absorbent gelling hydrocolloid pad.

Usage: Allow a margin of 2.5 cm around the wound. Clean and dry wound and surrounding skin before application. The dressing, together with a secondary plastic film carrier, is removed from the silicone release liner using the red tabs provided, and applied to the wound site. Can be left on wound for up to 7 days.

Precautions: Not for infected wounds, wounds with depth greater than 1.5 cm or undermining greater than 1 cm - 1.5 cm. Not for full thickness burns or wounds with muscle, bone or tendon exposed or with eschar covering more than 50% of the wound.

Size: 10cm x 12cm

Standard Box Quantity: Each



3M90007

Hydrocolloid Dressings

Tegaderm Sacral 12 x 14cm

Manufacturer: 3M **Manufacturer Code:** 90007

Features: This advanced formulation hydrocolloid creates ideal conditions for optimum wound healing. The polysaccharide particles in the dressing swell, absorbing exudate and producing a soft gel. A layer of 3M™ Tegaderm™ Film Dressing covers the hydrocolloid, extending to a border on oval dressings. This gives additional security against dressing edge roll or leakage.

Usage: For moderate to heavily exuding wounds. Gives up to a complete week's wear time.

Precautions: This product does not contain latex.

Size: 12x14cm

Standard Box Quantity: 1



CON187901

Hydrocolloid Dressings

Duoderm Sterile EXTRA Thin Control Gel Formula 7.5cm x 7.5cm

Manufacturer: Convatec **Manufacturer Code:** 187901

Features: Duoderm is a hydrocolloid dressing which interacts with the wound exudates to form a gel with water. This latex free dressing has a triple hydrocolloid matrix with a viral and bacterial barrier.

Usage: The wound should be cleaned using sterile techniques. The paper backing is removed from an appropriate size dressing for the wound and is applied gently over the wound. There should be at least a 2 cm margin over the wound edges to allow for a secure sealing of the dressing.

Precautions: Dressing not to be used in the treatment of clinically infected wounds, wounds with depth greater than 1.5 cm or undermining greater than 1 cm - 1.5 cm. The dressing is also not indicated for the management of full thickness burns or wounds with muscle, bone or tendon exposed or with eschar covering more than 50% of the wound.

Size: 7.5cm x 7.5cm

Standard Box Quantity: 1



CON187955

Hydrocolloid Dressings

Duoderm Sterile EXTRA Thin Control Gel Formula 10 x 10cm

Manufacturer: Convatec **Manufacturer Code:** 187955

Features: Duoderm is a hydrocolloid dressing which interacts with the wound exudates to form a gel with water. This latex free dressing has a triple hydrocolloid matrix with a viral and bacterial barrier.

Usage: The wound should be cleaned using sterile techniques. The paper backing is removed from an appropriate size dressing for the wound and is applied gently over the wound. There should be at least a 2 cm margin over the wound edges to allow for a secure sealing of the dressing.

Precautions: Dressing not to be used in the treatment of clinically infected wounds, wounds with depth greater than 1.5 cm or undermining greater than 1 cm - 1.5 cm. The dressing is also not indicated for the management of full thickness burns or wounds with muscle, bone or tendon exposed or with eschar covering more than 50% of the wound.

Size: 10cm x 10cm

Standard Box Quantity: 1



CON187957

Hydrocolloid Dressings

Duoderm Sterile EXTRA Thin Control Gel Formula 15 x 15cm

Manufacturer: Convatec **Manufacturer Code:** 187957

Features: Duoderm is a hydrocolloid dressing which interacts with the wound exudates to form a gel with water. This latex free dressing has a triple hydrocolloid matrix with a viral and bacterial barrier.

Usage: The wound should be cleaned using sterile techniques. The paper backing is removed from an appropriate size dressing for the wound and is applied gently over the wound. There should be at least a 2 cm margin over the wound edges to allow for a secure sealing of the dressing.

Precautions: Dressing not to be used in the treatment of clinically infected wounds, wounds with depth greater than 1.5 cm or undermining greater than 1 cm - 1.5 cm. The dressing is also not indicated for the management of full thickness burns or wounds with muscle, bone or tendon exposed or with eschar covering more than 50% of the wound.

Size: 15cm x 15cm

Standard Box Quantity: 1



CON403326

Hydrocolloid Dressings

Duoderm 10 x 10cm

Manufacturer: ConvaTec **Manufacturer Code:** 403326

Features: DuoDERM Signal® dressing is an adhesive, tapered edge, modern hydrocolloid wound dressing with a change indicator. It incorporates a unique ConvaTec hydrocolloid formulation that distinguishes it from other hydrocolloid dressings. The adhesive hydrocolloid matrix and the outer film on DuoDERM Signal® dressing are permeable by moisture vapour and provide an occlusive barrier against bacteria, viruses, dirt and water.

Usage: DuoDERM Signal® Dressing is indicated for the management of lightly to moderately exuding wounds.

Precautions: This product does not contain latex.

Size: 10x10cm

Standard Box Quantity: 1



CON410501

Hydrocolloid Dressings

Duoderm Sacral 20 x 23cm

Manufacturer: ConvaTec **Manufacturer Code:** 410501

Features: DuoDERM Signal® dressing is an adhesive, tapered edge, modern hydrocolloid wound dressing with a change indicator. It incorporates a unique ConvaTec hydrocolloid formulation that distinguishes it from other hydrocolloid dressings. The adhesive hydrocolloid matrix and the outer film on DuoDERM Signal® dressing are permeable by moisture vapour and provide an occlusive barrier against bacteria, viruses, dirt and water.

Usage: DuoDERM Signal® Dressing is indicated for the management of lightly to moderately exuding wounds.

Precautions: This product does not contain latex.

Size: 20x23cm

Standard Box Quantity: 1



CON410510

Hydrocolloid Dressings

Duoderm Oval 11 x 19cm

Manufacturer: ConvaTec **Manufacturer Code:** 410510

Features: DuoDERM Signal® dressing is an adhesive, tapered edge, modern hydrocolloid wound dressing with a change indicator. It incorporates a unique ConvaTec hydrocolloid formulation that distinguishes it from other hydrocolloid dressings. The adhesive hydrocolloid matrix and the outer film on DuoDERM Signal® dressing are permeable by moisture vapour and provide an occlusive barrier against bacteria, viruses, dirt and water.

Usage: DuoDERM Signal® Dressing is indicated for the management of lightly to moderately exuding wounds.

Precautions: This product does not contain latex.

Size: 11x19cm

Standard Box Quantity: 1



CON177902

Hydrofibre Dressings

Aquacel 10x10cm

Manufacturer: Convatec **Manufacturer Code:** 177902

Features: Aquacel is a primary wound dressing made from sodium carboxymethylcellulose. It is produced as a textile fibre and presented in the form of a fleece held together by a needle bonding process, and is available both as a 'ribbon' for packing cavities, and as a flat non-woven pad for application to larger open wounds.

Usage: Aquacel sheets may be applied to exuding lesions including leg ulcers, pressure areas, donor sites, and most other granulating wounds, but for deeper cavity wounds and sinuses the ribbon packing is generally preferred. The dressing is easy to remove without causing pain or trauma, and leaves minimal residue on the surface of the wound.

Precautions: Although there are no known contra-indications to the use of Aquacel, the dressing will be of little value if applied to wounds that are very dry, or covered with hard black necrotic tissue.

Size: 10cm x 10cm

Standard Box Quantity: Each



CON420127

Hydrofibre Dressings

Aquacel Rope 1 x 45cm

Manufacturer: ConvaTec **Manufacturer Code:** 420127

Features: AQUACEL® dressing incorporating unique Hydrofiber® Technology is a sterile, white, fibrous dressing derived from 100% sodium carboxymethylcellulose. It is a versatile, primary dressing indicated for use on moderately and highly exuding chronic and acute wounds.

Usage: Indicated for: excised abscesses, excised pilonidal cysts, plus tunnelling and undermining wounds, cavity wounds, sinus tracts (post-procedure), suture lines, fingers and toes of paediatric patients (burns, abrasions), orthopaedic pin insertion sites, tube insertion site (g-tube, percutaneous endoscopic gastronomy), finger or toe degloving, circumcision sites.

Precautions: This product does not contain latex.

Size: 1 x 45cm

Standard Box Quantity: 1



CON420671

Hydrofibre Dressings

Aquacel Extra dressing, 5 x 5cm

Manufacturer: ConvaTec **Manufacturer Code:** 420671

Features: AQUACEL® EXTRA™ Hydrofiber® Wound Dressing with Strengthening Fiber is composed of two layers of Hydrofiber® Technology. Cleared for the same indications as original AQUACEL® dressing. Specifically suitable for moderate to highly exuding wounds.

Usage: Leg ulcers, pressure ulcers (stage II-IV) and diabetic ulcers; surgical wounds; partial-thickness burns; traumatic or surgical wounds left to heal by secondary intention; surgical wounds that heal by primary intent; traumatic wounds; wounds that are prone to bleeding; painful wounds.

Precautions: Do not use AQUACEL® EXTRA™ Hydrofiber® Wound Dressing with Strengthening Fiber if you are sensitive to the dressing, or its components, or have had an allergic reaction to the dressing.

Size: 5 x 5cm

Standard Box Quantity: 1



CON420672

Hydrofibre Dressings

Aquacel Extra dressing, 10 x 10cm

Manufacturer: ConvaTec **Manufacturer Code:** 420672

Features: AQUACEL® EXTRA™ Hydrofiber® Wound Dressing with Strengthening Fiber is composed of two layers of Hydrofiber® Technology. Cleared for the same indications as original AQUACEL® dressing. Specifically suitable for moderate to highly exuding wounds.

Usage: Leg ulcers, pressure ulcers (stage II-IV) and diabetic ulcers; surgical wounds; partial-thickness burns; traumatic or surgical wounds left to heal by secondary intention; surgical wounds that heal by primary intent; traumatic wounds; wounds that are prone to bleeding; painful wounds.

Precautions: Do not use AQUACEL® EXTRA™ Hydrofiber® Wound Dressing with Strengthening Fiber if you are sensitive to the dressing, or its components, or have had an allergic reaction to the dressing.

Size: 10 x 10cm

Standard Box Quantity: 1



CON420673

Hydrofibre Dressings

Aquacel Extra dressing, 15 x 15cm

Manufacturer: ConvaTec **Manufacturer Code:** 420673

Features: AQUACEL® EXTRA™ Hydrofiber® Wound Dressing with Strengthening Fiber is composed of two layers of Hydrofiber® Technology. Cleared for the same indications as original AQUACEL® dressing. Specifically suitable for moderate to highly exuding wounds.

Usage: Leg ulcers, pressure ulcers (stage II-IV) and diabetic ulcers; surgical wounds; partial-thickness burns; traumatic or surgical wounds left to heal by secondary intention; surgical wounds that heal by primary intent; traumatic wounds; wounds that are prone to bleeding; painful wounds.

Precautions: Do not use AQUACEL® EXTRA™ Hydrofiber® Wound Dressing with Strengthening Fiber if you are sensitive to the dressing, or its components, or have had an allergic reaction to the dressing.

Size: 15 x 15cm

Standard Box Quantity: 1



SN66000324

Hydrogels

Intrasite Conformable 10 x 10cm

Manufacturer: Smith & Nephew **Manufacturer Code:** 66000324

Features: IntraSite® Conformable is a conformable hydrogel dressing which combines the advantages of IntraSite Gel with a non-woven dressing to aid in the gentle packing of deep shallow, open or undermined wounds. The non-woven substrate ensures the gel remains in close contact with the entire wound surface. It is suitable for use in necrotic, sloughy and granulating wounds, and helps create a moist wound environment to facilitate healing.

Usage: IntraSite Conformable is indicated for use in open wounds e.g. pressure sores, surgical wounds healing by secondary intention and malignant wounds. It is especially useful for the packing of deep cavity wounds. It can be used on infected wounds under medical supervision.

Precautions: Known sensitivity to IntraSite Conformable or any of its ingredients. IntraSite Conformable should be used with care in the vicinity of the eye. Avoid ingestion. Avoid use on cracked nipples.

Size: 10x10cm

Standard Box Quantity: 1



Intrasite Gel 8 gm**Manufacturer:** Smith & Nephew **Manufacturer Code:** 7308**Features:** IntraSite Gel is an amorphous hydrogel which promotes debridement of necrotic tissue. Loosens and absorbs slough and exudate. Provides a moist environment during later stages of wound healing. Non-adherent and does not harm viable or periwound tissue.**Usage:** For pressure sores, leg ulcers, surgical and malignant wounds, partial thickness burns, scalds, lacerations and grazes, and for treatment of granulating cavity wounds. Use on infected wounds under medical supervision.**Precautions:** Known sensitivity to IntraSite Gel or any of its ingredients. Use with care in the vicinity of the eyes and in deep wounds with narrow openings (e.g. fistulas) where removal of gel may be difficult. Should not be taken internally.**Size:** 8 g**Standard Box Quantity:** 1**Intrasite Gel 15gr****Manufacturer:** Smith & Nephew **Manufacturer Code:** 7311**Features:** IntraSite Gel is an amorphous hydrogel which promotes debridement of necrotic tissue. Loosens and absorbs slough and exudate. Provides a moist environment during later stages of wound healing. Non-adherent and does not harm viable or periwound tissue.**Usage:** For pressure sores, leg ulcers, surgical and malignant wounds, partial thickness burns, scalds, lacerations and grazes, and for treatment of granulating cavity wounds. Use on infected wounds under medical supervision.**Precautions:** Known sensitivity to IntraSite Gel or any of its ingredients. Use with care in the vicinity of the eyes and in deep wounds with narrow openings (e.g. fistulas) where removal of gel may be difficult. Should not be taken internally.**Size:** 15g**Standard Box Quantity:** Each**Intrasite Gel 25g****Manufacturer:** Smith & Nephew **Manufacturer Code:** 7313**Features:** IntraSite Gel is an amorphous hydrogel which promotes debridement of necrotic tissue. Loosens and absorbs slough and exudate. Provides a moist environment during later stages of wound healing. Non-adherent and does not harm viable or periwound tissue.**Usage:** For pressure sores, leg ulcers, surgical and malignant wounds, partial thickness burns, scalds, lacerations and grazes, and for treatment of granulating cavity wounds. Use on infected wounds under medical supervision.**Precautions:** Known sensitivity to IntraSite Gel or any of its ingredients. Use with care in the vicinity of the eyes and in deep wounds with narrow openings (e.g. fistulas) where removal of gel may be difficult. Should not be taken internally.**Size:** 25g**Standard Box Quantity:** 1

3-way Stop Cock**Manufacturer:** Med-Rx **Manufacturer Code:** 10-1006**Features:** Three way stopcock with Luer Lock connectors**Usage:****Precautions:****Size:****Standard Box Quantity:** Each**Opsite Iv3000 10cmx12cm****Manufacturer:** Smith & Nephew **Manufacturer Code:** 4008**Features:** IV3000 is a dedicated range of IV dressings designed to meet the demands of modern IV site care. IV3000 offers several important benefits: IV3000 REACTICT film and unique grid pattern adhesive allowing optimum moisture vapour transfer to reduce moisture accumulation and bacterial growth; Waterproof, allowing the patient to shower with the dressing in-situ.**Usage:** IV3000 is indicated for peripheral and central venous catheter fixation.**Precautions:** IV3000 should be used on IV sites only. Do not use for wounds such as pressure sores, burns or donor sites. Use on prepared dry skin only. Do not stack dressings or allow dressings to overlap. Periodically monitor the dressing and catheter site to confirm secure attachment and proper infusion, especially after bathing, showering, or if the dressing and catheter site becomes wet. Apply and remove carefully from sensitive or fragile skin.**Size:** 10 x 12cm**Standard Box Quantity:** Each**pk/2 Plastic Cotton Tipped Applicator Sterile 6 inch****Manufacturer:** AMG **Manufacturer Code:** 018-435**Features:** Sterile 6 inch plastic applicator stick with firmly wound cotton bud.**Usage:****Precautions:** Single use only. 100% Cotton tip.**Size:** 6 inch**Standard Box Quantity:** PKG/2

AMG118395

Miscellaneous

Silver Nitrate Sticks

Manufacturer: AMG **Manufacturer Code:** 118-395

Features: Flexible Caustic Applicators, 6 inch, Silver Nitrate 75% Potassium Nitrate 25%.

Usage: Silver nitrate applicators are useful for cauterisation of skin or mucous membrane and for the removal of granulation tissue, warts and verrucae.

Precautions: Exposure to light will cause the silver nitrate tip to turn black, but will not affect the product's potency. KEEP OUT OF REACH OF CHILDREN. DO NOT USE IN OR NEAR EYES. Silver nitrate applicators must not be used on genital warts. Continued application to mucous membranes and open wounds leads to argyria, a bluish-black discolouration of the skin due to depositions of granules of silver compounds in the connective tissues.

Size: 6 inches

Standard Box Quantity: 1



AMG500515

Miscellaneous

Scalpel #15 With Handle

Manufacturer: AMG **Manufacturer Code:** 500-515

Features: Disposable Scalpel #15 With Handle. Small, precision knives. Made of stainless steel.

Usage: The fine blade is used to make precise cuts.

Precautions: Sterile, the fine blade is very sharp.

Size: #15

Standard Box Quantity: 1



BAX25295A

Miscellaneous

Pipe Cleaners 125cm

Manufacturer: Misc **Manufacturer Code:**

Features: Pipe Cleaners 125cm.

Usage: General Use

Precautions: Non-sterile.

Size: 125cm

Standard Box Quantity: 1



CHC12240

Miscellaneous

Toothettes

Manufacturer: Kimberly Clark **Manufacturer Code:** KC12240

Features: Oral Care Swab with reinforced plastic handle.

Usage: Oral Swabs are designed to assist the practitioner in providing proper oral care.

Precautions: Latex Free

Size:

Standard Box Quantity: 1



CHC854201

Miscellaneous

Littauer Suture Scissors

Manufacturer: Med-Rx **Manufacturer Code:** 85-4201

Features: Disposable Littauer scissors with hook type blade for lifting and cutting sutures.

Usage: For suture removal.

Precautions: Sterile. Latex free.

Size:

Standard Box Quantity: Each



CHC854206

Miscellaneous

Trache Brush Large

Manufacturer: Med-Rx **Manufacturer Code:** 85-4206

Features: Trache Brush Large

Usage: Indicated for cleaning cannulae to prevent obstruction.

Precautions: Follow manufacturers cleaning instructions to avoid damaging cannulae.

Size: Large

Standard Box Quantity: 1



CHC9021

Miscellaneous

Specimen Container, C&S, Urine Or Sputum

Manufacturer: NCS Diagnostics **Manufacturer Code:** NCS902-1W


Features: Shatter-resistant containers feature a proprietary, all plastic, leak-resistant seal. They are tested at varying pressures, making them ideal for use in various transport systems.

Usage: For urine samples.

Precautions: Sterile

Size: 90mL

Standard Box Quantity: Each



CM0020

Miscellaneous

Paracentesis bag 4000ml

Manufacturer: Cardiomed **Manufacturer Code:** CM-0020

Features: Paracentesis bag 4000ml

Usage: Designed for patients who require ascites drainage.

Precautions: Sterile. Latex free.

Size: 4000ml

Standard Box Quantity: 1



COR207366

Miscellaneous

Nasogastric Tube 6FR w 1.5gr stylet 36 inch

Manufacturer: CorPak **Manufacturer Code:** COR20-7366

Features: Flow-through stylets to flush, auscultate and aspirate during the insertion procedure. Water activated lubricant for greater patient comfort during intubation and easy stylet removal. Radiopaque. "Y" ports minimize touch contamination and make flushing easy.

Usage: NG tubes can be used for enteral feeding.

Precautions: Complications may include aspiration and tissue trauma. Placement of the catheter can induce gagging or vomiting, therefore suction should always be ready to use. Gloves must be worn and if the risk of vomiting is high, the operator and team members should consider face and eye protection as well as a gown.

Size: 6FR 91CM/36 inches

Standard Box Quantity: 1



COR208366

Miscellaneous

Nasogastric Pediatric Tube 6FR w 1gr stylet 36 inch

Manufacturer: CorPak **Manufacturer Code:** 20-8366

Features: Flow-through stylets to flush, auscultate and aspirate during the insertion procedure. Water activated lubricant for greater patient comfort during intubation and easy stylet removal. Radiopaque. "Y" ports minimize touch contamination and make flushing easy.

Usage: NG tubes can be used for enteral feeding.

Precautions: Complications may include aspiration and tissue trauma. Placement of the catheter can induce gagging or vomiting, therefore suction should always be ready to use. Gloves must be worn and if the risk of vomiting is high, the operator and team members should consider face and eye protection as well as a gown.

Size: 6FR 91CM/36 inches

Standard Box Quantity: 1



COR209368

Miscellaneous

Nasogastric Tube 8FR w stylet 36 inch

Manufacturer: CorPak **Manufacturer Code:** COR20-9368

Features: Flow-through stylets to flush, auscultate and aspirate during the insertion procedure. Water activated lubricant for greater patient comfort during intubation and easy stylet removal. Radiopaque. "Y" ports minimize touch contamination and make flushing easy.

Usage: NG tubes can be used for enteral feeding.

Precautions: Complications may include aspiration and tissue trauma. Placement of the catheter can induce gagging or vomiting, therefore suction should always be ready to use. Gloves must be worn and if the risk of vomiting is high, the operator and team members should consider face and eye protection as well as a gown.

Size: 8FR 36 inches

Standard Box Quantity: 1



DS600

Miscellaneous

Nephrostomy bag 600ml

Manufacturer: Dartin Medical **Manufacturer Code:** DS600

Features: The Dartin Select 600 is one of the most patient friendly drainage bags available. It is a convenient 600 ml capacity, has a flock backing for enhanced comfort when worn next to the skin, has a 50cm silicone tube terminating in a male luer lock to allow a secure and kink free connection to a nephrostomy or other percutaneously inserted drainage catheter, has an easy to use tap at the bottom for emptying, and has a PVC extension below the tap to allow the connection of a larger bag for night time use.

Usage: Dartin Select bag is designed for the collection of fluid during nephrostomy, biliary, and abscess drainages, with particular attention being paid to achieving maximum procedural efficiency and patient comfort.

Precautions: DEHP & Latex Free

Size: 600ml

Standard Box Quantity: 1



DUM680

Miscellaneous

Niko Drain Fix 5-14FR

Manufacturer: Unomedical **Manufacturer Code:** MK-680M


Features: Drain-Fix is a sterile, securement device for percutaneous catheters and tubes. It is specifically designed to stop migration, movement and accidental removal of the catheter.

Usage: For Securement of Pigtail catheters, Balloon catheters, Suprapubic catheters, PEG feeding tubes and other percutaneous catheters and drains.

Precautions: Designed for 7 days in situ.

Size: 5-14FR

Standard Box Quantity: 1



IBM199129A

Miscellaneous

Sterile Safety Pin 1.5in

Manufacturer: Stevens **Manufacturer Code:** 347-30-091

Features: Sterile Safety Pin 2 per package.

Usage: General Use

Precautions: Sterile.

Size: 1.5 inches

Standard Box Quantity: PKG/2



ML1674

Miscellaneous

Lymphedema Drainage Bag 2L

Manufacturer: Medline **Manufacturer Code:** DYNC1674

Features: Drain bags come with a reinforced hanger with built-in anti-kink drain tube guide. 50" of 11/32" STAR tubing facilitates drainage. Adjustable clamp secures tubing and needle-less sample port that is accessible with either a blunt cannula or needle. 2000mL capacity.

Usage: These units are latex free and in sterile packaging. They are considered one time use or disposable.

Precautions: Monitor bag and empty to avoid overflow.

Size: 2000 cc

Standard Box Quantity: 1



Tongue Depressors**Manufacturer:** Medline **Manufacturer Code:** MDS202073**Features:** Wooden Tongue Depressor.**Usage:** Useful for spreading ointments and for stirring and mixing medications.**Precautions:** Sterile.**Size:** 5-1/2 inch**Standard Box Quantity:** 1**Mouth Spray (Moi-Stir)****Manufacturer:** Pharma Science **Manufacturer Code:** 1439**Features:** Relieves dry mouth caused by certain medication, therapeutic treatments, diseases, and advanced age.**Usage:** Spray 1 or 2 times into mouth whenever it feels uncomfortably dry.**Precautions:** Keep out of reach of children.**Size:** 120mL**Standard Box Quantity:** Each**Non-Sterile plastic disposable A-clamp (CVL Hickman Clamp)****Manufacturer:** AMG **Manufacturer Code:** MSA-110**Features:** Reusable locking plastic clamp.**Usage:** For closing Hickman catheter lines to prevent entry of air.**Precautions:** Apply clamp between skin and break or disconnection point.**Size:****Standard Box Quantity:** 1

MSTR012RSC

Miscellaneous

Nephrostomy 14FR Connecting Tube

Manufacturer: Medical Systems **Manufacturer Code:** TR-012RSC

Features: Nephrostomy 14FR Connecting Tube

Usage: Connecting Tube 30cm, with Blue Lip and Removable Low Pressure Stop Cock.

Precautions: Latex free. Non DEHP.

Size: 14FR

Standard Box Quantity: 1



POS8197L

Miscellaneous

Large Foam Trache Ties 33-48cm

Manufacturer: Posey **Manufacturer Code:** 8197L

Features: Soft foam padding with hook & loop "double track" closure. Helps prevent extubation.

Usage: A convenient, comfortable collar to secure tracheostomy tubes. May be used with most tracheostomy tubes.

Precautions: Do not attach too tightly. Foam should not be compressed when applied.

Size: Large

Standard Box Quantity: 1



POS8197M

Miscellaneous

Medium Trache Ties

Manufacturer: Posey **Manufacturer Code:** 8197M

Features: Soft foam padding with hook & loop "double track" closure. Helps prevent extubation.

Usage: A convenient, comfortable collar to secure tracheostomy tubes. May be used with most tracheostomy tubes.

Precautions: Do not attach too tightly. Foam should not be compressed when applied.

Size: Medium

Standard Box Quantity: 1



Phydrion Brilliant PH stik Ranges 1-14**Manufacturer:** Mikro Essential Labs **Manufacturer Code:** 7800**Features:** A simple single colour match for all pH. Distinct clear colour separations. Immediate response.**Usage:** Immediate pH testing is available in brilliant, long lasting colors. These direct-reading pH papers give a simple, single color match for all pH values in distinct, bright color separations.**Precautions:** Use one strip per each solution.**Size:****Standard Box Quantity:** 1**Lymphedema Catheter Butterfly 19G Cadium****Manufacturer:** Terumo **Manufacturer Code:** TC-SVS19BLS**Features:** Lymphedema Catheter Butterfly 19G Cadium**Usage:****Precautions:** Single use.**Size:** 19G**Standard Box Quantity:** 1**Scissors - Sterile****Manufacturer:** Med-Rx **Manufacturer Code:** 85-4202**Features:** Scissors - Sterile sharp/blunt, disposable.**Usage:** Clinical Instrument, for wound care treatment.**Precautions:** Sterile.**Size:****Standard Box Quantity:** 1

ML250317

Non Adherent Dressing

Curad 7.5 x 40cm Non-Adherent Oil Emulsion Dressing

Manufacturer: Curad **Manufacturer Code:** CUR250317

Features: This sterile, rayon/polyester non adhesive petrolatum pad will not interfere with granulating tissue by sticking to the wound. The soft perforated film allows air to circulate, and the absorbent pad allows the fluid to collect. It can easily be cut to fit any shape without separating.

Usage: Place pads over lightly bleeding post-surgical wounds.

Precautions: Latex-free.

Size: 7.5 x 40cm

Standard Box Quantity: 1



ML250330

Non Adherent Dressing

Curad 7.5 x 7.5cm Non-Adherent Oil Emulsion Dressing

Manufacturer: Curad **Manufacturer Code:** CUR250330

Features: This sterile, rayon/polyester non adhesive petrolatum pad will not interfere with granulating tissue by sticking to the wound. The soft perforated film allows air to circulate, and the absorbent pad allows the fluid to collect. It can easily be cut to fit any shape without separating.

Usage: Place pads over lightly bleeding post-surgical wounds.

Precautions: Latex-free.

Size: 7.5cm x 7.5cm

Standard Box Quantity: 1



ML250381

Non Adherent Dressing

Curad 7.5 x 20cm Non-Adherent Oil Emulsion Dressing

Manufacturer: Curad **Manufacturer Code:** CUR250381

Features: This sterile, rayon/polyester non adhesive petrolatum pad will not interfere with granulating tissue by sticking to the wound. The soft perforated film allows air to circulate, and the absorbent pad allows the fluid to collect. It can easily be cut to fit any shape without separating.

Usage: Place pads over lightly bleeding post-surgical wounds.

Precautions: Latex-free.

Size: 20cm x 7.5cm

Standard Box Quantity: 1



Adaptic 7.5x7.5cm**Manufacturer:** Acelity **Manufacturer Code:** 2012**Features:** Ideal for most draining wounds. Coated with a petrolatum emulsion to prevent dressing adherence while allowing exudate to easily pass through the absorbent layer. Sterile .**Usage:** ADAPTIC® Non-adhering Dressing Features: Sterile, Protects regenerating tissue and minimizes patient pain during dressing changes, Exudate easily flows through to second absorbent dressing, Can be cut to fit the wound size .**Precautions:** Non adhering.**Size:** 7.5x7.5cm**Standard Box Quantity:** Each**Adaptic 7.5x20cm****Manufacturer:** Acelity **Manufacturer Code:** 2013**Features:** Ideal for most draining wounds. Coated with a petrolatum emulsion to prevent dressing adherence while allowing exudate to easily pass through the absorbent layer. Sterile .**Usage:** ADAPTIC® Non-adhering Dressing Features: Sterile, Protects regenerating tissue and minimizes patient pain during dressing changes, Exudate easily flows through to second absorbent dressing, Can be cut to fit the wound size .**Precautions:** Non adhering.**Size:** 7.6 X 20.3 cm**Standard Box Quantity:** Pack/3**Adaptic Touch 12.5cm x 15cm****Manufacturer:** Acelity **Manufacturer Code:** TCH503**Features:** ADAPTIC TOUCH® is a flexible, open-mesh primary wound contact layer comprised of cellulose acetate coated with a soft tack silicone designed to facilitate fluid transfer to secondary dressing and minimize adherence and pain at dressing change.**Usage:** For use in the management of dry to heavily exuding, partial and full-thickness chronic wounds including venous ulcers, decubitus (pressure) ulcers and diabetic ulcers, and for traumatic and surgical wounds, donor sites and 1st and 2nd degree burns.**Precautions:** May be left in place for several days depending upon amount of exudate.**Size:** 12.7x15cm**Standard Box Quantity:** Each

Speedicath Pediatric Intermittent Catheters Female 8FR (2.7mm)**Manufacturer:** Coloplast **Manufacturer Code:** 27708**Features:** Speedicath pediatric Nelaton intermittent catheter with polished eyelets. Ships packaged in saline for immediate use.**Usage:** For pediatric use.**Precautions:** Latex, PVC and phthalate free.**Size:** 8FR 8 inches**Standard Box Quantity:** 1

ANE507235

Pleural Drainage

PleurX Replacement Valve Cap**Manufacturer:** Care Fusion **Manufacturer Code:** 50-7235**Features:** PleurX replacement valve cap.**Usage:** Used with PleurX kits (ordered separately)**Precautions:****Size:****Standard Box Quantity:** Each

ANE507510

Pleural Drainage

PleurX Drainage Kit With 1000mL Vacuum Bottle**Manufacturer:** Care Fusion **Manufacturer Code:** 50-7510**Features:** PleurX drainage kit with 1000ml vacuum bottle, drainage line, Self-adhesive dressing, Blue wrapping, Alcohol pads (qty: 3), Pair of gloves, Valve cap, Blue emergency clamp, 4 Gauze pads, Foam catheter pad**Usage:** The PleurX catheter system is for home management of recurrent pleural effusions and malignant ascites. With the PleurX in-dwelling catheter and vacuum bottles, patients can manage fluid accumulation through intermittent drainages at home, with minimal physician intervention.**Precautions:** Consult frequently asked questions brochure. If drainage slows or declines, call your doctor. Redness, warmth to touch, swelling, fever or drainage from around the catheter site may indicate catheter site is infected. If catheter pulls out or is cut, clamp shut and seek immediate medical attention.**Size:** 1000ml**Standard Box Quantity:** Each

Pleural Drainage Kit (Adaptor/Dressing Kit & 1000ml btl)**Manufacturer:** Medquest **Manufacturer Code:** KT5000

Features: Bottle includes anti-reflux valve and green accordion indicator to show whether the bottle has vacuum pressure. Bottle includes graduated markings and a pour spout to allow disposal of fluid by snipping off the accordion indicator. Adaptor line allows attachment to PleurX catheter. Roller clamp allows adjustment of drainage rate.

Usage: For pleural drainage. Manufacturer training required for proper use of this product. Follow attachment and clamping instructions closely to avoid accidental loss of vacuum pressure, or introduction of air into the drainage catheter. For supporting documents please refer to the CCAC Partner Site and the last page of this catalogue.

Precautions: If drainage slows or declines, call your doctor. Redness, warmth to touch, swelling, fever or drainage from the catheter site may indicate infection. If catheter pulls out or is cut, clamp shut and seek immediate medical attention.

Size: 1000ml**Standard Box Quantity:** 1**Pleural Procedure Dressing Kit****Manufacturer:** Medquest **Manufacturer Code:** 2976C

Features: MedQuest pleural procedure dressing kit contains: gloves, fenestrated foam pad, Tegaderm dressing, emergency clamp, gauze pads, alcohol wipes, Cavilon barrier wipe, valve cap for MedQuest catheter.

Usage: For pleural drainage procedures.

Precautions: Use caution to protect the sterility of kit contents after opening.

Size: Procedure Kit**Standard Box Quantity:** 1**MedQuest 600ml Pleural Drainage Bottle with Vacuum****Manufacturer:** Medquest **Manufacturer Code:** S55120

Features: Bottle includes anti-reflux valve and green accordion indicator to show whether the bottle has vacuum pressure. Bottle includes graduated markings and a pour spout to allow disposal of fluid by snipping off the accordion indicator. Roller clamp allows adjustment of drainage rate.

Usage: For pleural drainage. Manufacturer training required for proper use of this product. Follow attachment and clamping instructions closely to avoid accidental loss of vacuum pressure, or introduction of air into the drainage catheter. For supporting documents please refer to the CCAC Partner Site and the last page of this catalogue.

Precautions: If drainage slows or declines, call your doctor. Redness, warmth to touch, swelling, fever or drainage from the catheter site may indicate infection. If catheter pulls out or is cut, clamp shut and seek immediate medical attention.

Size: 600mL**Standard Box Quantity:** 1

SYSM772028

Protease Inhibitors

Promogran 28cm²

Manufacturer: Acelity **Manufacturer Code:** M772028

Features: PROMOGRAN™ matrix is a topically applied interactive wound therapy. It is a sterile, freeze dried composite of oxidised regenerated cellulose (ORC) and collagen. PROMOGRAN™ matrix modulates and re-balances the wound environment by binding and inactivation of proteases and protection of naturally occurring growth factors against degradation by these excess proteases.

Usage: Indicated for the management of wounds healing by secondary intent which are clear of necrotic tissue, including: Diabetic ulcers, Venous ulcers, Pressure ulcers, ulcers caused by mixed vascular aetiologies, Traumatic and surgical wounds.

Precautions: PROMOGRAN™ matrix is contraindicated in patients with known hypersensitivity to ORC and Collagen. Discontinue use if signs of sensitivity appear. If infection is suspected an appropriate antimicrobial dressing or systemic therapy should be used.

Size: 28cm

Standard Box Quantity: 1



SYSM772123

Protease Inhibitors

Promogran 123cm²

Manufacturer: Acelity **Manufacturer Code:** M772123

Features: PROMOGRAN™ matrix is a topically applied interactive wound therapy. It is a sterile, freeze dried composite of oxidised regenerated cellulose (ORC) and collagen. PROMOGRAN™ matrix modulates and re-balances the wound environment by binding and inactivation of proteases and protection of naturally occurring growth factors against degradation by these excess proteases.

Usage: Indicated for the management of wounds healing by secondary intent which are clear of necrotic tissue, including: Diabetic ulcers, Venous ulcers, Pressure ulcers, ulcers caused by mixed vascular aetiologies, Traumatic and surgical wounds.

Precautions: PROMOGRAN™ matrix is contraindicated in patients with known hypersensitivity to ORC and Collagen. Discontinue use if signs of sensitivity appear. If infection is suspected an appropriate antimicrobial dressing or systemic therapy should be used.

Size: 123cm

Standard Box Quantity: 1



SYSPS2028

Protease Inhibitors

Prisma Silver Dressing 28cm²

Manufacturer: Acelity **Manufacturer Code:** PS2028

Features: Promogran Prisma Wound Balancing Matrix is comprised of a sterile, freeze dried composite of oxidized regenerated cellulose (ORC), collagen and 1% silver-ORC which contains 25% w/w ionically bound antimicrobial silver.

Usage: Prisma is indicated for the management of all wounds healing by secondary intent which are clear of necrotic tissue including: Diabetic ulcers; Venous ulcers; Pressure ulcers; Ulcers caused by mixed vascular aetiologies; Traumatic wounds and surgical wounds. Prisma has known haemostatic properties and can be used under compression therapy.

Precautions: Contraindicated in patients with extensive burns or known hypersensitivity to ORC, Collagen and Silver. Discontinue use if signs of sensitivity appear.

Size: 28cm²

Standard Box Quantity: 1



Prisma Silver Dressing 123cm²**Manufacturer:** Acelity **Manufacturer Code:** PS2123**Features:** Promogran Prisma Wound Balancing Matrix is comprised of a sterile, freeze dried composite of oxidized regenerated cellulose (ORC), collagen and 1% silver-ORC which contains 25% w/w ionically bound antimicrobial silver.**Usage:** Prisma is indicated for the management of all wounds healing by secondary intent which are clear of necrotic tissue including: Diabetic ulcers; Venous ulcers; Pressure ulcers; Ulcers caused by mixed vascular aetiologies; Traumatic wounds and surgical wounds. Prisma has known haemostatic properties and can be used under compression therapy.**Precautions:** Contraindicated in patients with extensive burns or known hypersensitivity to ORC, Collagen and Silver. Discontinue use if signs of sensitivity appear.**Size:** 123cm²**Standard Box Quantity:** 1

DUP82910

Respiratory Supplies - Tracheotomy

Drain Sponge, 10x10cm Sterile 12-Ply Woven**Manufacturer:** Derma Science **Manufacturer Code:** DUP82910**Features:** An added fenestration makes these sponges ideal for dressing around a tracheostomy tube or any type of tube on the body.

The woven sponges are ideal for when natural breathability is desired.

Usage: For dressing around a tracheostomy tube or any type of tube on the body.**Precautions:** Sterile. Latex-free.**Size:** 10 x 10cm**Standard Box Quantity:** Each

DYND73043

Respiratory Supplies - Tracheotomy

Tracheostomy, Brush**Manufacturer:** Medline **Manufacturer Code:** DYND73043**Features:** Tracheostomy Sterile Brushes come in a set of 2**Usage:** For cleaning trach tubes**Precautions:** Sterile**Size:** 5in x 1/2in**Standard Box Quantity:** Pack/2

Trach Ties, Foam, Small 23cmx3cm - Neck 18-23cm**Manufacturer:** Posey **Manufacturer Code:** 8197S**Features:** Soft Foam Padding With Hook & Loop "double Track" Closure. Helps Prevent Extubation.**Usage:** A Convenient, Comfortable Collar to Secure Tracheostomy Tubes. May Be Used With Most Tracheostomy Tubes.**Precautions:** Do not attach too tightly. Foam should not be compressed when applied.**Size:** Small**Standard Box Quantity:** Each

BD300450

Safety Needles & Syringes

Sharps Container, L (Only 1 Per Client)**Manufacturer:** BD **Manufacturer Code:** 300450**Features:** 3 Litre sharps collector.**Usage:** Specially designed for safe sharps disposal.**Precautions:** Does not contain Latex.**Size:** 3.0L**Standard Box Quantity:** Each

BD300460

Safety Needles & Syringes

Sharps Container 1.4L**Manufacturer:** BD **Manufacturer Code:** BD300460**Features:** 1.4 Litre Tray sharps collector.**Usage:** Specially designed for safe sharps disposal.**Precautions:** Product & product labelling does not contain Latex.**Size:** 1.4L**Standard Box Quantity:** 1

BD300466**Safety Needles & Syringes****Sharps Container 3.1L****Manufacturer:** BD **Manufacturer Code:** BD300466**Features:** 3.1 Litre One-piece sharps collector. Uniquely designed with point-of-use convenience and numerous built-in features. Yellow with vertical entry.**Usage:** Specially designed for safe sharps disposal.**Precautions:** Product & product labelling does not contain Latex.**Size:** 3.1L**Standard Box Quantity:** 1**BD301604****Safety Needles & Syringes****Syringe, 10mL, Slip Tip****Manufacturer:** BD **Manufacturer Code:** 301604**Features:** 10 mL BD™ slip-tip disposable syringe. Has 1/5 mL graduation.**Usage:** General Use Syringes.**Precautions:** Sterile, single use**Size:** 10cc**Standard Box Quantity:** Each**BD302830****Safety Needles & Syringes****Luer Lock Tip Syringe 20cc****Manufacturer:** BD **Manufacturer Code:** BD302830**Features:** 20 mL BD™ syringe with BD Luer-Lok™ tip. Graduation - 1 mL, 3/4 oz in 1/8 oz.**Usage:** BD™ syringes feature a clear barrel with bold scale markings, tapered plunger rod for ease of aspiration, positive plunger rod stop, and an added BD Luer-Lok™ thread for increased secure connection.**Precautions:** Sterile, single use**Size:** 20cc**Standard Box Quantity:** 1

BD302831

Safety Needles & Syringes

Slip Tip Syringe 20cc

Manufacturer: BD **Manufacturer Code:** BD302831

Features: 20 mL BD™ slip-tip disposable syringe with slip tip. Graduation 1 mL, 3/4 oz in 1/8 oz.

Usage: General Use Syringes.

Precautions: Sterile, single use

Size: 20cc

Standard Box Quantity: 1



BD302832

Safety Needles & Syringes

Syringe, 30mL, Luer Lock

Manufacturer: BD **Manufacturer Code:** 302832

Features: 30 mL BD™ syringe with BD Luer-Lok™ tip. Clear barrel with bold scale markings, tapered plunger rod for ease of aspiration, positive plunger rod stop, and an added BD Luer-Lok™ thread for increased secure connection. Packaging is clearly labeled latex free.

Usage: General Use Syringes.

Precautions: Sterile, single use

Size: 30mL

Standard Box Quantity: Each



BD302833

Safety Needles & Syringes

Slip Tip Syringe 30cc

Manufacturer: BD **Manufacturer Code:** BD302833

Features: 30 mL BD™ slip-tip disposable syringe with slip tip. Graduation 1 mL, 3/4 oz in 1/8 oz.

Usage: General Use Syringes.

Precautions: Sterile, single use

Size: 30cc

Standard Box Quantity: 1



BD303134

Safety Needles & Syringes

Slip Tip Syringe 10cc

Manufacturer: BD **Manufacturer Code:** BD303134

Features: 10 mL BD™ slip-tip disposable syringe. Has 1/5 mL graduation.

Usage: General Use Syringes.

Precautions: Cannot be re-used.

Size: 10cc

Standard Box Quantity: 1



BD305180

Safety Needles & Syringes

18 G Blunt Fill Needle

Manufacturer: BD **Manufacturer Code:** BD305180

Features: 18 G Blunt Fill Safety Needle.

Usage: General Use Needles.

Precautions: Standard hygienic practices are needed to use this product safely. Contaminated syringes must be disposed of according to the nature of the contaminate.

Size: 18 G

Standard Box Quantity: 1



BD305758

Safety Needles & Syringes

27 G x 1.27cm (.5 inch) Safety Needle

Manufacturer: BD **Manufacturer Code:** BD305758

Features: 27 G x 1.27cm (.5") Safety Needle

Usage: General Use Needles.

Precautions: Standard hygienic practices are needed to use this product safely. Contaminated syringes must be disposed of according to the nature of the contaminate.

Size: 27G x 1.27cm

Standard Box Quantity: 1



BD305762

Safety Needles & Syringes

Needle, Luer Lock, 23g x 1" (Safety)

Manufacturer: BD **Manufacturer Code:** 305762

Features: 23 G x 1 in. BD Eclipse™ Hypodermic needle, for use with BD Luer-Lok™ syringe. Single-handed activation. No hard surface required for activation. Bevel oriented needle allows for low-angle injections.

Usage: Hypodermic needle with pivoting shielding mechanism. BD PrecisionGlide™ needles minimize injection pain and patient discomfort. Needle-based design facilitates standardization and SKU reduction.

Precautions: Sterile, single use

Size: 23G 1 inch

Standard Box Quantity: Each



BD305763

Safety Needles & Syringes

Needle, Luer Lock, 22gx1 1/2" (Safety)

Manufacturer: BD **Manufacturer Code:** 305763

Features: 22G x 3.81cm (1.5") Safety Needle

Usage: General Use Needles.

Precautions: Sterile, single use

Size: 22G x 3.81cm

Standard Box Quantity: Each



BD305764

Safety Needles & Syringes

Safety Needle 21g x 1"

Manufacturer: BD **Manufacturer Code:** 305764

Features: 21 G x 2.54cm (1.0") Safety Needle

Usage: General Use Needles.

Precautions: Sterile, single use

Size: 21 G x 2.54cm

Standard Box Quantity: Each



BD305781

Safety Needles & Syringes

Syringe W/Needle, 3mL, 25g x 5/8" (Safety)

Manufacturer: BD **Manufacturer Code:** 305781

Features: 3mL luer-lock syringe with 25G x 5/8in safety needle

Usage: Hypodermic needle and syringe with shielding mechanism.

Precautions: Sterile, single use

Size: 3mL 25G 5/8 inch

Standard Box Quantity: Each



BD305783

Safety Needles & Syringes

Syringe W/Needle , 3mL, 22g x 1-1/2" (Safety)

Manufacturer: BD **Manufacturer Code:** 305783

Features: 3mL luer-lock syringe with 22G x 1-1/2in safety needle

Usage: Hypodermic needle and syringe with shielding mechanism.

Precautions: Sterile, single use

Size: 3mL 22G 1-1/2 inch

Standard Box Quantity: Each



BD305789

Safety Needles & Syringes

1 ml/27G combo pack syringe

Manufacturer: BD **Manufacturer Code:** BD305789

Features: 1 ml/27G 1/2in TB safety syringe

Usage: Hypodermic needle and syringe with shielding mechanism.

Precautions: Standard hygienic practices are needed to use this product safely. Contaminated syringes must be disposed of according to the nature of the contaminate.

Size: 1ml/27G 1/2 inch

Standard Box Quantity: 1



BD305932**Safety Needles & Syringes****Insulin Syringe 1/2cc 29G x 1/2 inch****Manufacturer:** BD **Manufacturer Code:** BD305932**Features:** 1/2mL insulin safety syringe with 29 G x 1/2 in needle.**Usage:** Hypodermic needle and syringe with shielding mechanism.**Precautions:** Standard hygienic practices are needed to use this product safely. Contaminated syringes must be disposed of according to the nature of the contaminate.**Size:** 29G x 1/2 inch 1/2cc**Standard Box Quantity:** 1**BD305945****Safety Needles & Syringes****TB Syringe 1 cc 27G x 1.27cm (.5 inch)****Manufacturer:** BD **Manufacturer Code:** BD305945**Features:** 1 mL BD™ Tuberculin syringe with 27 G x 1/2 in. BD PrecisionGlide™ detachable needle, slip tip, regular bevel, regular wall.**Usage:** The BD™ Tuberculin syringe with BD PrecisionGlide™ Needle helps improve dosage accuracy and virtually eliminates needle pop-off. Plus, bold, clear scale markings make easy and accurate dosage measurements. The thinner and smaller syringe size is easier to control and is less threatening to patients.**Precautions:** Cannot be re-used.**Size:** 1cc 27Gx1.27cm**Standard Box Quantity:** 1**BD309604****Safety Needles & Syringes****Luer Lock Tip Syringe 10 or 12cc****Manufacturer:** BD **Manufacturer Code:** BD309604**Features:** 10 mL BD Luer-Lok™ tip syringe, disposable. Has 1/5 mL graduation. BD™ syringes feature a clear barrel with bold scale markings, tapered plunger rod for ease of aspiration, positive plunger rod stop, and an added BD Luer-Lok™ thread for increased secure connection. Packaging is clearly labeled latex free.**Usage:** General Use Syringes.**Precautions:** Cannot be re-used.**Size:** 10 or 12cc**Standard Box Quantity:** 1

BD309620

Safety Needles & Syringes

Catheter Tip Syringe 60 cc

Manufacturer: BD **Manufacturer Code:** BD309620

Features: 60 mL BD™ syringe catheter tip 2 oz.

Usage: General Use Syringes.

Precautions: Standard hygienic practices are needed to use this product safely. Contaminated syringes must be disposed of according to the nature of the contaminate.

Size: 60cc

Standard Box Quantity: 1



BD309646

Safety Needles & Syringes

Luer Lock Tip Syringe 5 or 6cc

Manufacturer: BD **Manufacturer Code:** BD309646

Features: 5 mL BD Luer-Lok™ tip syringe, disposable. Has 1/5 mL graduation. BD™ syringes feature a clear barrel with bold scale markings, tapered plunger rod for ease of aspiration, positive plunger rod stop, and an added BD Luer-Lok™ thread for increased secure connection. Packaging is clearly labeled latex free.

Usage: General Use Syringes.

Precautions: Cannot be re-used.

Size: 5 or 6cc

Standard Box Quantity: 1



BD309653

Safety Needles & Syringes

Luer Lock Tip Syringe 60cc

Manufacturer: BD **Manufacturer Code:** BD309653

Features: 60 mL BD™ syringe with BD Luer-Lok™ tip. Graduation - 1 mL, 2 oz. in 1/4 oz. BD™ syringes feature a clear barrel with bold scale markings, tapered plunger rod for ease of aspiration, positive plunger rod stop, and an added BD Luer-Lok™ thread for increased secure connection. Packaging is clearly labeled latex free.

Usage: General Use Syringes.

Precautions: Cannot be re-used.

Size: 60cc

Standard Box Quantity: 1



BD309654

Safety Needles & Syringes

Syringe, 60mL, Slip Tip

Manufacturer: BD **Manufacturer Code:** 309654

Features: 60 mL syringe with slip tip.

Usage: General Use Syringes.

Precautions: Sterile, single use

Size: 60mL

Standard Box Quantity: Each



BD309656

Safety Needles & Syringes

Slip Tip Syringe 3cc

Manufacturer: BD **Manufacturer Code:** BD309656

Features: 3 mL BD™ slip-tip disposable syringe. Has 1/10 mL graduation.

Usage: BD™ syringes feature a clear barrel with bold scale markings, tapered plunger rod for ease of aspiration, and a positive plunger rod stop.

Precautions: Sterile, single use

Size: 3cc

Standard Box Quantity: 1



BD309657

Safety Needles & Syringes

Luer Lock Tip Syringe - 3cc

Manufacturer: BD **Manufacturer Code:** BD309657

Features: 3 mL BD Luer-Lok™ disposable syringe with BD Luer-Lok™ tip. Has 1/10 mL graduation.

Usage: BD™ syringes feature a clear barrel with bold scale markings, tapered plunger rod for ease of aspiration, positive plunger rod stop, and an added BD Luer-Lok™ thread for increased secure connection.

Precautions: Sterile, single use

Size: 3cc

Standard Box Quantity: 1



BD309659

Safety Needles & Syringes

Slip Tip Syringe 1cc

Manufacturer: BD **Manufacturer Code:** BD309659

Features: 1mL BD™ slip-tip disposable syringe. Has 1/10 mL graduation.

Usage: BD™ syringes feature a clear barrel with bold scale markings, tapered plunger rod for ease of aspiration, and a positive plunger rod stop.

Precautions: Sterile, single use

Size: 1cc

Standard Box Quantity: 1



BD320440

Safety Needles & Syringes

Insulin Syringe Ultra-Fine 11, 31G Short Needle 5/16"

Manufacturer: BD **Manufacturer Code:** BD320440

Features: 0.3mL ultra fine insulin syringe with 31 G x 8mm needle.

Usage: Hypodermic needle and syringe with shielding mechanism.

Precautions: Standard hygienic practices are needed to use this product safely. Contaminated syringes must be disposed of according to the nature of the contaminate.

Size: 31G x 8mm 3/10cc

Standard Box Quantity:



CHC402115

Safety Needles & Syringes

21 G x 3.8cm (1.5 inch) Safety Needle

Manufacturer: Smiths Medical **Manufacturer Code:** S402115

Features: 21 G x 3.8cm (1.5") Safety Needle

Usage: General Use Needles.

Precautions: Standard hygienic practices are needed to use this product safely. Contaminated syringes must be disposed of according to the nature of the contaminate.

Size: 21 G x 3.8cm

Standard Box Quantity: 1



CHC402558

Safety Needles & Syringes

25 G x 1.58cm (.675 inch) Safety Needle

Manufacturer: Smiths Medical **Manufacturer Code:** S402558

Features: 25 G x 1.58cm (5/8in) Safety Needle

Usage: General Use Needles.

Precautions: Standard hygienic practices are needed to use this product safely. Contaminated syringes must be disposed of according to the nature of the contaminate.

Size: 25 G x 1.58cm

Standard Box Quantity: 1



CHC4027125

Safety Needles & Syringes

27 G x 3.17cm (1.25 inch) Safety Needle

Manufacturer: Smiths Medical **Manufacturer Code:** S4027125

Features: 27 G x 3.17cm Safety Needle.

Usage: General Use Needles.

Precautions: Standard hygienic practices are needed to use this product safely. Contaminated syringes must be disposed of according to the nature of the contaminate.

Size: 27 Gx3.17cm

Standard Box Quantity: 1



CHC4319

Safety Needles & Syringes

Insulin Syringe 1 cc 27G x 1.27cm (.5 inch)

Manufacturer: Smiths Medical **Manufacturer Code:** 4319

Features: Insulin Syringe 1 cc 27G x 1.27cm (.5")

Usage: Hypodermic needle and syringe with shielding mechanism.

Precautions: Standard hygienic practices are needed to use this product safely. Contaminated syringes must be disposed of according to the nature of the contaminate.

Size: 1cc 27G x 1.27cm

Standard Box Quantity: 1



CHC44285

Safety Needles & Syringes

Insulin Syringe 100U (.5ml) 28G x 1/2 inch

Manufacturer: Smiths Medical **Manufacturer Code:** 4428-5

Features: 100U insulin safety syringe with 28 G x 1/2 in needle.

Usage: Hypodermic needle and syringe with shielding mechanism.

Precautions: Standard hygienic practices are needed to use this product safely. Contaminated syringes must be disposed of according to the nature of the contaminate.

Size: 28G x 1/2 inch 100U

Standard Box Quantity: 1



COV1180100777

Safety Needles & Syringes

Luer Lock Tuberculin Syringe 1ml

Manufacturer: Covidien **Manufacturer Code:** 1180100777

Features: 1 mL Luer-Lok disposable Tuberculin syringe. Has 1/100 mL graduation.

Usage: General Use Syringes.

Precautions: Cannot be re-used.

Size: 1ml

Standard Box Quantity: 1



COV8881850210

Safety Needles & Syringes

22 G x 2.54cm (1 inch) Safety Needle

Manufacturer: Covidien **Manufacturer Code:** 8881850210

Features: 22 G x 2.54cm Safety Needle.

Usage: General Use Needles.

Precautions: Standard hygienic practices are needed to use this product safely. Contaminated syringes must be disposed of according to the nature of the contaminate.

Size: 22 Gx2.54cm

Standard Box Quantity: 1



COV8881882712

Safety Needles & Syringes

Syringe W/Needle 1mL Tb, 27g x 1/2" (Safety)

Manufacturer: Covidien **Manufacturer Code:** 8881-882712

Features: 1 ml/27G 1/2in TB safety syringe

Usage: Hypodermic needle and syringe with shielding mechanism.

Precautions: Sterile, single use

Size: 27G x1/2 inch x 1cc

Standard Box Quantity: Each



SL201104

Safety Needles & Syringes

19G Flexible Wound Irrigation Tip F/Luer Lock

Manufacturer: Med-Rx **Manufacturer Code:** 20-1104

Features: Fits Luer Lock syringes.

Usage: Indwelling catheter.

Precautions:

Size: 19G

Standard Box Quantity: 1



BD306575

Solutions

Normal Saline 10cc Prefilled Syringe

Manufacturer: BD **Manufacturer Code:** BD306575

Features: 10 mL BD PosiFlush™ S.P. saline syringe

Usage: Pre-filled flush syringes specifically designed to reduce blood reflux that may be introduced by the plunger rod of standard syringes.

Precautions: Preservative free

Size: 10ml

Standard Box Quantity: 1



Sterile Water 10cc Vial**Manufacturer:** Hospira **Manufacturer Code:** 4887-10**Features:** Single Dose Plastic Fliptop Vial**Usage:****Precautions:** Sterile .**Size:** 10mL**Standard Box Quantity:** Each**Betadine 250ml (500ml)****Manufacturer:** LernaPharm **Manufacturer Code:** LP109-08**Features:** Aqueous solution of 10% povidone-iodine. It is a fast-acting, broad-spectrum anti-septic that kills gram-positive and gram-negative bacteria (including antibiotic resistant organisms), as well as most fungi/yeasts, viruses and protozoa.**Usage:** It is indicated for degerming skin, wounds and mucous membranes, and for preoperative skin preparation of patients. Helps reduce bacteria that potentially can cause skin infection.**Precautions:** May cause burning, redness, or a rash. Notify your doctor if you develop any of these symptoms. If you notice other effects not listed above, contact your doctor or pharmacist.**Size:** 500ml bottle**Standard Box Quantity:** 1**Normal Saline 15ml pink snap vial****Manufacturer:** CHS **Manufacturer Code:** 116**Features:** Normal saline solution dispenses by the drop or quickly for procedural efficiency in tracheal cleansing. Saline 0.9%, 15 ml Unit Dose.**Usage:** Used as a cleansing agent.**Precautions:** Single use only.**Size:** 15ml**Standard Box Quantity:** 1

Normal Saline 4oz bottle for flushing 118ml**Manufacturer:** Derma Science **Manufacturer Code:** DUP48501**Features:** Dual-top bottle with twist-open nozzle for controlled pressure dispensing.**Usage:** Precision spray nozzle provides irrigation pressure of <15 PSI. Spray nozzle can be easily removed to convert to a pour bottle.**Precautions:** Single use product. Do not re-use.**Size:** 118 ml bottle**Standard Box Quantity:** 1**Steri Stat 2% Chlorhexidine gluconate skin cleanser 450ml pump****Manufacturer:** Germiphene **Manufacturer Code:** SS-XM-P-2**Features:** 2% Chlorhexidine gluconate antiseptic skin cleanser in an easy to use 450ml pump dispenser.**Usage:** For general hand washing.**Precautions:** External use only.**Size:** 450ml pump bottle**Standard Box Quantity:** 1**Hibidil I:2000 Solution****Manufacturer:** Glebe Pharmacy **Manufacturer Code:** 450ML 2% HIBIDIL**Features:** HIBIDIL® 1:2000 Chlorhexidine Gluconate solution.**Usage:** Antibacterial agent active against a wide range of organisms. Recommended for swabbing wounds and burns and in obstetrics.**Precautions:** Not for patients who have previously shown a hypersensitivity reaction to chlorhexidine. For external use only. Keep out of the eyes and avoid contact with brain, meninges, or middle ear. Do not inject. Do not use in body cavities. If chlorhexidine solutions come into contact with eyes, wash out promptly and thoroughly with water.**Size:** 450ml**Standard Box Quantity:** 1

Sterile Water Bag For Irrigation (1000ml)**Manufacturer:** Baxter **Manufacturer Code:** JB7114**Features:** Baxter 1000ml sterile water for irrigation.**Usage:** For urethral irrigation given via an indwelling catheter.**Precautions:** Sterile. Not for injection.**Size:** 1000mL**Standard Box Quantity:** 1**Normal Saline, 3L bag****Manufacturer:** Baxter **Manufacturer Code:** JB7127**Features:** Normal Saline 0.9% NaCl.**Usage:** For Continuous Bladder Irrigation 3000 ml.**Precautions:****Size:** 3000ml**Standard Box Quantity:** 1**Sterile Water 500ml (btl)****Manufacturer:** Baxter **Manufacturer Code:** BAXJF7623P**Features:** Plastic pour bottle of sterile water.**Usage:** For irrigation and device care.**Precautions:** Do not reseal after opening.**Size:** 500 ml**Standard Box Quantity:** 1

JF7633

Solutions

Normal Saline 500 ml (btl)

Manufacturer: Baxter **Manufacturer Code:** BAXJF7633P

Features: Plastic pour bottle of sterile 0.9% saline.

Usage: For irrigation and device care.

Precautions: Do not reseal after opening.

Size: 500 ml

Standard Box Quantity: 1



JF7634

Solutions

Saline, Normal, 1000mL

Manufacturer: Baxter **Manufacturer Code:** JF7634

Features: Plastic pour bottle of sterile 0.9% saline.

Usage: For irrigation and device care.

Precautions: Do not reseal after opening.

Size: 1000mL

Standard Box Quantity: Each



SHP500

Solutions

Hydrogen Peroxide 3% 500ml (btl)

Manufacturer: Misc **Manufacturer Code:** SHP500


Features: Hydrogen Peroxide 3% 500 ml.

Usage: Used as a cleaning agent.

Precautions:

Size: 500 ml

Standard Box Quantity: 1



Chlorhexidine Gluconate Solution 500mL (Et Only)**Manufacturer:** Omega **Manufacturer Code:** L0000009**Features:** Chlorhexidine Gluconate 2% w/v aqueous with Isopropyl Alcohol 4% v/v.**Usage:** Antibacterial for the skin. Treatment of superficial skin infections.**Precautions:****Size:** 450mL**Standard Box Quantity:** Each

3M10106

Swabs and Wipes

Swabs - Chlorhexidine (Packet)**Manufacturer:** Solu Med **Manufacturer Code:** 101.06**Features:** Individually sealed wipes for cleaning wounds when water is not readily available.**Usage:** For wound cleansing and pre-dressing preparation. Ensures greater adhesion for dressings.**Precautions:** Should redness or other signs of irritation appear, discontinue use.**Size:** Small 0.65mL**Standard Box Quantity:** Each

3M10207

Swabs and Wipes

Chlorhexidine Swab Stick, Non-Alcohol**Manufacturer:** Solu Med **Manufacturer Code:** 102.07**Features:** Individually sealed swapsticks for cleaning wounds when water is not readily available.**Usage:** For wound cleansing and pre-dressing preparation. Ensures greater adhesion for dressings.**Precautions:** Should redness or other signs of irritation appear, discontinue use.**Size:** 1.6mL 4.75in**Standard Box Quantity:** Each

Cavilon No Sting Barrier Wipes**Manufacturer:** 3M **Manufacturer Code:** 3344E


Features: A fast-drying, non-sticky, alcohol-free liquid barrier film that forms a breathable, transparent coating on the skin. Protects damaged or intact skin from body fluids, adhesive trauma, friction, and incontinence. The film is hypoallergenic, non-cytotoxic and will not sting even when applied to damaged or denuded skin.

Usage: A hypoallergenic, non-cytotoxic liquid film which will not sting even when applied to damaged or denuded skin.

Precautions: If Cavilon No Sting Barrier Film is applied to an area with skin folds or other skin-to-skin contact, ensure that these areas of skin are separated to allow coating to dry thoroughly before allowing skin to return to normal position.

Size: 1 ml**Standard Box Quantity:** 1**Clinic Only - Specimen Swab Wire, Micropur, 5.5"- Pkg Of 100****Manufacturer:** AMD-Ritmed **Manufacturer Code:** 1001D

Features: Aluminum shaft specimen swab.

Usage:**Precautions:****Size:** 5.5 inch, pack of 100**Standard Box Quantity:** Pack/100**Swab Container, C&S****Manufacturer:** COPAN **Manufacturer Code:** 408C

Features: M40 Transystem Amies Agar Gel without charcoal – single plastic swab - red cap.

Usage: Throat, vaginal, wound and skin swabs.

Precautions:**Size:** 6 inch**Standard Box Quantity:** Each

Betadine Swabsticks**Manufacturer:** PDI **Manufacturer Code:** S42050**Features:** Povidone iodine swabstick**Usage:****Precautions:****Size:** 4 inch**Standard Box Quantity:** Each**Alcohol Swabs****Manufacturer:** Covidien **Manufacturer Code:** 63001**Features:** Used to wipe needles, thermometers, syringes and smaller body parts like finger tips, nails etc.**Usage:** Used to wipe needles, thermometers, syringes and smaller body parts like finger tips, nails etc.**Precautions:** Should redness or other signs of irritation appear, discontinue use.**Size:** Medium**Standard Box Quantity:** 1**Skin Prep Wipe Packet Pkg****Manufacturer:** Covidien **Manufacturer Code:** 6560**Features:** Forms protective barrier between skin and adhesive products. Helps reduce the chance of excoriation. Suitable for use in both hospital and long-term care. One wipe per packet helps reduce waste by using only what is needed.**Usage:** Helps protect skin from feces, urine, chemical irritants and gastrointestinal fluids.**Precautions:** Non-sterile.**Size:** Large**Standard Box Quantity:** Each

Lemon Glycerine Swabs (3/pkg)**Manufacturer:** PDI **Manufacturer Code:** S16225**Features:** Provides oral lubrication and stimulates saliva. Lint free swabs are impregnated with Lemon Glycerin solution.**Usage:** Effective for patients with parched lips, on oxygen, after general anesthesia and drug therapy.**Precautions:****Size:****Standard Box Quantity:** PKG/3**Skin Prep Wipe****Manufacturer:** Medline **Manufacturer Code:** MSC1500**Features:** This is a liquid film-forming wipe that forms a protective film to help reduce friction during removal of tapes and films.**Usage:** Wipes are fast and easy to use, only one coat is required, so less product is used per application. Removes easily using skin cleanser or soap and water, so unsanitary residue isn't left in or around the wound area.**Precautions:** For external use only. Keep out of the reach of children. Avoid contact with eyes. In the case of accidental contact, flush eyes well with water. Do not apply directly to open wounds or in deep puncture wounds. Should redness or other signs of irritation appear, discontinue use.**Size:****Standard Box Quantity:** 1**Povidone-iodone Swab Sticks****Manufacturer:** Dynarex **Manufacturer Code:** S42325**Features:** Aqueous solution of 10% povidone-iodine. It is a fast-acting, broad-spectrum antiseptic that kills gram-positive and gram-negative bacteria (including antibiotic resistant organisms), as well as most fungi/yeasts, viruses and protozoa. Use swab to coat area with iodine solution.**Usage:** It is indicated for degerming skin, wounds and mucous membranes, and for preoperative skin preparation of patients. Helps reduce bacteria that potentially can cause skin infection.**Precautions:** May cause burning, redness, or a rash. Notify your doctor if you develop any of these symptoms. If you notice other effects not listed above, contact your doctor or pharmacist.**Size:****Standard Box Quantity:** PKG/3

3M1527/1

Tapes

Transpore Non-Allergic 2.5cm roll

Manufacturer: 3M **Manufacturer Code:** 1527-1

Features: A latex-free, hypoallergenic, transparent and perforated plastic tape that offers strong adhesion. Transpore tape is easy to tear into very thin strips, and is easy to use even with gloves.

Usage: Securing tubing and devices (e.g., catheters, I.V. tubing). Anchoring bulky dressings.

Precautions: Apply tape without tension. Assure that skin is clean and dry before applying tape. Protect skin with 3M™ Cavilon™ No Sting Barrier Film allowing it to dry before applying tape. If possible, avoid using tackifiers.

Size: 2.5cm

Standard Box Quantity: 1



3M15272

Tapes

Transpore Non-Allergic 5cm roll

Manufacturer: 3M **Manufacturer Code:** 1527-2

Features: A latex-free, hypoallergenic, transparent and perforated plastic tape that offers strong adhesion. Transpore tape is easy to tear into very thin strips, and is easy to use even with gloves.

Usage: Securing tubing and devices (e.g., catheters, I.V. tubing). Anchoring bulky dressings.

Precautions: Apply tape without tension. Assure that skin is clean and dry before applying tape. Protect skin with 3M™ Cavilon™ No Sting Barrier Film allowing it to dry before applying tape. If possible, avoid using tackifiers.

Size: 5cm

Standard Box Quantity: 1



3M1530/1

Tapes

Micropore (Paper) 2.5cm roll

Manufacturer: 3M **Manufacturer Code:** 1530-1

Features: A latex-free, hypoallergenic paper tape that is gentle to the skin yet adheres well and leaves minimal adhesive residue upon removal. General purpose, breathable surgical tape.

Usage: Securing small to medium dressings. Securing lightweight tubing. Securing ostomy pouches. Taping fragile, at-risk skin when repeated taping is needed.

Precautions: Apply tape without tension to clean, dry skin. Protect at-risk skin with 3M™ Cavilon™ No Sting Barrier Film allowing it to dry before applying tape. In general, do not encircle a limb completely with tape. Tape should extend at least one-half inch beyond the dressing. Tape should not be pulled or stretched when applied.

Size: 2.5cm

Standard Box Quantity: 1



3M15302

Tapes

Micropore (Paper) 5cm roll

Manufacturer: 3M **Manufacturer Code:** 1530-2

Features: A latex-free, hypoallergenic paper tape that is gentle to the skin yet adheres well and leaves minimal adhesive residue upon removal. General purpose, breathable surgical tape.

Usage: Securing small to medium dressings. Securing lightweight tubing. Securing ostomy pouches. Taping fragile, at-risk skin when repeated taping is needed.

Precautions: Apply tape without tension to clean, dry skin. Protect at-risk skin with 3M™ Cavilon™ No Sting Barrier Film allowing it to dry before applying tape. In general, do not encircle a limb completely with tape. Tape should extend at least one-half inch beyond the dressing. Tape should not be pulled or stretched when applied.

Size: 5cm

Standard Box Quantity: 1



3M1541

Tapes

Steri-strips 6-7mm

Manufacturer: 3M **Manufacturer Code:** R1541

Features: 3M Steri-Strips are flexible skin closure strips made of a soft, non-woven elastic backing coated with a pressure-sensitive, hypoallergenic adhesive.

Usage: Sterile skin closure strips that are made of a porous, non-woven backing coated with a pressure-sensitive, hypoallergenic adhesive and reinforced with polyester filaments for added strength.

Precautions: Care should be taken to avoid skin damage by repeated applications on patients with thin or fragile skin.

Size: 6mm X 75mm

Standard Box Quantity: Pk/3



3M1547

Tapes

Steri-strips 1.25cmx 10cm Adhesive Skin Closures (6pkg)

Manufacturer: 3M **Manufacturer Code:** R1547

Features: 3M Steri-Strips are flexible skin closure strips made of a soft, non-woven elastic backing coated with a pressure-sensitive, hypoallergenic adhesive.

Usage: Sterile skin closure strips that are made of a porous, non-woven backing coated with a pressure-sensitive, hypoallergenic adhesive and reinforced with polyester filaments for added strength.

Precautions: Care should be taken to avoid skin damage by repeated applications on patients with thin or fragile skin.

Size: 1.25cmx 10cm

Standard Box Quantity: Pk/6



Steri-strips 2.5cmx12.5cm Adhesive Skin Closures (4pkg)**Manufacturer:** 3M **Manufacturer Code:** R1548**Features:** 3M Steri-Strips are flexible skin closure strips made of a soft, non-woven elastic backing coated with a pressure-sensitive, hypoallergenic adhesive.**Usage:** Sterile skin closure strips that are made of a porous, non-woven backing coated with a pressure-sensitive, hypoallergenic adhesive and reinforced with polyester filaments for added strength.**Precautions:** Care should be taken to avoid skin damage by repeated applications on patients with thin or fragile skin.**Size:** 2.5cm x 12.5cm**Standard Box Quantity:** Pk/4**3M Soft Cloth Medical Tape 5cm x 10m****Manufacturer:** 3M **Manufacturer Code:** 2762**Features:** Soft Cloth Tape with Liner features an S-slit liner with a printed grid. It can easily be cut to a desired size and is easy to apply. Good adhesion and excellent cross and diagonal stretch to accommodate movement or swelling.**Usage:** Breathable and gentle on the skin, supporting a broad range of applications including the securement of tubing, devices, and dressings. Gentle enough for at-risk skin, strong enough for critical securement, and stretches to accommodate movement or swelling.**Precautions:** Hypoallergenic and latex-free**Size:** 5cm x 10m**Standard Box Quantity:** 1**3M Soft Cloth Medical Tape 10cm x 10m****Manufacturer:** 3M **Manufacturer Code:** 2764**Features:** Soft Cloth Tape with Liner features an S-slit liner with a printed grid. It can easily be cut to a desired size and is easy to apply. Good adhesion and excellent cross and diagonal stretch to accommodate movement or swelling.**Usage:** Breathable and gentle on the skin, supporting a broad range of applications including the securement of tubing, devices, and dressings. Gentle enough for at-risk skin, strong enough for critical securement, and stretches to accommodate movement or swelling.**Precautions:** Hypoallergenic and latex-free**Size:** 10cm x 10m**Standard Box Quantity:** 1

3M2862

Tapes

Tape, Medipore, Soft Cloth Roll 5cm x 9.14m

Manufacturer: 3M **Manufacturer Code:** 2862

Features: Medipore tape is a non-woven polyester, soft, highly conformable tape. It has excellent adhesion yet is gentle to skin, reducing the risk of skin irritation. Medipore tape is comfortable, porous, breathable, and hypoallergenic. It has moderate stretch to accommodate swelling and movement. The tape is perforated so that no scissors are required. It is easy to handle, with or without gloves, without curling edges, or sticking to itself.

Usage: Medipore tape is particularly useful for securing wound dressings, tubing and catheters, and repeated applications where preservation of "at risk" skin is critical.

Precautions: Hypoallergenic and latex free.

Size: 2 inch

Standard Box Quantity: Each



3M2864

Tapes

Tape, Medipore, Soft Cloth Roll 10cmx9.14m

Manufacturer: 3M **Manufacturer Code:** 2864

Features: Medipore tape is a non-woven polyester, soft, highly conformable tape. It has excellent adhesion yet is gentle to skin, reducing the risk of skin irritation. Medipore tape is comfortable, porous, breathable, and hypoallergenic. It has moderate stretch to accommodate swelling and movement. The tape is perforated so that no scissors are required. It is easy to handle, with or without gloves, without curling edges, or sticking to itself.

Usage: Medipore tape is particularly useful for securing wound dressings, tubing and catheters, and repeated applications where preservation of "at risk" skin is critical.

Precautions: Hypoallergenic and latex free.

Size: 4 inch

Standard Box Quantity: Each



BSN7235907

Tapes

Leukoplast Sleek LF 5cm x 3M waterproof pink tape, latex free

Manufacturer: Smith & Nephew **Manufacturer Code:** 7235907

Features: Plastic Waterproof Adhesive Tape. Often referred to as sleek or pink tape. Leukoplast® Sleek is easy to unroll, easy to apply and can stay on for hours and days. This tape resists sweat, grease, oil and water.


Usage: Indicated for any situation where a reliable and conformable waterproof tape is required (for holding tubes or tracheotomy).

Precautions: In general, do not encircle a limb completely with tape. Tape should extend at least one-half inch beyond the dressing. Tape should not be pulled or stretched when applied.

Size: 5cm x 3m

Standard Box Quantity: Each



Leukoplast Sleek 5cm x 3M waterproof pink tape**Manufacturer:** BSN **Manufacturer Code:** BSN7237906**Features:** Plastic Waterproof Adhesive Tape. Often referred to as sleek or pink tape. Leukoplast® Sleek is easy to unroll, easy to apply and can stay on for hours and days. This tape resists sweat, grease, oil and water.**Usage:** Indicated for any situation where a reliable and conformable waterproof tape is required (for holding tubes or tracheotomy).**Precautions:** In general, do not encircle a limb completely with tape. Tape should extend at least one-half inch beyond the dressing. Tape should not be pulled or stretched when applied.**Size:** 5cm x 3m**Standard Box Quantity:** 1**Adhesive Montgomery Straps 18cm x 28cm****Manufacturer:** Medline **Manufacturer Code:** NON5129**Features:** Montgomery Straps are hypoallergenic adhesive straps used to facilitate frequent dressing changes without having to remove and reapply tape. Breathable, strong, and comfortable, they help prevent patient skin trauma associated with frequent taping.**Usage:** Used to facilitate frequent dressing changes without having to remove and reapply tape. Sheet sizes can be cut to any length.**Precautions:** Do not tie a Montgomery strap too tightly or you could restrict blood flow to the wound, which will slow healing and possibly cause further damage. Latex-free.**Size:** 18cm x 28cm**Standard Box Quantity:** Pair**Fabric Tie Tapes 130cm****Manufacturer:** Misc **Manufacturer Code:** 48200-010**Features:** This product is the twill tape used in conjunction with Montgomery straps.**Usage:** Use these twill tape ties to secure Montgomery straps for frequent dressing changes.**Precautions:** Do not tie a Montgomery strap too tightly or you could restrict blood flow to the wound, which will slow healing and possibly cause further damage.**Size:** 130cm**Standard Box Quantity:** 1

DERTCC21100

Total Contact Casting (TCC)

TCC-EZ Boot Regular

Manufacturer: Derma Science **Manufacturer Code:** TCC2-1100

Features: TCC-EZ® lightweight, customized boots simplify the process and provide additional support and stability.

Usage: For use with TCC-EZ® Total Contact Cast System.

Precautions: Patients who wear wide shoes may need to be evaluated for a Large or Extra Large Boot regardless of foot length.

Size: Male Shoe: 6-10/Female Shoe: 5-11.5

Standard Box Quantity: 1



DERTCC21114

Total Contact Casting (TCC)

TCC-EZ Boot Large

Manufacturer: Derma Science **Manufacturer Code:** TCC2-1114

Features: TCC-EZ® lightweight, customized boots simplify the process and provide additional support and stability.

Usage: For use with TCC-EZ® Total Contact Cast System.

Precautions: Patients who wear wide shoes may need to be evaluated for a Large or Extra Large Boot regardless of foot length.

Size: Male Shoe: 10.5-13W/Fem. Shoe: 11.5+ W

Standard Box Quantity: 1



DERTCC21116

Total Contact Casting (TCC)

TCC - Charcot Boot

Manufacturer: Derma Science **Manufacturer Code:** TCC21116

Features: TCC-EZ® offers Large and Extra Large Charcot Boots with innovative cushion insert specially designed to accommodate patients who have developed rocker bottom deformities.

Usage: For use with TCC-EZ® Total Contact Cast System.

Precautions: Specially designed and customizable to make dealing with complicated wounds easier.

Size: Large

Standard Box Quantity: 1



DERTCC21124

Total Contact Casting (TCC)

TCC-EZ Boot XL

Manufacturer: Derma Science **Manufacturer Code:** TCC2-1124

Features: TCC-EZ® lightweight, customized boots simplify the process and provide additional support and stability.

Usage: For use with TCC-EZ® Total Contact Cast System.

Precautions: Patients who wear wide shoes may need to be evaluated for a Large or Extra Large Boot regardless of foot length.

Size: Male Shoe: 13.5-18 XW/Fem. Shoe: 11.5+ XW **Standard Box Quantity:** 1



DERTCC21126

Total Contact Casting (TCC)

TCC - Charcot Boot

Manufacturer: Derma Science **Manufacturer Code:** TCC21126

Features: TCC-EZ® offers Large and Extra Large Charcot Boots with innovative cushion insert specially designed to accommodate patients who have developed rocker bottom deformities.

Usage: For use with TCC-EZ® Total Contact Cast System.

Precautions: Specially designed and customizable to make dealing with complicated wounds easier.

Size: Xlarge **Standard Box Quantity:** 1



DERTCC23002

Total Contact Casting (TCC)

TCC EZ Casting System Kit

Manufacturer: Derma Science **Manufacturer Code:** TCC2-3002

Features: The TCC-EZ® Total Contact Cast System provides the healthcare provider with a primary dressing and the necessary components to apply a total contact cast. The TCC-EZ® system includes instructional materials for application and removal.

Usage: Indicated for the treatment of: non-infected neuropathic foot ulcers without involvement of deeper structures, post-operative care, Charcot Neuroarthropathy, Pre-ulcerative conditions. Intended to be used in patients that are not infected and have adequate blood supply to the involved foot, in conjunction with standard wound care regimens.

Precautions: Inappropriate use or removal of the total contact cast could result in serious injury. A 4-prong cane or walker may be needed. **SPECIAL HANDLING:** store kit contents at room temperature. Keep upright at all times

Size: Calf Size: 11-17 inches **Standard Box Quantity:** 1



TCC-EZ Casting System Kit

Manufacturer: Derma Science **Manufacturer Code:** TCC2-4002

Features: The TCC-EZ® Total Contact Cast System provides the healthcare provider with a primary dressing and the necessary components to apply a total contact cast. The TCC-EZ® system includes instructional materials for application and removal.

Usage: Indicated for the treatment of: non-infected neuropathic foot ulcers without involvement of deeper structures, post-operative care, Charcot Neuroarthropathy, Pre-ulcerative conditions. Intended to be used in patients that are not infected and have adequate blood supply to the involved foot, in conjunction with standard wound care regimens.

Precautions: Inappropriate use or removal of the total contact cast could result in serious injury. A 4-prong cane or walker may be needed. **SPECIAL HANDLING:** store kit contents at room temperature. Keep upright at all times

Size: Calf Size: 15-24 inches

Standard Box Quantity: 1



TCC - Even-Up Shoe Balancer

Manufacturer: Derma Science **Manufacturer Code:** TCCEVENUPL

Features: Evenup allows patients to walk more evenly and comfortably. Adjustable straps for any enclosed flat soled shoe and fits right or left foot.

Usage: Used on opposite foot from TCC-EZ® cast for patients who have hip or back issues. Use with laced oxfords, orthopedic shoes or sneakers.

Precautions: Not a replacement for TCC-EZ® boot. Fasten securely to shoe and remove before operating a motor vehicle. Cannot be used with sandals. Balance may be compromised. Walk slowly and avoid uneven or soft ground to prevent falls. Use of a four prong cane or walker recommended.

Size: Large

Standard Box Quantity: 1



TCC - Even-Up Shoe Balancer

Manufacturer: Derma Science **Manufacturer Code:** TCC2EVENUPM

Features: Evenup allows patients to walk more evenly and comfortably. Adjustable straps for any enclosed flat soled shoe and fits right or left foot.

Usage: Used on opposite foot from TCC-EZ® cast for patients who have hip or back issues. Use with laced oxfords, orthopedic shoes or sneakers.

Precautions: Not a replacement for TCC-EZ® boot. Fasten securely to shoe and remove before operating a motor vehicle. Cannot be used with sandals. Balance may be compromised. Walk slowly and avoid uneven or soft ground to prevent falls. Use of a four prong cane or walker recommended.

Size: Medium

Standard Box Quantity: 1



TCC - Even-Up Shoe Balancer

Manufacturer: Derma Science **Manufacturer Code:** TCC2EVENUPS

Features: Evenup allows patients to walk more evenly and comfortably. Adjustable straps for any enclosed flat soled shoe and fits right or left foot.

Usage: Used on opposite foot from TCC-EZ® cast for patients who have hip or back issues. Use with laced oxfords, orthopedic shoes or sneakers.

Precautions: Not a replacement for TCC-EZ® boot. Fasten securely to shoe and remove before operating a motor vehicle. Cannot be used with sandals. Balance may be compromised. Walk slowly and avoid uneven or soft ground to prevent falls. Use of a four prong cane or walker recommended.

Size: Small

Standard Box Quantity: 1



3M16004

Transparent Films

Tegaderm Roll 10cm x 10M

Manufacturer: 3M **Manufacturer Code:** 16004

Features: This non-sterile roll version of 3M(tm) Tegaderm(tm) Transparent Film Dressing provides all the benefits of surgical tape while providing a waterproof, bacterial and viral barrier. Can be easily cut to desired shape or size. Conformable flexes with skin for greater patient comfort.

Usage: Can be used for protection of at-risk skin from friction and body fluids, as a waterproof fixation of primary dressings, and for fixation of medical devices or infusion tubing.

Precautions: Leave on as long as possible until the dressing becomes compromised (i.e. edge lift, rolled or soiled) to prevent disruption of the healing process. Dressing generally changed when fluid is "excessive," dressing is loose or routinely at 7 days.

Size: 10cm x 10M

Standard Box Quantity: Each



3M1624

Transparent Films

Tegaderm Film 6x7cm

Manufacturer: 3M **Manufacturer Code:** 1624WB

Features: Transparent film dressings designed for protecting skin and wound sites.

Usage: Superficial partial thickness burns, Stage I or II pressure ulcers, Protective Eye coverings, IV sites. Soak with Saline if difficult to remove.

Precautions: It is recommended that hemostasis be achieved prior to applying this dressing.

Size: 6cm x 7cm

Standard Box Quantity: Each



3M1626

Transparent Films

Tegaderm 10x12cm

Manufacturer: 3M **Manufacturer Code:** 1626WB

Features: Transparent film dressings designed for protecting skin and wound sites.

Usage: Superficial partial thickness burns, Stage I or II pressure ulcers, Protective Eye coverings, IV sites. Soak with Saline if difficult to remove.

Precautions: It is recommended that hemostasis be achieved prior to applying this dressing.

Size: 10cm x 12cm

Standard Box Quantity: Each



3M1628

Transparent Films

Tegaderm Flat Film 15 x 20cm

Manufacturer: 3M **Manufacturer Code:** 1628

Features: A sterile, waterproof, thin, transparent film dressing. Designed to be applied over I.V. devices, clean skin or wounds, providing a breathable bacterial and viral barrier to outside contaminants.

Usage: The hydrophilic nature of Tegaderm™ HP (holding power) Film adhesive on the sacral shaped dressing makes it more adherent in moist conditions. Tegaderm™ Film, with the "frame" delivery system, allows maximum versatility as it can be cut to any shape and size. Dressings conform to the body and flex with skin for greater patient comfort. Transparent dressings allow easy monitoring of wounds without the need to remove the dressing.

Precautions: Remember to peel the dressing back, rather than pulling it up from the skin - To aid in lifting a dressing edge, secure a piece of surgical tape to one corner and rub firmly. Use the tape as a tab to help slowly peel back the dressing.

Size: 15x20cm

Standard Box Quantity: 1



3M1629

Transparent Films

Tegaderm Flat Film 20 x 30cm

Manufacturer: 3M **Manufacturer Code:** 1629

Features: A sterile, waterproof, thin, transparent film dressing. Designed to be applied over I.V. devices, clean skin or wounds, providing a breathable bacterial and viral barrier to outside contaminants.

Usage: The hydrophilic nature of Tegaderm™ HP (holding power) Film adhesive on the sacral shaped dressing makes it more adherent in moist conditions. Tegaderm™ Film, with the "frame" delivery system, allows maximum versatility as it can be cut to any shape and size. Dressings conform to the body and flex with skin for greater patient comfort. Transparent dressings allow easy monitoring of wounds without the need to remove the dressing.

Precautions: Remember to peel the dressing back, rather than pulling it up from the skin - To aid in lifting a dressing edge, secure a piece of surgical tape to one corner and rub firmly. Use the tape as a tab to help slowly peel back the dressing.

Size: 20x30cm

Standard Box Quantity: 1



3M90800

Transparent Films

Tegaderm Clear Acrylic Sm Oval 7.6 x 9.5cm

Manufacturer: 3M **Manufacturer Code:** 90800

Features: Clear Acrylic Dressing allows wound exudate to move through perforations in the bottom layer which is coated with a moist-skin adhesive. Exudate is absorbed into the patented, clear, acrylic polymer pad. Moisture vapor is released through the top, breathable, waterproof film layer which also provides a barrier to outside contaminants. Easy to center on the wound and apply. Gentle removal from the skin. Promotes autolytic debridement while allowing visualization.

Usage: Enables clinicians to monitor pressure ulcers, skin tears and other moderately exuding wounds without changing the dressing.

Precautions: Low potential for skin maceration. Thin profile, conformable. Low friction surface. Washable. Hypoallergenic and latex-free.

Size: 7.6 x 9.5cm

Standard Box Quantity: 1



3M90801

Transparent Films

Tegaderm Clear Acrylic Med Oval 11.1 x 12.7cm

Manufacturer: 3M **Manufacturer Code:** 90801

Features: Clear Acrylic Dressing allows wound exudate to move through perforations in the bottom layer which is coated with a moist-skin adhesive. Exudate is absorbed into the patented, clear, acrylic polymer pad. Moisture vapor is released through the top, breathable, waterproof film layer which also provides a barrier to outside contaminants. Easy to center on the wound and apply. Gentle removal from the skin. Promotes autolytic debridement while allowing visualization.

Usage: Enables clinicians to monitor pressure ulcers, skin tears and other moderately exuding wounds without changing the dressing.

Precautions: Low potential for skin maceration. Thin profile, conformable. Low friction surface. Washable. Hypoallergenic and latex-free.

Size: 11.1 x 12.7cm

Standard Box Quantity: 1



3M90803

Transparent Films

Tegaderm Absorbent - 14.2 cm x 15.8 cm

Manufacturer: 3M **Manufacturer Code:** 90803

Features: Absorbent Clear Acrylic Dressing has advanced 3M technology where wound exudate moves through perforations in the bottom layer of Tegaderm™ Transparent Film, which is coated with a moist-skin adhesive. Through the process of diffusion, wound exudate is absorbed into the patented, clear, acrylic polymer pad.

Usage: Ideal for Pressure ulcers Skin tears, Abrasions and superficial wound Donor sites, Superficial partial thickness burns, Secondary (cover) dressing over wound fillers, Clean approximated surgical incisions.

Precautions: Discontinue use if the following occur. Redness, Fever, increased pain, Bleeding, Swelling, Unusual odour, cloudy or foul drainage.

Size: Large (3.4" X 4.25")

Standard Box Quantity: Each



Opsite Flexigrid 10cm x 12cm**Manufacturer:** Smith & Nephew **Manufacturer Code:** 4630**Features:** OpSite Flexigrid is a transparent, adhesive film, with a unique wound mapping grid. The film is moisture vapour permeable, conformable and extensible.**Usage:** Indicated for superficial pressure sores, minor burns, scalds abrasions and lacerations, leg ulcers in the final stage of healing, pressure sore prophylaxis and skin protection around stoma and under leg bags, to relieve pain in painful diabetic neuropathy, to provide catheter fixation, and to protect the skin from friction and external contamination.**Precautions:** Care should be taken to avoid skin damage by repeated applications on patients with thin or fragile skin.**Size:** 10x12cm**Standard Box Quantity:** 1**Opsite Flexifix Roll 5cm x 10m****Manufacturer:** Smith & Nephew **Manufacturer Code:** 66000040**Features:** The film is waterproof and a barrier to bacteria. Its high moisture vapour permeability allows the skin to breathe. The acrylic adhesive minimises the risk of skin damage on removal, especially after long periods of wear. Allows easy application of the film to awkward areas of the body and over dressings and tubes. Highly conformable and extensible to increase patient comfort.**Usage:** It can be used for :Retention of low adherent primary dressings Fixation of tubing. Treatment of painful peripheral neuropathy. Skin protection under leg bags, stoma devices etc. Reduction of shearing forces on unbroken skin, e.g. in Pressure Sore Prophylaxis.**Precautions:** Care should be taken to avoid skin damage by repeated applications on patients with thin or fragile skin.**Size:** 5cmx10m**Standard Box Quantity:** 1**Opsite Flexifix Roll 10cm x 10m****Manufacturer:** Smith & Nephew **Manufacturer Code:** 66000041**Features:** The film is waterproof and a barrier to bacteria. Its high moisture vapour permeability allows the skin to breathe. The acrylic adhesive minimises the risk of skin damage on removal, especially after long periods of wear. Allows easy application of the film to awkward areas of the body and over dressings and tubes. Highly conformable and extensible to increase patient comfort.**Usage:** It can be used for :Retention of low adherent primary dressings Fixation of tubing. Treatment of painful peripheral neuropathy. Skin protection under leg bags, stoma devices etc. Reduction of shearing forces on unbroken skin, e.g. in Pressure Sore Prophylaxis.**Precautions:** Care should be taken to avoid skin damage by repeated applications on patients with thin or fragile skin.**Size:** 10cmx10m**Standard Box Quantity:** 1

Dressing Tray**Manufacturer:** Med-Rx **Manufacturer Code:** 85-4009**Features:** 1-three compartment tray, 3-plastic transfer forceps, 6 absorbent rayon balls, 4-4x4 gauze sponges, 1-waste bag w/tie, 1-towel/underpad, 1-CSR wrap.**Usage:** Used for dressing changes.**Precautions:** Open package using aseptic technique.**Size:****Standard Box Quantity:** 1**Tracheotomy Cleaning Tray****Manufacturer:** Med-Rx **Manufacturer Code:** 85-4040**Features:** 1-three compartment tray, 1-plastic transfer forceps, 2-vinyl gloves, 5-cotton tipped applicator, 1-tracheostomy brush, 1-twill tape, 2-4x4 gauze sponges, 1-towel/underpad, 1-CSR wrap.**Usage:** Trach care tray.**Precautions:** Open package using aseptic technique.**Size:****Standard Box Quantity:** 1**Staple Removal Set****Manufacturer:** Medline **Manufacturer Code:** DYNJ04058Z**Features:** Disposable Skin Staple Remover. Sterile, Metal.**Usage:** Used to remove staples.**Precautions:** Open package using aseptic technique.**Size:****Standard Box Quantity:** 1

Suture Removal Set**Manufacturer:** Medline **Manufacturer Code:** MDS707555**Features:** Components are packaged in sequence of use, allowing the clinician to progress efficiently through procedures. Suture Removal Tray With Metal Littauer Scissors Contains: Scissors (Wire Littauer), Forceps (Plastic), Gauze (3" x 3").**Usage:** Used to remove sutures.**Precautions:** Open package using aseptic technique.**Size:****Standard Box Quantity:** 1

MLNET01

Tubular Bandages

Stretch Mesh - Finger 1.2cm #1 (per mtr)**Manufacturer:** Medline **Manufacturer Code:** NONNET01**Features:** Rayon/polyester knitted elastic net tube holds dressings firmly in place and reduces dressing time. Soft, open weave does not constrict. Easily cut to customize length. Accommodates virtually any part of the body from finger to large torso.**Usage:** For securing wound dressings. An essential product for use with fragile, geriatric, pediatric, and at-risk skin, or those allergic to tape.**Precautions:** Latex-free and neoprene free.**Size:** Finger 1.2**Standard Box Quantity:** METRE

MLNET03

Tubular Bandages

Bandage, Surgilast #3 - 1 Ft Length**Manufacturer:** Medline **Manufacturer Code:** NONNET03**Features:** Tubular bandage provides an extremely reliable method of securing wound dressings, particularly for fingers and hands. There are sizes for all parts of the body.**Usage:** For securing wound dressings on fingers and hands.**Precautions:** Non sterile.**Size:** 1ft (stretched) size #3**Standard Box Quantity:** Ft

MLNET04

Tubular Bandages

Stretch Mesh-Leg/Arm/Hand 2.5cm #4(per mtr)

Manufacturer: Medline **Manufacturer Code:** NONNET04

Features: Rayon/polyester knitted elastic net tube holds dressings firmly in place and reduces dressing time. Soft, open weave does not constrict. Easily cut to customize length. Accommodates virtually any part of the body from finger to large torso.

Usage: For securing wound dressings. An essential product for use with fragile, geriatric, pediatric, and at-risk skin, or those allergic to tape.

Precautions: Latex-free and neoprene free.

Size: Leg/Arm/Hand 2.5

Standard Box Quantity: METRE



MLNET06

Tubular Bandages

Bandage, Surgilast #6 - 1ft Length

Manufacturer: Medline **Manufacturer Code:** NONNET06

Features: Tubular bandage provides an extremely reliable method of securing wound dressings, particularly for fingers and hands. There are sizes for all parts of the body.

Usage: For securing wound dressings on fingers and hands.

Precautions: Non sterile.

Size: 1ft (stretched) size #6

Standard Box Quantity: Ft



MLNET10

Tubular Bandages

Stretch Mesh-Chest/Abdomen 6cm #10 (per mtr)

Manufacturer: Medline **Manufacturer Code:** NONNET10

Features: Rayon/polyester knitted elastic net tube holds dressings firmly in place and reduces dressing time. Soft, open weave does not constrict. Easily cut to customize length. Accommodates virtually any part of the body from finger to large torso.

Usage: For securing wound dressings. An essential product for use with fragile, geriatric, pediatric, and at-risk skin, or those allergic to tape.

Precautions: Latex-free and neoprene free.

Size: Chest/Abdomen 6

Standard Box Quantity: METRE



MM707LF

Tubular Bandages

Surgilast #6 (per mtr)

Manufacturer: Glenwood Laboratories **Manufacturer Code:** GL707LF

Features: Elasticated mesh for securement of dressings. A tubular elastic stretch net. Easy to customize. Allows maximum air-flow. Available in a wide range of sizes.

Usage: Designed to serve as a secondary dressing, applying gentle pressure to keep bandages or other medical devices securely in place without adhesive tape. Ideal for securing burn or post-op dressings, hot or cold packs, IV tubes, incontinence pads and decubitus ulcer dressings. Large size for head shoulder and or thigh.

Precautions: Latex free.

Size: 6

Standard Box Quantity: METRE



MM709LF

Tubular Bandages

Surgilast #8 (per mtr)

Manufacturer: Glenwood Laboratories **Manufacturer Code:** GL709LF

Features: Elasticated mesh for securement of dressings. A tubular elastic stretch net. Easy to customize. Allows maximum air-flow. Available in a wide range of sizes.

Usage: Designed to serve as a secondary dressing, applying gentle pressure to keep bandages or other medical devices securely in place without adhesive tape. Ideal for securing burn or post-op dressings, hot or cold packs, IV tubes, incontinence pads and decubitus ulcer dressings. Large size for head shoulder and or thigh.

Precautions: Latex free.

Size: 8

Standard Box Quantity: METRE



MOL1434

Tubular Bandages

Tubigrip E 24-28cm stretched fits large ankles, medium knees, and small thighs

Manufacturer: Molnlycke **Manufacturer Code:** 1434

Features: Tubigrip provides tissue support in the treatment of strains and sprains, soft tissue injuries, joint effusions, general oedema, post- burn scarring and ribcage injuries and is also used for pressure dressings and arm fixation.

Usage: It is easy to apply and reapply. Delivers consistent support regardless of the variation in application technique, and it is not bulky, allowing for freedom of movement and comfort.

Precautions: Product contains natural rubber latex.

Size: E

Standard Box Quantity: METRE



MOL1436

Tubular Bandages

Tubigrip B 13-16cm stretched fits sm hands & arms

Manufacturer: Molnlyke **Manufacturer Code:** 1436

Features: Tubigrip provides tissue support in the treatment of strains and sprains, soft tissue injuries, joint effusions, general oedema, post- burn scarring and ribcage injuries and is also used for pressure dressings and arm fixation.

Usage: It is easy to apply and reapply. Delivers consistent support regardless of the variation in application technique, and it is not bulky, allowing for freedom of movement and comfort.

Precautions: Product contains natural rubber latex.

Size: B

Standard Box Quantity: METRE



MOL1437

Tubular Bandages

Tubigrip D 20-24cm stretched fits large arms, med ankles, sm knees

Manufacturer: Molnlyke **Manufacturer Code:** 1437

Features: Tubigrip provides tissue support in the treatment of strains and sprains, soft tissue injuries, joint effusions, general oedema, post- burn scarring and ribcage injuries and is also used for pressure dressings and arm fixation.

Usage: It is easy to apply and reapply. Delivers consistent support regardless of the variation in application technique, and it is not bulky, allowing for freedom of movement and comfort.

Precautions: Product contains natural rubber latex.

Size: D

Standard Box Quantity: METRE



MOL1438

Tubular Bandages

Tubigrip F 28-36cm stretched fits large knees and med thighs

Manufacturer: Molnlyke **Manufacturer Code:** 1438

Features: Tubigrip provides tissue support in the treatment of strains and sprains, soft tissue injuries, joint effusions, general oedema, post- burn scarring and ribcage injuries and is also used for pressure dressings and arm fixation.

Usage: It is easy to apply and reapply. Delivers consistent support regardless of the variation in application technique, and it is not bulky, allowing for freedom of movement and comfort.

Precautions: Product contains natural rubber latex.

Size: F

Standard Box Quantity: METRE



MOL1439

Tubular Bandages

Tubigrip G 10m/Box

Manufacturer: Molnlyke **Manufacturer Code:** 1439

Features: Tubigrip provides tissue support in the treatment of strains and sprains, soft tissue injuries, joint effusions, general oedema, post- burn scarring and ribcage injuries and is also used for pressure dressings and arm fixation.

Usage: It is easy to apply and reapply. Delivers consistent support regardless of the variation in application technique, and it is not bulky, allowing for freedom of movement and comfort.

Precautions: Product contains natural rubber latex.

Size: G 36-48cm

Standard Box Quantity: Box/10M



MOL1440

Tubular Bandages

Tubigrip J 46-60cm stretched fits small trunks

Manufacturer: Molnlyke **Manufacturer Code:** 1440

Features: Tubigrip provides tissue support in the treatment of strains and sprains, soft tissue injuries, joint effusions, general oedema, post- burn scarring and ribcage injuries and is also used for pressure dressings and arm fixation.

Usage: It is easy to apply and reapply. Delivers consistent support regardless of the variation in application technique, and it is not bulky, allowing for freedom of movement and comfort.

Precautions: Product contains natural rubber latex.

Size: J

Standard Box Quantity: METRE



MOL2438

Tubular Bandages

Tubifast (65cm) Blue

Manufacturer: Molnlyke **Manufacturer Code:** 2438

Features: Unique technology holds dressings securely, without constriction or compression and without the need for tapes, pins or ties. Manufactured from rayon with very fine elastane threads knitted into the fabric to provide light elasticity. It holds dressings securely, without constriction or compression. Its light elasticity allows patients complete freedom of movement. The light, open weave fabric allows air to circulate over the skin, helping to avoid maceration.

Usage: Ideal for dressing fixation and skin covering for any part of the body. Can also be used for patch wrapping and as an under-cast stockinette.

Precautions: Latex Free.

Size: Large limbs 7.5cm wide

Standard Box Quantity: Metre



MOL2440

Tubular Bandages

Tubifast (65cm) Yellow

Manufacturer: Molnlyke **Manufacturer Code:** 2440

Features: Unique technology holds dressings securely, without constriction or compression and without the need for tapes, pins or ties. Manufactured from rayon with very fine elastane threads knitted into the fabric to provide light elasticity. It holds dressings securely, without constriction or compression. Its light elasticity allows patients complete freedom of movement. The light, open weave fabric allows air to circulate over the skin, helping to avoid maceration.

Usage: Ideal for dressing fixation and skin covering for any part of the body. Can also be used for patch wrapping and as an under-cast stockinette.

Precautions: Latex Free.

Size: Heads, children trunks 10.75cm wide

Standard Box Quantity: Metre



MOL2444

Tubular Bandages

Tubifast (65cm) Purple

Manufacturer: Molnlyke **Manufacturer Code:** 2444

Features: Unique technology holds dressings securely, without constriction or compression and without the need for tapes, pins or ties. Manufactured from rayon with very fine elastane threads knitted into the fabric to provide light elasticity. It holds dressings securely, without constriction or compression. Its light elasticity allows patients complete freedom of movement. The light, open weave fabric allows air to circulate over the skin, helping to avoid maceration.

Usage: Ideal for dressing fixation and skin covering for any part of the body. Can also be used for patch wrapping and as an under-cast stockinette.

Precautions: Latex Free.

Size: Adult trunk 20cm wide

Standard Box Quantity: Metre



SM133104LF

Tubular Bandages

Surgilast #4 (per mtr)

Manufacturer: Glenwood Laboratories **Manufacturer Code:** 13-3104LF

Features: Elasticated mesh for securement of dressings. A tubular elastic stretch net. Easy to customize. Allows maximum air-flow. Available in a wide range of sizes.

Usage: Designed to serve as a secondary dressing, applying gentle pressure to keep bandages or other medical devices securely in place without adhesive tape. Ideal for securing burn or post-op dressings, hot or cold packs, IV tubes, incontinence pads and decubitus ulcer dressings.

Precautions: Latex free.

Size: 4


Standard Box Quantity: METRE



ZipZoc**Manufacturer:** Smith & Nephew **Manufacturer Code:** 66000747**Features:** Zipzoc is a sterile rayon stocking impregnated with an ointment containing 20% Zinc Oxide Ph. Eur. Each stocking contains a mean value of 41.5g ointment.**Usage:** Stocking presentation means it is easy to apply and remove with less mess therefore saving time, compared with zinc oxide bandages. Contains zinc oxide which is soothing and has a role in wound healing. Preservative free formulation means skin reactions are less likely for patients with sensitive skin. Zipzoc and one bandage are usually thin enough to allow patients to wear their normal shoes.**Precautions:** Arterial leg ulcers. Known hypersensitivity to zinc oxide or the excipients in the ointment.**Size:** 82cm x 14cm**Standard Box Quantity:** Each

ALMP600

Wound Miscellaneous

Probe Double Ended, Non-Sterile 5"**Manufacturer:** Almedic **Manufacturer Code:** P600**Features:** Made of steel with two blunt ends.**Usage:** A blunt-ended surgical instrument used for exploring a wound or body cavity.**Precautions:** Non-sterile.**Size:** 5 inch**Standard Box Quantity:** Each

AMG118421

Wound Miscellaneous

Tongue Depressor (1/Pk) Sterile**Manufacturer:** AMG **Manufacturer Code:** 118-421**Features:** Individually wrapped tongue depressors. Sterilized and wrapped in peel open envelope.**Usage:** Single use only.**Precautions:** watch for wood splinters.**Size:** 15 cm x 2m**Standard Box Quantity:** Each

AMG500510

Wound Miscellaneous

Scalpel, Size 10, Sterile

Manufacturer: Lance Paragon **Manufacturer Code:** 500-510

Features: Disposable Scalpel #10 With Handle.

Made of stainless steel.

Usage: The fine blade is used to make precise cuts.

Precautions: Sterile.

Size: Size #10

Standard Box Quantity: Each



COL2655

Wound Miscellaneous

Coloplast, Strip Paste, 57gm

Manufacturer: Coloplast **Manufacturer Code:** 2655

Features: Coloplast Ostomy Paste is a reliable, high quality ostomy paste with less alcohol content than competitive brands. The paste provides a secure seal yet remains soft and flexible in use.

Usage: Shape the strip around the hole in your flange prior to applying it to the skin. Run tube under hot water to soften it prior to application.

Precautions: Keep in package until ready to use to avoid drying of product.

Size: 57gm

Standard Box Quantity: 1



COL3210

Wound Miscellaneous

Coloplast Skin Barrier 10cm x 10cm

Manufacturer: Coloplast **Manufacturer Code:** 3210

Features: Skin Barrier Protective Sheets are secure and flexible for total skin protection around stomas, draining wounds and other areas subject to irritation. Skin Barrier Protective Sheets are non-pectin based.

Usage: Medicated to heal and protect broken down skin due to stool or urine.

Precautions: Latex free.

Size: 10 cm X 10 cm

Standard Box Quantity: 1




CON839262

Wound Miscellaneous

Convatec Eakin Wound Pouch 175mm x 110mm (839262)**Manufacturer:** ConvaTec **Manufacturer Code:** 839262**Features:** Eakin Wound Pouches provide a comfortable, skin friendly, cost effective and time saving procedure for managing fistulas, multiple or irregular stomas and high-output wounds.**Usage:** Large pouch with remote drainage attachment for wounds up to 175 X 110mm.**Precautions:****Size:** 175mm x 110mm**Standard Box Quantity:** 1

HOL9775

Wound Miscellaneous

Hollister Wound Drainage Collector with Adhesive Tape Border Sterile 9.5cm (9775)**Manufacturer:** Hollister **Manufacturer Code:** 9775**Features:** Protects patient by containing corrosive exudate while reducing exposure to bacteria and contaminants. Keeps patient dry and comfortable by collecting exudate, eliminating odour and conforming to the wound site. Easy to use and cost effective.**Usage:** Cost effective alternatives for the management of draining wounds, fistulae or draining wound sites. All wound drainage collectors are made with odour-barrier plastic, a convenience drain and a unique reclosable cap that permits access to the wound site without removing the pouch.**Precautions:** Sterile. Stays in-situ up to 7 days.**Size:** Medium 9.5cm**Standard Box Quantity:** 1

HOL9777

Wound Miscellaneous

Hollister Wound Drainage Collector with Adhesive Tape Border Sterile 10cm x 20cm (9777)**Manufacturer:** Hollister **Manufacturer Code:** 9777**Features:** Protects patient by containing corrosive exudate while reducing exposure to bacteria and contaminants. Keeps patient dry and comfortable by collecting exudate, eliminating odour and conforming to the wound site. Easy to use and cost effective.**Usage:** Cost effective alternatives for the management of draining wounds, fistulae or draining wound sites. All wound drainage collectors are made with odour-barrier plastic, a convenience drain and a unique reclosable cap that permits access to the wound site without removing the pouch.**Precautions:** Sterile. Stays in-situ up to 7 days.**Size:** Large 10cm x 20cm**Standard Box Quantity:** 1

Dermal Curette, 4mm Disposable

Manufacturer: Integra York PA Inc. **Manufacturer Code:** 33-54

Features: Dermal Curettes feature a unique symmetrically ground beveled edge. Manufactured of surgical grade stainless steel providing a clean, precise edge every time guaranteed.

Usage:

Precautions: Individually packaged in sterile pouches.

Size: 4mm

Standard Box Quantity: Each

