

The City model makes for a comfortable and smooth ride at a price point unmatched for its capability.

Technical specs

Max load:	180 kg/396 lbs
Tires:	8" PU solid tire (Front) 12" PU solid tire (Back)
Max slope capacity:	12°
Max speed:	6 km/3.7 miles per hour
Max distance (full charge):	25 km/15.5 miles
Steering:	360° joystick controller
Motor:	2 x 250W brushless motors
Seat depth:	43 cm/17"
Seat width:	49 cm/19.25"
Seat height:	46 cm/18"
Braking distance:	1.5 m (when at top speed)
Turning radius:	85 cm/33"
Battery:	2 x Lithium 24V/6AH
Charger:	24V/2A
Time required to fully charge battery:	6-8 h
Net weight (including battery):	26.5 kg/58.5 lbs
Net weight (without battery):	23.5 kg/52 lbs
Frame material:	Aluminum alloy
Unfolded (L x W x H):	38" x 24" x 37"
Folded (L x W x H):	23" x 13" x 31"