



Continence Supplies Catalogue

2018-10-26

Category Table of Contents

Catheter Trays	Solutions
Catheter Trays 3	Solutions 43
Catheters - External	Swabs/Wipes/Creams/Moisturizers
Catheters - External 3	Swabs/Wipes/Creams/ Moisturizers 44
Catheters - Foley	Tapes
Catheters - Foley..... 6	Tapes..... 50
Catheters - Intermittent	Urinals/Bedpans
Catheters - Intermittent..... 19	Urinals/Bedpans..... 52
Catheters - Rubber	
Catheters - Rubber..... 24	
Drainage Bags/Tubing	
Drainage Bags/Tubing..... 26	
Gloves	
Gloves 34	
Pads/Briefs	
Pads/Briefs..... 37	
Pediatric Catheters	
Pediatric Catheters 40	
Safety Needles & Syringes	
Safety Needles & Syringes 41	
Securements	
Securements 42	

Tray, Catheter Add-A-Foley**Manufacturer:** Amsino **Manufacturer Code:** AS890**Features:** Contains: 1000ml Outer Basin Tray; Syringe Prefilled with 30cc Sterile Water; 2 Vinyl Gloves; Waterproof Drape; Lubricating Jelly; Fenestrated Drape; 10% Povidone Iodine Swabsticks**Usage:** For catheterization procedures.**Precautions:** Sterile. Latex free.**Size:** Tray**Standard Box Quantity:** Each**Catheter Tray****Manufacturer:** Med-Rx **Manufacturer Code:** 85-5027**Features:** Contains 1250ml graduated basin, 10cc pre-filled syringe sterile water, 2 plastic transfer forceps, prep cup, 5 absorbent rayon ball, lubricating jelly, PVP iodine, 2 non-latex gloves, waterproof underpad, fenestrated drape, specimen container with ID label, CSR wrap.**Usage:** Latex Free.**Precautions:** Sterile unless package is damaged or opened.**Size:** 1250ml basin, 10cc pre-filled syringe**Standard Box Quantity:** 1**Condom Catheter Non-Latex - Small 25mm****Manufacturer:** Coloplast **Manufacturer Code:** 5225H**Features:** Conveen Security+ self-sealing male external catheter.**Usage:** This external catheter needs to be changed every 24 to 48 hours. After removal, the penis and surrounding area should be cleansed thoroughly with soap and water. The area should then be dried completely before applying a new catheter.**Precautions:** Latex free**Size:** 25mm small**Standard Box Quantity:** 1

COL6100

Catheters - External

Condom Catheter Non-Latex 100% Silicone 23mm SM 6100

Manufacturer: Coloplast **Manufacturer Code:** 504311

Features: The Condom Catheter is a flexible sheath that fits over the penis just like a condom. The condom part is then attached to a tube that drains the urine into a urinary storage bag. These types of catheters are usually good for those individuals who suffer from male urinary incontinence due to a medical condition or advanced Alzheimer's disease.

Usage: This external catheter needs to be changed every 24 to 48 hours. After removal, the penis and surrounding area should be cleansed thoroughly with soap and water. The area should then be dried completely before applying a new catheter.

Precautions: Should be changed frequently, either every day or on alternate days. If there's any leakage of urine, it should be changed promptly as urine is an irritant.

Size: 23mm small

Standard Box Quantity: 1



COL6200

Catheters - External

Condom Catheter Non-Latex 100% Silicone 28mm MED 6200

Manufacturer: Coloplast **Manufacturer Code:** 504331

Features: The Condom Catheter is a flexible sheath that fits over the penis just like a condom. The condom part is then attached to a tube that drains the urine into a urinary storage bag. These types of catheters are usually good for those individuals who suffer from male urinary incontinence due to a medical condition or advanced Alzheimer's disease.

Usage: This external catheter needs to be changed every 24 to 48 hours. After removal, the penis and surrounding area should be cleansed thoroughly with soap and water. The area should then be dried completely before applying a new catheter.

Precautions: Should be changed frequently, either every day or on alternate days. If there's any leakage of urine, it should be changed promptly as urine is an irritant.

Size: 28mm medium

Standard Box Quantity: 1



COL6300

Catheters - External

Condom Catheter Non-Latex 100% Silicone 31mm INTER 6300

Manufacturer: Coloplast **Manufacturer Code:** 504341

Features: The Condom Catheter is a flexible sheath that fits over the penis just like a condom. The condom part is then attached to a tube that drains the urine into a urinary storage bag. These types of catheters are usually good for those individuals who suffer from male urinary incontinence due to a medical condition or advanced Alzheimer's disease.

Usage: This external catheter needs to be changed every 24 to 48 hours. After removal, the penis and surrounding area should be cleansed thoroughly with soap and water. The area should then be dried completely before applying a new catheter.

Precautions: Should be changed frequently, either every day or on alternate days. If there's any leakage of urine, it should be changed promptly as urine is an irritant.

Size: 31mm intermediate

Standard Box Quantity: 1



Condom Catheter Non-Latex 100% Silicone 35mm LG 6400**Manufacturer:** Coloplast **Manufacturer Code:** 504351

Features: The Condom Catheter is a flexible sheath that fits over the penis just like a condom. The condom part is then attached to a tube that drains the urine into a urinary storage bag. These types of catheters are usually good for those individuals who suffer from male urinary incontinence due to a medical condition or advanced Alzheimer's disease.

Usage: This external catheter needs to be changed every 24 to 48 hours. After removal, the penis and surrounding area should be cleansed thoroughly with soap and water. The area should then be dried completely before applying a new catheter.

Precautions: Should be changed frequently, either every day or on alternate days. If there's any leakage of urine, it should be changed promptly as urine is an irritant.

Size: 35mm large**Standard Box Quantity:** 1**Condom Catheter Latex w/tape - Large 35mm****Manufacturer:** Medline **Manufacturer Code:** DYND12301

Features: The Condom Catheter is a flexible sheath that fits over the penis just like a condom. The condom part is then attached to a tube that drains the urine into a urinary storage bag. These types of catheters are usually good for those individuals who suffer from male urinary incontinence due to a medical condition or advanced Alzheimer's disease.

Usage: This external catheter needs to be changed every 24 to 48 hours. After removal, the penis and surrounding area should be cleansed thoroughly with soap and water. The area should then be dried completely before applying a new catheter.

Precautions: Should be changed frequently, either every day or on alternate days. If there's any leakage of urine, it should be changed promptly as urine is an irritant. Contains latex.



Size: Large**Standard Box Quantity:** 1**Condom Catheter Latex w/tape - Medium 30mm****Manufacturer:** Medline **Manufacturer Code:** DYND12302

Features: The Condom Catheter is a flexible sheath that fits over the penis just like a condom. The condom part is then attached to a tube that drains the urine into a urinary storage bag. These types of catheters are usually good for those individuals who suffer from male urinary incontinence due to a medical condition or advanced Alzheimer's disease.

Usage: This external catheter needs to be changed every 24 to 48 hours. After removal, the penis and surrounding area should be cleansed thoroughly with soap and water. The area should then be dried completely before applying a new catheter.

Precautions: Should be changed frequently, either every day or on alternate days. If there's any leakage of urine, it should be changed promptly as urine is an irritant. Contains latex.

Size: Medium**Standard Box Quantity:** 1

Catheter, Coude Tip Foley, 16F 5-10cc -Lubricious coated**Manufacturer:** Bard **Manufacturer Code:** 0102L16**Features:** The Olive coude of this model features a slightly larger bulb at the tip to assist in negotiation of restrictions.**Usage:** Coude-tipped catheters are used to navigate the male urethra, particularly when a stricture is present.**Precautions:** Sterile.**Size:** 16FR 5CC**Standard Box Quantity:** Each**Lubricath Latex 2 way 5cc Foley 12FR****Manufacturer:** Bard **Manufacturer Code:** 0165L12**Features:** Bard Medical Bardex® Lubricath® 2-Way Foley Catheter 12Fr 5cc Balloon Capacity, Hydrogel Coating, Sterile.**Usage:** The BARDEX® LUBRICATH® latex Foley catheter provides for smoother insertion and enhanced patient comfort through the hydrogel coating's unique ability to absorb fluid, thus creating a hydrophilic "cushion" between the catheter surface and the urethra.**Precautions:** Contains latex. Do not over inflate balloon past recommended volume. Use sterile water or normal saline.**Size:** 12FR 5cc**Standard Box Quantity:** 1**Silicone Coated 2 way 30cc 24FR****Manufacturer:** Bard **Manufacturer Code:** 33424**Features:** This urinary catheter is latex, but it is coated internally and externally with a silicone elastomer coating. This enables the catheter to resist encrustation and bacterial colonization.**Usage:** A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment.**Precautions:** Do not over inflate balloon past recommended volume. Use sterile water or normal saline.**Size:** 24FR 30 cc**Standard Box Quantity:** 1

Catheter, Foley, 24 Fr 30cc-Lubricious Coated**Manufacturer:** Bard **Manufacturer Code:** 366724

Features: The exclusive hydrogel coating interlocks with the inner latex substrate to create a lubricious barrier that helps reduce the risk of infection. -Absorbs mucosal fluid to form a "cushion" between the catheter surface and delicate urethral tissue, reducing friction and irritation that can lead to infection. -Provides a surface which resists encrustation.

Usage:


Precautions: Do not over inflate balloon past recommended volume. Use sterile water to inflate balloon. Contains latex.

Size: 30cc 24FR**Standard Box Quantity:** Each**Catheter, Foley 2-way, 28F 5cc balloon-Silicone Coated****Manufacturer:** Amsino **Manufacturer Code:** AS41028

Features: The Silicone-Elastomer-Coated Foley Catheter is an indwelling catheter featuring a new and innovative technology. This catheter undergoes a manufacturing process that provides the maximum amount of silicone coating available. This urinary catheter is latex, but it is coated internally and externally with a silicone elastomer coating. This enable the catheter to resist encrustation and bacterial colonization.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.

Precautions: Do not over inflate balloon past recommended volume. Use sterile water to inflate balloon.

Size: 5cc 28FR**Standard Box Quantity:** Each**Catheter, Foley, 24F 5-10cc balloon-Silicone Coated****Manufacturer:** Covidien **Manufacturer Code:** 402724

Features: The Silicone-Elastomer-Coated Foley Catheter is a latex indwelling catheter, but it is coated with a silicone elastomer coating for superior lubricity. A ribbed balloon provides increased strength and symmetry. Bspak spring-loaded valve ensures proper inflation/deflation of Foley catheter balloon.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.

Precautions: Do not over inflate balloon past recommended volume. Use sterile water to inflate balloon. Contains latex.

Size: 5cc 24FR**Standard Box Quantity:** Each

Silicone Coated 2 way 30cc 12FR**Manufacturer:** Covidien **Manufacturer Code:** 403712

Features: The Silicone-Elastomer-Coated Foley Catheter is an indwelling catheter featuring a new and innovative technology. This urinary catheter is latex, but it is coated internally and externally with a silicone elastomer coating. This enables the catheter to resist encrustation and bacterial colonization.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.


Precautions: Do not over inflate balloon past recommended volume. Use sterile water or normal saline.

Size: 12FR 30cc**Standard Box Quantity:** 1**Catheter, 100% Silicone, 22F 5cc****Manufacturer:** Covidien **Manufacturer Code:** 8887-605221

Features: 100% Silicone Foley Catheters: Pliable catheter with a firm tip eases insertion. Features large, smooth eyes for maximum comfort and drainage and symmetrical balloon for added reliability. 2-way Foley, 5cc balloon, 20Fr.

Usage: Silicone Foley catheters are widely accepted as having better biocompatibility characteristics and lower encrustation and blockage rates.

Precautions: Do not over inflate balloon past recommended volume. Use sterile water to inflate balloon.

Size: 5cc 22FR**Standard Box Quantity:** Each**Catheter, 3-way Foley, 16F 5-10cc-Lubricious Coated****Manufacturer:** Covidien **Manufacturer Code:** 8887688169

Features: The Silicone-Elastomer-Coated Foley Catheter is latex, but it is coated internally and externally with a silicone elastomer coating. This enables the catheter to resist encrustation and bacterial colonization. Silicone-Elastomer coating provides superior lubricity. Large, smooth drain eyes for maximum drainage. Ribbed balloon increases strength and symmetry. Bspak spring-loaded valve ensures proper inflation/deflation of balloon.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.

Precautions: Do not over inflate balloon past recommended volume. Use sterile water to inflate balloon. Contains latex.

Size: 16FR 5CC**Standard Box Quantity:** Each

DYND11504

Catheters - Foley

All Silicone (100%) 2 way 5cc 20FR

Manufacturer: Medline **Manufacturer Code:** DYND11504

Features: 100% Silicone Foley Catheters: Pliable catheter with a firm tip eases insertion. Features large, smooth eyes for maximum comfort and drainage and symmetrical balloon for added reliability. 2-way Foley, 5cc balloon, 20Fr.

Usage: Silicone Foley catheters are widely accepted as having better biocompatibility characteristics and lower encrustation and blockage rates.

Precautions: Do not over inflate balloon past recommended volume. Use sterile water or normal saline.

Size: 5cc 20FR

Standard Box Quantity: 1



DYND11505

Catheters - Foley

All Silicone (100%) 2 way 5cc 22FR

Manufacturer: Medline **Manufacturer Code:** DYND11505

Features: 100% Silicone Foley Catheters: Pliable catheter with a firm tip eases insertion. Features large, smooth eyes for maximum comfort and drainage and symmetrical balloon for added reliability. 2-way Foley, 5cc balloon, 22Fr.

Usage: Silicone Foley catheters are widely accepted as having better biocompatibility characteristics and lower encrustation and blockage rates.

Precautions: Do not over inflate balloon past recommended volume. Use sterile water or normal saline.

Size: 5cc 22FR

Standard Box Quantity: 1



DYND11532

Catheters - Foley

All Silicone (100%) 2 way 30cc 16FR

Manufacturer: Medline **Manufacturer Code:** DYND11532

Features: 100% Silicone Foley Catheters: Pliable catheter with a firm tip eases insertion. Features large, smooth eyes for maximum comfort and drainage and symmetrical balloon for added reliability. 2-way Foley, 30cc balloon, 16Fr.

Usage: Silicone Foley catheters are widely accepted as having better biocompatibility characteristics and lower encrustation and blockage rates.

Precautions: Do not over inflate balloon past recommended volume. Use sterile water or normal saline.

Size: 30cc 16FR

Standard Box Quantity: 1



DYND11533

Catheters - Foley

All Silicone (100%) 2 way 30cc 18FR

Manufacturer: Medline **Manufacturer Code:** DYND11533

Features: 100% Silicone Foley Catheters: Pliable catheter with a firm tip eases insertion. Features large, smooth eyes for maximum comfort and drainage and symmetrical balloon for added reliability. 2-way Foley, 30cc balloon, 18Fr.

Usage: Silicone Foley catheters are widely accepted as having better biocompatibility characteristics and lower encrustation and blockage rates.

Precautions: Do not over inflate balloon past recommended volume. Use sterile water or normal saline.

Size: 30cc 18FR

Standard Box Quantity: 1



ML11500

Catheters - Foley

All Silicone (100%) 2 way 5cc 12FR

Manufacturer: Medline **Manufacturer Code:** DYND11500

Features: 100% Silicone Foley Catheters: Pliable catheter with a firm tip eases insertion. Features large, smooth eyes for maximum comfort and drainage and symmetrical balloon for added reliability. 2-way Foley, 5cc balloon, 12Fr.

Usage: Silicone Foley catheters are widely accepted as having better biocompatibility characteristics and lower encrustation and blockage rates.

Precautions: Do not over inflate balloon past recommended volume. Use sterile water or normal saline.

Size: 5cc 12FR

Standard Box Quantity: 1



ML11501

Catheters - Foley

All Silicone (100%) 2 way 5cc 14FR

Manufacturer: Medline **Manufacturer Code:** DYND11501

Features: 100% Silicone Foley Catheters: Pliable catheter with a firm tip eases insertion. Features large, smooth eyes for maximum comfort and drainage and symmetrical balloon for added reliability. 2-way Foley, 5cc balloon, 14Fr.

Usage: Silicone Foley catheters are widely accepted as having better biocompatibility characteristics and lower encrustation and blockage rates.

Precautions: Do not over inflate balloon past recommended volume. Use sterile water or normal saline.

Size: 5cc 14FR

Standard Box Quantity: 1



ML11502

Catheters - Foley

All Silicone (100%) 2 way 5cc 16FR

Manufacturer: Medline **Manufacturer Code:** DYND11502

Features: 100% Silicone Foley Catheters: Pliable catheter with a firm tip eases insertion. Features large, smooth eyes for maximum comfort and drainage and symmetrical balloon for added reliability. 2-way Foley, 5cc balloon, 16Fr.

Usage: Silicone Foley catheters are widely accepted as having better biocompatibility characteristics and lower encrustation and blockage rates.

Precautions: Do not over inflate balloon past recommended volume. Use sterile water or normal saline.

Size: 5cc 16FR

Standard Box Quantity: 1



ML11503

Catheters - Foley

All Silicone (100%) 2 way 5cc 18FR

Manufacturer: Medline **Manufacturer Code:** DYND11503

Features: 100% Silicone Foley Catheters: Pliable catheter with a firm tip eases insertion. Features large, smooth eyes for maximum comfort and drainage and symmetrical balloon for added reliability. 2-way Foley, 5cc balloon, 18Fr.

Usage: Silicone Foley catheters are widely accepted as having better biocompatibility characteristics and lower encrustation and blockage rates.

Precautions: Do not over inflate balloon past recommended volume. Use sterile water or normal saline.

Size: 5cc 18FR

Standard Box Quantity: 1



ML11506

Catheters - Foley

Catheter, Silicone, 24Fr 5-10cc

Manufacturer: Medline **Manufacturer Code:** DYND11506

Features: 100% Silicone Foley Catheters: Pliable catheter with a firm tip eases insertion. Features large, smooth eyes for maximum comfort and drainage and symmetrical balloon for added reliability. 2-way Foley, 5cc balloon, 18Fr.

Usage: Silicone Foley catheters are widely accepted as having better biocompatibility characteristics and lower encrustation and blockage rates.

Precautions: Do not over inflate balloon past recommended volume. Use sterile water to inflate balloon.

Size: 24Fr 5-10cc

Standard Box Quantity: Each



ML11534

Catheters - Foley

Catheter, Foley, 100% Silicone, 20F 30cc

Manufacturer: Medline **Manufacturer Code:** DYND11534

Features: Pliable catheter with a firm tip eases insertion. Features large, smooth eyes for maximum comfort and drainage and symmetrical balloon for added reliability. Clear silicone color makes it easy for visualization of clots, mucus and urine flow through the catheter. Blue radiopaque strip makes the catheter easily visible during scans.

Usage: Silicone material causes less tissue irritation and encrustation during extended periods of indwelling use.

Precautions: Latex free.

Size: 20FR 30CC

Standard Box Quantity: Each



ML11593

Catheters - Foley

Catheter, Coude Tip, 100% Silicone, 18F 5-10cc

Manufacturer: Medline **Manufacturer Code:** DYND11593

Features: 2-Way, 100% silicone coude tip Foley catheter.

Usage: For difficult insertions, coude-tip Foley catheters feature a bent tip that maneuvers around obstructions leading to the bladder. Silicone material causes less tissue irritation and encrustation during extended periods of indwelling use.

Precautions: Latex free.

Size: 18FR 10CC

Standard Box Quantity: Each



ML11752

Catheters - Foley

Silicone Coated 2 way 5cc 12FR

Manufacturer: Medline **Manufacturer Code:** DYND11752

Features: The Silicone-Elastomer-Coated Foley Catheter is an indwelling catheter featuring a new and innovative technology. This catheter undergoes a manufacturing process that provides the maximum amount of silicone coating available. This urinary catheter is latex, but it is coated internally and externally with a silicone elastomer coating. This enable the catheter to resist encrustation and bacterial colonization.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.

Precautions: Do not over inflate balloon past recommended volume. Use sterile water or normal saline.

Size: 5cc 12FR

Standard Box Quantity: 1



ML11754

Catheters - Foley

Silicone Coated 2 way 5cc 14FR

Manufacturer: Medline **Manufacturer Code:** DYND11754

Features: The Silicone-Elastomer-Coated Foley Catheter is an indwelling catheter featuring a new and innovative technology. This catheter undergoes a manufacturing process that provides the maximum amount of silicone coating available. This urinary catheter is latex, but it is coated internally and externally with a silicone elastomer coating. This enables the catheter to resist encrustation and bacterial colonization.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.

Precautions: Do not over inflate balloon past recommended volume. Use sterile water or normal saline.

Size: 5cc 14FR

Standard Box Quantity: 1



ML11756

Catheters - Foley

Silicone Coated 2 way 5cc 16FR

Manufacturer: Medline **Manufacturer Code:** DYND11756

Features: The Silicone-Elastomer-Coated Foley Catheter is an indwelling catheter featuring a new and innovative technology. This catheter undergoes a manufacturing process that provides the maximum amount of silicone coating available. This urinary catheter is latex, but it is coated internally and externally with a silicone elastomer coating. This enable the catheter to resist encrustation and bacterial colonization.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.

Precautions: Do not over inflate balloon past recommended volume. Use sterile water or normal saline.

Size: 5cc 16FR

Standard Box Quantity: 1



ML11758

Catheters - Foley

Silicone Coated 2 way 5cc 18FR

Manufacturer: Medline **Manufacturer Code:** DYND11758

Features: The Silicone-Elastomer-Coated Foley Catheter is an indwelling catheter featuring a new and innovative technology. This catheter undergoes a manufacturing process that provides the maximum amount of silicone coating available. This urinary catheter is latex, but it is coated internally and externally with a silicone elastomer coating. This enable the catheter to resist encrustation and bacterial colonization.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.

Precautions: Do not over inflate balloon past recommended volume. Use sterile water or normal saline.

Size: 5cc 18FR

Standard Box Quantity: 1



ML11760

Catheters - Foley

Silicone Coated 2 way 5cc 20FR

Manufacturer: Medline **Manufacturer Code:** DYND11760

Features: The Silicone-Elastomer-Coated Foley Catheter is an indwelling catheter featuring a new and innovative technology. This catheter undergoes a manufacturing process that provides the maximum amount of silicone coating available. This urinary catheter is latex, but it is coated internally and externally with a silicone elastomer coating. This enable the catheter to resist encrustation and bacterial colonization.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.

Precautions: Do not over inflate balloon past recommended volume. Use sterile water or normal saline.

Size: 5cc 20FR

Standard Box Quantity: 1



ML11762

Catheters - Foley

Silicone Coated 2 way 5cc 22FR

Manufacturer: Medline **Manufacturer Code:** DYND11762

Features: The Silicone-Elastomer-Coated Foley Catheter is an indwelling catheter featuring a new and innovative technology. This catheter undergoes a manufacturing process that provides the maximum amount of silicone coating available. This urinary catheter is latex, but it is coated internally and externally with a silicone elastomer coating. This enable the catheter to resist encrustation and bacterial colonization.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.

Precautions: Do not over inflate balloon past recommended volume. Use sterile water or normal saline.

Size: 5cc 22FR

Standard Box Quantity: 1



ML11774

Catheters - Foley

Silicone Coated 2 way 30cc 14FR

Manufacturer: Medline **Manufacturer Code:** DYND11774

Features: The Silicone-Elastomer-Coated Foley Catheter is an indwelling catheter featuring a new and innovative technology. This catheter undergoes a manufacturing process that provides the maximum amount of silicone coating available. This urinary catheter is latex, but it is coated internally and externally with a silicone elastomer coating. This enable the catheter to resist encrustation and bacterial colonization.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.

Precautions: Do not over inflate balloon past recommended volume. Use sterile water or normal saline.

Size: 30cc 14FR

Standard Box Quantity: 1



ML11776

Catheters - Foley

Silicone Coated 2 way 30cc 16FR

Manufacturer: Medline **Manufacturer Code:** DYND11776

Features: The Silicone-Elastomer-Coated Foley Catheter is an indwelling catheter featuring a new and innovative technology. This catheter undergoes a manufacturing process that provides the maximum amount of silicone coating available. This urinary catheter is latex, but it is coated internally and externally with a silicone elastomer coating. This enable the catheter to resist encrustation and bacterial colonization.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.

Precautions: Do not over inflate balloon past recommended volume. Use sterile water or normal saline.

Size: 30cc 16FR

Standard Box Quantity: 1



ML11778

Catheters - Foley

Silicone Coated 2 way 30cc 18FR

Manufacturer: Medline **Manufacturer Code:** DYND11778

Features: The Silicone-Elastomer-Coated Foley Catheter is an indwelling catheter featuring a new and innovative technology. This catheter undergoes a manufacturing process that provides the maximum amount of silicone coating available. This urinary catheter is latex, but it is coated internally and externally with a silicone elastomer coating. This enable the catheter to resist encrustation and bacterial colonization.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.

Precautions: Do not over inflate balloon past recommended volume. Use sterile water or normal saline.

Size: 30cc 18FR

Standard Box Quantity: 1



ML11780

Catheters - Foley

Silicone Coated 2 way 30cc 20FR

Manufacturer: Medline **Manufacturer Code:** DYND11780

Features: The Silicone-Elastomer-Coated Foley Catheter is an indwelling catheter featuring a new and innovative technology. This catheter undergoes a manufacturing process that provides the maximum amount of silicone coating available. This urinary catheter is latex, but it is coated internally and externally with a silicone elastomer coating. This enable the catheter to resist encrustation and bacterial colonization.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.

Precautions: Do not over inflate balloon past recommended volume. Use sterile water or normal saline.

Size: 30cc 20FR

Standard Box Quantity: 1



ML11782

Catheters - Foley

Silicone Coated 2 way 30cc 22FR

Manufacturer: Medline **Manufacturer Code:** DYND11782

Features: The Silicone-Elastomer-Coated Foley Catheter is an indwelling catheter featuring a new and innovative technology. This catheter undergoes a manufacturing process that provides the maximum amount of silicone coating available. This urinary catheter is latex, but it is coated internally and externally with a silicone elastomer coating. This enable the catheter to resist encrustation and bacterial colonization.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.

Precautions: Do not over inflate balloon past recommended volume. Use sterile water or normal saline.

Size: 30cc 22FR

Standard Box Quantity: 1



ML11800

Catheters - Foley

Silicone Coated 3 way 30cc 16FR

Manufacturer: Medline **Manufacturer Code:** DYND11800

Features: The Silicone-Elastomer-Coated Foley Catheter is an indwelling catheter featuring a new and innovative technology. This catheter undergoes a manufacturing process that provides the maximum amount of silicone coating available. This urinary catheter is latex, but it is coated internally and externally with a silicone elastomer coating. This enable the catheter to resist encrustation and bacterial colonization.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.

Precautions: Do not over inflate balloon past recommended volume. Use sterile water or normal saline.

Size: 30cc 16FR

Standard Box Quantity: 1



ML11801

Catheters - Foley

Silicone Coated 3 way 30cc 18FR

Manufacturer: Medline **Manufacturer Code:** DYND11801

Features: The Silicone-Elastomer-Coated Foley Catheter by Medline is an indwelling catheter featuring a new and innovative technology. This catheter undergoes a manufacturing process that provides the maximum amount of silicone coating available. This urinary catheter is latex, but it is coated internally and externally with a silicone elastomer coating. This enable the catheter to resist encrustation and bacterial colonization. Additional features include a balloon, 2-way, and sterile packaging. A variety of popular sizes are available.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.

Precautions: Do not over inflate balloon past recommended volume. Use sterile water or normal saline.

Size: 30cc

Standard Box Quantity: 1




R173830220

Catheters - Foley

Catheter, 3-way, 22 F 30cc balloon - 100% Silicone**Manufacturer:** Rusch **Manufacturer Code:** 173830220**Features:** This is a 100% Silicone, 3 way Foley with 2 drainage eyes and 1 irrigation eye. It has a 30cc balloon and is 40 cm (16") long. It is sterile for single use and is color coded.**Usage:** Silicone Foley catheters are widely accepted as having better biocompatibility characteristics and lower encrustation and blockage rates.**Precautions:** Do not over inflate balloon past recommended volume. Use sterile water to inflate balloon.**Size:** 30cc Balloon**Standard Box Quantity:** Each


R318114


Catheters - Foley

Coated Latex 2way Foley Coude Tip 14FR 5cc**Manufacturer:** Rusch **Manufacturer Code:** RUS3181-14**Features:** Rusch pure Gold PTFE Coated Latex 2 way has 1 eye, and a Coude/Tiemann Tip. It has a 5cc balloon, is 40 cm (16") long and color coded.**Usage:** Coude tip simplifies insertion.**Precautions:** Sterile for single use. PTFE coated latex.**Size:** 14FR 5cc 16in**Standard Box Quantity:** 1

R318116

Catheters - Foley

Coated Latex 2way Foley Coude Tip 16FR 5cc**Manufacturer:** Rusch **Manufacturer Code:** RUS3181-16**Features:** Rusch pure Gold PTFE Coated Latex 2 way has 1 eye, and a Coude/Tiemann Tip. It has a 5cc balloon, is 40 cm (16") long and color coded.**Usage:** Coude tip simplifies insertion.**Precautions:** Sterile for single use. PTFE coated latex.**Size:** 16FR 5cc 16in**Standard Box Quantity:** 1

R318118**Catheters - Foley****Coated Latex 2way Foley Coude Tip 18FR 5cc****Manufacturer:** Rusch **Manufacturer Code:** RUS3181-18**Features:** Rusch pure Gold PTFE Coated Latex 2 way has 1 eye, and a Coude/Tiemann Tip. It has a 5cc balloon, is 40 cm (16") long and color coded.**Usage:** Coude tip simplifies insertion.**Precautions:** Sterile for single use. PTFE coated latex.**Size:** 18FR 5cc 16in**Standard Box Quantity:** 1**RUS183430200****Catheters - Foley****Silicone Coated 3 way 30cc 20FR****Manufacturer:** Rusch **Manufacturer Code:** RUS183430200**Features:** The Silicone-Elastomer-Coated Foley Catheter is an indwelling catheter featuring a new and innovative technology. This catheter undergoes a manufacturing process that provides the maximum amount of silicone coating available. This urinary catheter is latex, but it is coated internally and externally with a silicone elastomer coating. This enable the catheter to resist encrustation and bacterial colonization.**Usage:** A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.**Precautions:** Do not over inflate balloon past recommended volume. Use sterile water or normal saline.**Size:** 30cc 20FR**Standard Box Quantity:** 1**RUS183430220****Catheters - Foley****Silicone Coated 3 way 30cc 22FR****Manufacturer:** Rusch **Manufacturer Code:** RUS183430220**Features:** The Silicone-Elastomer-Coated Foley Catheter is an indwelling catheter featuring a new and innovative technology. This catheter undergoes a manufacturing process that provides the maximum amount of silicone coating available. This urinary catheter is latex, but it is coated internally and externally with a silicone elastomer coating. This enable the catheter to resist encrustation and bacterial colonization.**Usage:** A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.**Precautions:** Do not over inflate balloon past recommended volume. Use sterile water or normal saline.**Size:** 30cc 22FR**Standard Box Quantity:** 1

Catheter, #12F Short Disposable for Intermittent Catheterization (Female)**Manufacturer:** Amici **Manufacturer Code:** 5612**Features:** Intermittent catheters feature specialized low-friction coatings and polished eyes that help improve clinical outcomes and patient comfort.**Usage:** A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.**Precautions:** Precautions should be used to keep the catheter clean and free of bacteria. Patients requiring intermittent catheterization should be well trained in the technique by a qualified health care professional.**Size:** 6"**Standard Box Quantity:** Each**Catheter, #14F Short Disposable for Intermittent Catheterization (Female)****Manufacturer:** Amici **Manufacturer Code:** 5614**Features:** Intermittent catheters feature specialized low-friction coatings and polished eyes that help improve clinical outcomes and patient comfort.**Usage:** A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.**Precautions:** Precautions should be used to keep the catheter clean and free of bacteria. Patients requiring intermittent catheterization should be well trained in the technique by a qualified health care professional.**Size:** 6"**Standard Box Quantity:** Each**Catheter, Intermittent 8F Male****Manufacturer:** Amici **Manufacturer Code:** 5908**Features:** Intermittent catheters feature specialized low-friction coatings and polished eyes that help improve clinical outcomes and patient comfort. They are available in 16" length in 8fr, 10fr, 12fr, 14fr, 16fr and 18fr sizes.**Usage:** A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.**Precautions:** Precautions should be used to keep the catheter clean and free of bacteria. Patients requiring intermittent catheterization should be well trained in the technique by a qualified health care professional.**Size:** 8FR 16 inch**Standard Box Quantity:** Each

Catheter, Intermittent, 10F Male**Manufacturer:** Amici **Manufacturer Code:** 5910

Features: Intermittent catheters feature specialized low-friction coatings and polished eyes that help improve clinical outcomes and patient comfort. They are available in 16" length in 8fr, 10fr, 12fr, 14fr, 16fr and 18fr sizes.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.

Precautions: Precautions should be used to keep the catheter clean and free of bacteria. Patients requiring intermittent catheterization should be well trained in the technique by a qualified health care professional.

Size: 10FR 16 inch**Standard Box Quantity:** Each**Catheter, #12F Long Disposable For Intermittent Catheterization Male****Manufacturer:** Amici **Manufacturer Code:** 5912

Features: Intermittent catheters feature specialized low-friction coatings and polished eyes that help improve clinical outcomes and patient comfort. They are available in 16" length in 8fr, 10fr, 12fr, 14fr, 16fr and 18fr sizes.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.

Precautions: Precautions should be used to keep the catheter clean and free of bacteria. Patients requiring intermittent catheterization should be well trained in the technique by a qualified health care professional.

Size: 12FR 16 inch**Standard Box Quantity:** Each**Catheter, #14F Long Disposable For Intermittent Catheterization Male****Manufacturer:** Amici **Manufacturer Code:** 5914

Features: Intermittent catheters feature specialized low-friction coatings and polished eyes that help improve clinical outcomes and patient comfort. They are available in 16" length in 8fr, 10fr, 12fr, 14fr, 16fr and 18fr sizes.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.

Precautions: Precautions should be used to keep the catheter clean and free of bacteria. Patients requiring intermittent catheterization should be well trained in the technique by a qualified health care professional.

Size: 14FR 16 inch**Standard Box Quantity:** Each

COL210

Catheters - Intermittent

Female 10FR 6 inch

Manufacturer: Coloplast **Manufacturer Code:** 504610

Features: The female intermittent catheter is pre-lubricated - ready to use. Easy-to-open sterile foil container. Extremely smooth, safe and comfortable.

Usage: Urinary catheter for intermittent use. The Self-Cath female catheter has been designed for the shorter female urethra.

Precautions: The product is for single use only. It contains phthalates and should not be used by either children or pregnant or nursing women without consulting a physician. Keep away from sunlight.

Size: 10FR 6 inches

Standard Box Quantity: 1



COL212

Catheters - Intermittent

Female 12FR 6 inch

Manufacturer: Coloplast **Manufacturer Code:** 504620

Features: The female intermittent catheter is pre-lubricated - ready to use. Easy-to-open sterile foil container. Extremely smooth, safe and comfortable.

Usage: Urinary catheter for intermittent use. The Self-Cath female catheter has been designed for the shorter female urethra.

Precautions: The product is for single use only. It contains phthalates and should not be used by either children or pregnant or nursing women without consulting a physician. Keep away from sunlight.

Size: 12FR 6 inches

Standard Box Quantity: 1



COL214

Catheters - Intermittent

Female 14FR 6 inch

Manufacturer: Coloplast **Manufacturer Code:** 504630

Features: The female intermittent catheter is pre-lubricated - ready to use. Easy-to-open sterile foil container. Extremely smooth, safe and comfortable.

Usage: Urinary catheter for intermittent use. The Self-Cath female catheter has been designed for the shorter female urethra.

Precautions: The product is for single use only. It contains phthalates and should not be used by either children or pregnant or nursing women without consulting a physician. Keep away from sunlight.

Size: 14FR 6 inches

Standard Box Quantity: 1



COL410

Catheters - Intermittent

Mentor Male Catheter 10FR 16 inch

Manufacturer: Coloplast **Manufacturer Code:** 504510

Features: Popular straight tip intermittent catheter. Straight tip fire polished eyelets. Siliconized surface for smooth insertion.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra.

Precautions: 100% Silicone - latex free.

Size: 10FR

Standard Box Quantity: 1



COL412

Catheters - Intermittent

Mentor Male Catheter 12FR 16 inch

Manufacturer: Coloplast **Manufacturer Code:** 504520

Features: Popular straight tip intermittent catheter. Straight tip fire polished eyelets. Siliconized surface for smooth insertion.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra.

Precautions: 100% Silicone - latex free.

Size: 12FR

Standard Box Quantity: 1



COL414

Catheters - Intermittent

Mentor Male Catheter 14FR 16 inch

Manufacturer: Coloplast **Manufacturer Code:** 504530


Features: Popular straight tip intermittent catheter. Straight tip fire polished eyelets. Siliconized surface for smooth insertion.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra.

Precautions: 100% Silicone - latex free.

Size: 14FR

Standard Box Quantity: 1



COL416

Catheters - Intermittent

Mentor Male Catheter 16FR 16 inch

Manufacturer: Coloplast **Manufacturer Code:** 504540

Features: Popular straight tip intermittent catheter. Straight tip fire polished eyelets. Siliconized surface for smooth insertion.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra.

Precautions: 100% Silicone - latex free.

Size: 16FR

Standard Box Quantity: 1



COL814

Catheters - Intermittent

Catheter, Intermittent, Coude Olive Tip, 14fr Male guide stripe

Manufacturer: Coloplast **Manufacturer Code:** 814

Features: Intermittent catheters feature specialized low-friction coatings and polished eyes that help improve clinical outcomes and patient comfort. They are available in 16" length in 8fr, 10fr, 12fr, 14fr, 16fr and 18fr sizes.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.

Precautions: Precautions should be used to keep the catheter clean and free of bacteria. Patients requiring intermittent catheterization should be well trained in the technique by a qualified health care professional.

Size: 16"

Standard Box Quantity: Each



SL605016

Catheters - Intermittent

Catheter, 16F Male Intermittent

Manufacturer: Med-Rx **Manufacturer Code:** 60-5016


Features: Intermittent catheters feature specialized low-friction coatings and polished eyes that help improve clinical outcomes and patient comfort.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra.

Precautions: Precautions should be used to keep the catheter clean and free of bacteria. Patients requiring intermittent catheterization should be well trained in the technique by a qualified health care professional.

Size: 16FR 16 inch

Standard Box Quantity: Each



SL605018

Catheters - Intermittent

Catheter, Long Intermittent, 18F Male

Manufacturer: Med-Rx **Manufacturer Code:** 605018


Features: Intermittent catheters feature specialized low-friction coatings and polished eyes that help improve clinical outcomes and patient comfort.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra.

Precautions: Precautions should be used to keep the catheter clean and free of bacteria. Patients requiring intermittent catheterization should be well trained in the technique by a qualified health care professional.

Size: 18FR 16 inch

Standard Box Quantity: Each



COV8888492017

Catheters - Rubber

All Purpose (Red Rubber) 8FR

Manufacturer: Covidien **Manufacturer Code:** 8888492017

Features: Urethral Catheter. Features two staggered eyes. Maximizes Urine Flow. Made from thermosensitive polyvinyl chloride (PVC) material.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.

Precautions: Latex Free. Radiopaque. Visible in X-Rays.

Size: 8FR

Standard Box Quantity: 1



COV8888492025

Catheters - Rubber

All Purpose (Red Rubber) 10FR

Manufacturer: Covidien **Manufacturer Code:** 8888492025

Features: Urethral Catheter. Features two staggered eyes. Maximizes Urine Flow. Made from thermosensitive polyvinyl chloride (PVC) material.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.

Precautions: Latex Free. Radiopaque. Visible in X-Rays.

Size: 10FR

Standard Box Quantity: 1



COV8888492033

Catheters - Rubber

All Purpose (Red Rubber) 12FR

Manufacturer: Covidien **Manufacturer Code:** 8888492033

Features: Urethral Catheter. Features two staggered eyes. Maximizes Urine Flow. Made from thermosensitive polyvinyl chloride (PVC) material.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.

Precautions: Latex Free. Radiopaque. Visible in X-Rays.

Size: 12FR

Standard Box Quantity: 1



COV8888492041

Catheters - Rubber

All Purpose (Red Rubber) 14FR

Manufacturer: Covidien **Manufacturer Code:** 8888492041

Features: Urethral Catheter. Features two staggered eyes. Maximizes Urine Flow. Made from thermosensitive polyvinyl chloride (PVC) material.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.

Precautions: Latex Free. Radiopaque. Visible in X-Rays.

Size: 14FR

Standard Box Quantity: 1



COV8888492058

Catheters - Rubber

All Purpose (Red Rubber) 16FR

Manufacturer: Covidien **Manufacturer Code:** 8888492058

Features: Urethral Catheter. Features two staggered eyes. Maximizes Urine Flow. Made from thermosensitive polyvinyl chloride (PVC) material.

Usage: A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.

Precautions: Latex Free. Radiopaque. Visible in X-Rays.

Size: 16FR

Standard Box Quantity: 1



All Purpose (Red Rubber) 18FR**Manufacturer:** Covidien **Manufacturer Code:** 8888492066**Features:** Urethral Catheter. Features two staggered eyes. Maximizes Urine Flow. Made from thermosensitive polyvinyl chloride (PVC) material.**Usage:** A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.**Precautions:** Latex Free. Radiopaque. Visible in X-Rays.**Size:** 18FR**Standard Box Quantity:** 1**All Purpose (Red Rubber) 20FR****Manufacturer:** Medline **Manufacturer Code:** DYND13520**Features:** Urethral Catheter. Features two staggered eyes. Maximizes Urine Flow. Made from thermosensitive polyvinyl chloride (PVC) material.**Usage:** A larger french (FR) size can become necessary when the urine is thick, bloody or contains large amounts of sediment. However, larger may cause damage to the urethra. Most common sizes are 14 or 16 fr.**Precautions:** Latex Free. Radiopaque. Visible in X-Rays.**Size:** 20FR**Standard Box Quantity:** 1**Urinary Drainage Bag, Leg, Flip Flow Valve, 32 Oz****Manufacturer:** Bard **Manufacturer Code:** 1501LF**Features:** Made from heavy duty vinyl for extended wear, odor containment and leak resistance. Each leg bag contains an anti-reflux valve. FLIP-FLO™ drainage valve opens and closes with a simple flip.**Usage:** For urinary drainage.**Precautions:** Monitor bag and empty to avoid overflow. Replace every 5-7 days or if discoloration or odour is present.**Size:** 1000mL**Standard Box Quantity:** Each

CHS101006

Drainage Bags/Tubing

3-way Stop Cock

Manufacturer: Med-Rx **Manufacturer Code:** 10-1006

Features: Three way stopcock with Luer Lock connectors

Usage:

Precautions:

Size:

Standard Box Quantity: Each



COV145514

Drainage Bags/Tubing

Latex Free Leg Bag 500ml W Straps w/o tubing

Manufacturer: Covidien **Manufacturer Code:** 145514

Features: 500cc/17oz Drainage Leg Bag With Straps (No Tubing).

Usage: A drain bag is a urine collection device that is connected to a catheter in order to drain the contents of the bladder.

Precautions: Do not use if package is opened or damaged. Latex free. Single use.

Size: 500cc

Standard Box Quantity: 1



COV145516

Drainage Bags/Tubing

Leg Bag 730/1000ml Latex Free Kendall

Manufacturer: Covidien **Manufacturer Code:** 145516

Features: Anti-Reflux valve. Longer catheter adapter for greater security with all types of catheters. Twist valve for easier emptying and less spill; no bottom cap to lose or misplace.

Usage: Includes leg straps.

Precautions: Latex-Free.

Size: 1000ml

Standard Box Quantity: 1



COV30510

Drainage Bags/Tubing

Urinary Drainage Bag - Bedside, 2000mL

Manufacturer: Covidien **Manufacturer Code:** 30510

Features: Dover™ Urine Drainage Bag with needleless sampling port, Anti-Reflux device, drain tube, poly bag, 2000 mL

Usage: A drain bag is a urine collection device that is connected to a catheter in order to drain the contents of the bladder. Fits securely on bedrails and wheelchairs.

Precautions: Monitor bag and empty to avoid overflow. Replace every 5-7 days or if discoloration or odour is present.

Size: 2000 cc

Standard Box Quantity: Each



COV3512

Drainage Bags/Tubing

Drainage Bed Bag 2000cc/Tubing/Connector

Manufacturer: Covidien **Manufacturer Code:** 3512

Features: The Kendall line of drain bags have a 2000ml capacity and have available features such as the SPLASHGUARD II™ Drain Spout, SAFEGUARD™ Needleless Sampling Port, and MONO-FLO™ Anti-Reflux device.

Usage: A drain bag is a urine collection device that is connected to a catheter in order to drain the contents of the bladder. Fits securely on bedrails and wheelchairs.

Precautions: Latex Free.

Size: 2000cc

Standard Box Quantity: 1



DS600

Drainage Bags/Tubing

Nephrostomy bag 600ml

Manufacturer: Dartin Medical **Manufacturer Code:** DS600

Features: The Dartin Select 600 is one of the most patient friendly drainage bags available. It is a convenient 600 ml capacity, has a flock backing for enhanced comfort when worn next to the skin, has a 50cm silicone tube terminating in a male luer lock to allow a secure and kink free connection to a nephrostomy or other percutaneously inserted drainage catheter, has an easy to use tap at the bottom for emptying, and has a PVC extension below the tap to allow the connection of a larger bag for night time use.

Usage: Dartin Select bag is designed for the collection of fluid during nephrostomy, biliary, and abscess drainages, with particular attention being paid to achieving maximum procedural efficiency and patient comfort.

Precautions: DEHP & Latex Free

Size: 600ml

Standard Box Quantity: 1



HOL7331

Drainage Bags/Tubing

Hollister Drain Tube Adapter

Manufacturer: Hollister **Manufacturer Code:** 7331

Features: Convenient and easy to attach to urostomy pouch. Safe to use.

Usage: Allows easy attachment of a night drainage collector to the pouch. For use with Series 14xx, 390x, 1840x, and 1841x Urostomy Pouches.

Precautions: Latex free.

Size:

Standard Box Quantity: 10



HOL9343

Drainage Bags/Tubing

Hollister, Leg Strap, Medium

Manufacturer: Hollister **Manufacturer Code:** 9343


Features: Wide, lightweight, woven fabric allows skin to breathe and resists rolling or binding. Smooth plastic connectors. Machine washable.

Usage: The Hollister fabric leg bag strap is wide to help hold the bag comfortably in place.

Precautions: Latex Free.

Size: Medium

Standard Box Quantity: 10



HOL9346

Drainage Bags/Tubing

Urostomy Connecting Tubing

Manufacturer: Hollister **Manufacturer Code:** 9346

Features: Latex-free, kink-resistant extension tubing with connector.


Usage: For use with urinary catheters and drainage bags.

Precautions: Sterile

Size: 46cm

Standard Box Quantity: 1



Hollister, Horizontal Drain/Tube Attachment**Manufacturer:** Hollister **Manufacturer Code:** 9781**Features:** Sterile, flexible support secured to a synthetic skin barrier with FlexWear backing. The device may stay in place up to 10 days.**Usage:** Used to secure, stabilize and support a variety of catheters, drains and tubes. Eliminates the need for taping or suturing tubes to skin. Protects the surrounding skin and minimizes tube movement.**Precautions:** Sterile, Latex Free.**Size:** Fits 5-40 FR**Standard Box Quantity:** 1**Urinary Drainage Bag - Leg, Medium 500mL****Manufacturer:** Medline **Manufacturer Code:** DYND12574**Features:** Leg Bag W/ Comfort Strap and Twist Valve Drainage Port: These sturdy bags are latex-free and have a sterile fluid pathway.**Usage:** Leg bag should be attached to a tube and/or catheter. To drain, open the spigot at the bottom of the leg bag when it becomes one-third to one-half full. Either empty it into the toilet or into a sterile container, if doctor wants you to measure urine volume. Clean the spigot with rubbing alcohol and a cotton ball or an alcohol pad.**Precautions:** Leg bags are Sterile and standard precautions should be used when applying them. Ensure that straps are not too tight around leg. Monitor bag and empty to avoid overflow. Replace every 5-7 days or if discoloration or odour is present.**Size:** 550 ml**Standard Box Quantity:** Each**Urinary Drainage Bag - Leg, Large 1000mL****Manufacturer:** Medline **Manufacturer Code:** 12578**Features:** Leg Bags with comfort strap and twist valve drainage port. These sturdy bags are latex-free and have a sterile fluid pathway.**Usage:** Leg bag should be attached to a tube and/or catheter. To drain, open the spigot at the bottom of the leg bag when it becomes one-third to one-half full. Either empty it into the toilet or into a sterile container, if doctor wants you to measure urine volume. Clean the spigot with rubbing alcohol and a cotton ball or an alcohol pad.**Precautions:** Monitor bag and empty to avoid overflow. Replace every 5-7 days or if discoloration or odour is present. Latex-free.**Size:** 1000 ml**Standard Box Quantity:** Each

Bedside Bag 4000cc**Manufacturer:** Medline **Manufacturer Code:** DYND15405**Features:** Anti-reflux design minimizes back flow into the drainage tube and features a clear one way mylar valve for visualization of urine flow. Includes metal clamp, needleless sample port, and slide tap for easy one-handed drainage. Reinforced hanger with built-in anti-kink drain tube guide and adjustable bed sheet clamp.**Usage:** Designed for extra capacity during continuous irrigation.**Precautions:** Latex-free.**Size:** 4000ml**Standard Box Quantity:** 1**Nephrostomy 14FR Connecting Tube****Manufacturer:** Medical Systems **Manufacturer Code:** TR-012RSC**Features:** Nephrostomy 14FR Connecting Tube**Usage:** Connecting Tube 30cm, with Blue Lip and Removable Low Pressure Stop Cock.**Precautions:** Latex free. Non DEHP.**Size:** 14FR**Standard Box Quantity:** 1**Seamless Enema Bag Set****Manufacturer:** Med-Rx **Manufacturer Code:** 50-800**Features:** Features the Touch N Seal valve atop a value oriented bag and administration tubing set. Prelubricated tip is hygienically protected by a plastic sheath. Includes a package of castile soap, plastic-backed underpad, 1500cc bag and 60" administration tube."**Usage:** Enema kits for colon cleansing.**Precautions:** Sterile technique required.**Size:** 1500cc**Standard Box Quantity:** 1

SL588018

Drainage Bags/Tubing

18FR Rectal Tube, Sterile, Disposable

Manufacturer: Med-Rx **Manufacturer Code:** 58-8018

Features: Plastic rectal tube with lube packet

Usage: Used for the administration of enemas and certain medications. Also for relief of flatulence.

Precautions: Single use.

Size: 18Fr

Standard Box Quantity: 1



SL588024

Drainage Bags/Tubing

24FR Rectal Tube-Length 50cm

Manufacturer: Med-Rx **Manufacturer Code:** 58-8024

Features: Plastic rectal tube with lube packet

Usage: Used for the administration of enemas and certain medications. Also for relief of flatulence.

Precautions: Single use.

Size: 24FR 30 cc

Standard Box Quantity: 1



SL651234

Drainage Bags/Tubing

Extension Tubing + Connector (For Leg Bag) 45cm

Manufacturer: Sherwood Medical **Manufacturer Code:** 65-1234

Features: Urinary Extension tube Latex Free 45cm.

Usage: For connecting a male external catheter or Foley catheter to a drainage bag.

Precautions: Latex-free.

Size: 45cm

Standard Box Quantity: 1



Catheter Clamp

Manufacturer: Urocare **Manufacturer Code:** 7000-00

Features: Urocare Adjustable Thumb Clamp. Two position.

Usage: May be used with clear-vinyl, latex or white rubber drainage or extension tubing as a flow control device.

Precautions:

Size: Universal

Standard Box Quantity: 1

**Cook 14FR Urological Connecting Tube**

Manufacturer: Cook **Manufacturer Code:** 50040

Features: Standard with male Luer Lock and Drainage Bag Connector; Female Luer Lock to Tuohy-Borst Adapter Attached.

Usage: Used to provide connection to a drainage bag.

Precautions:

Size: 14FR

Standard Box Quantity: 1

**Lymphedema Tubing 50cm 14Fr**

Manufacturer: Cook **Manufacturer Code:** 50041

Features: Male Luer Lock and Drainage Bag Connector; Female Luer Lock to Tuohy-Borst Adapter Attached.

Usage: Used to provide connection to a drainage bag.

Precautions:

Size: 50cm

Standard Box Quantity: 1



Cysto bladder irrigation set**Manufacturer:** Baxter **Manufacturer Code:** JC4040**Features:** Cysto/Bladder Irrigation Set includes a regulating clamp, drip chamber, catheter adapter, 0.195"ID tubing sealed in a sterile peel pouch.**Usage:** Used for bladder irrigation. Set plugs directly into 3-way catheter or catheter adapter unit when 2-way catheter is used. .**Precautions:** Sterile.**Size:****Standard Box Quantity:** 1

MAX484405P

Gloves

Gloves Sterile Small (pr)**Manufacturer:** Misc **Manufacturer Code:** 484405**Features:** Examination Gloves, Sterile, Latex Free Vinyl Small.**Usage:** For examination procedures that require latex free gloves.**Precautions:** Vinyl gloves should not be used when handling chemotherapy drugs.**Size:** Small**Standard Box Quantity:** PAIR

MAX484406

Gloves

Gloves Sterile Medium (pr)**Manufacturer:** Misc **Manufacturer Code:** 484406**Features:** Examination Gloves, Sterile, Latex Free Vinyl Medium.**Usage:** For examination procedures that require latex free gloves.**Precautions:** Vinyl gloves should not be used when handling chemotherapy drugs.**Size:** Medium**Standard Box Quantity:** PAIR

MAX484407

Gloves

Gloves Sterile Large (pr)

Manufacturer: Misc **Manufacturer Code:** 484407

Features: Examination Gloves, Sterile, Latex Free Vinyl Large.

Usage: For examination procedures that require latex free gloves.

Precautions: Vinyl gloves should not be used when handling chemotherapy drugs.

Size: Large

Standard Box Quantity: PAIR



MAX484802

Gloves

Gloves Non Sterile Powder Free Chemo - Med (pr)

Manufacturer: Misc **Manufacturer Code:** 484802

Features: Non-Sterile Powder-Free Nitrile Exam Gloves have durable thin mil nitrile material that provides maximum protection against chemicals, harmful abrasions, Type I allergic reactions to NRL and viral exposure. Comfortable and secure even during the most difficult procedures.

Usage: Glove is approved for chemotherapy treatments.

Precautions: Latex free, non-sterile, single use.

Size: Medium

Standard Box Quantity: PAIR



MAX484803

Gloves

Gloves Non Sterile Powder Free Chemo - Lg (pr)

Manufacturer: Misc **Manufacturer Code:** 484803

Features: Non-Sterile Powder-Free Nitrile Exam Gloves have durable thin mil nitrile material that provides maximum protection against chemicals, harmful abrasions, Type I allergic reactions to NRL and viral exposure. Comfortable and secure even during the most difficult procedures.

Usage: Glove is approved for chemotherapy treatments.

Precautions: Latex free, non-sterile, single use.

Size: Large

Standard Box Quantity: PAIR



MED1155B

Gloves

Gloves Non Sterile Small (pr)

Manufacturer: Misc **Manufacturer Code:** 1155B

Features: Examination Gloves, Non-Sterile, Latex Free Vinyl Small.

Usage: For examination procedures that require latex free gloves.

Precautions: Vinyl gloves should not be used when handling chemotherapy drugs.

Size: Small

Standard Box Quantity: PAIR



MED1195CP

Gloves

Gloves Non Sterile Medium (pr)

Manufacturer: Misc **Manufacturer Code:** 1195C


Features: Examination Gloves, Non-Sterile, Latex Free Vinyl Medium.

Usage: For examination procedures that require latex free gloves.

Precautions: Vinyl gloves should not be used when handling chemotherapy drugs.

Size: Medium

Standard Box Quantity: PAIR



MED1195DP

Gloves

Gloves Non Sterile Large (pr)

Manufacturer: Misc **Manufacturer Code:** 1195D


Features: Examination Gloves, Non-Sterile, Latex Free Vinyl Large.


Usage: For examination procedures that require latex free gloves.

Precautions: Vinyl gloves should not be used when handling chemotherapy drugs.

Size: Large

Standard Box Quantity: PAIR



Reusable Incontinence Pads 34x36 inch**Manufacturer:** Medline **Manufacturer Code:** MDTIU4TEFGRE**Features:** Extra Absorbent Reusable Underpads. Heavy soaker pad is a blend of rayon and polyester and absorbs more fluid.**Usage:** Designed for moderate to heavy incontinence. Best protection against moisture damage to mattresses and bedding.**Precautions:** Perineal skin is particularly susceptible to damage commonly associated with incontinence.**Size:** 34x36 inches**Standard Box Quantity:** 1**Adult Brief - Small 22 - 35 inch****Manufacturer:** Covidien **Manufacturer Code:** 60032**Features:** Wings™ adult briefs have a full fit design to provide complete coverage and maximum absorbency. The design includes an extra wide crotch for added protection, four tape tabs and a blue dryness strip for optimum skin protection.**Usage:** Super absorbent core is designed to quickly lock in fluid and neutralize unpleasant odors. Fastening system features four sturdy tape tabs to ensure a secure fit.**Precautions:** Latex free.**Size:** Small 20 - 31 inches Peach**Standard Box Quantity:** 1**Incontinence Pad 58 x 91cm****Manufacturer:** Covidien **Manufacturer Code:** 7174**Features:** Curity Disposable Underpads with non-woven facing are soft against skin. Absorbent material helps contain fluid, and the polypropylene backing helps keep underpad in place.**Usage:** Absorbent material for fluid containment.**Precautions:** Perineal skin is particularly susceptible to damage commonly associated with incontinence.**Size:** 58x91cm**Standard Box Quantity:** 1

Day Liner 7 x 24 inch**Manufacturer:** Covidien **Manufacturer Code:** LSF5141BR-05**Features:** Durable, discreet non-rustle polyethylene backing. Have a special double pad acquisition layer-maximising absorbency where it is most needed in the central area.**Usage:** Have a time/wetness indicator strip which changes colour when the pad is wet instead of "saturated".**Precautions:** Perineal skin is particularly susceptible to damage commonly associated with incontinence.**Size:** 7x24 inches**Standard Box Quantity:** 1

FITEXTRA-LG

Pads/Briefs

Adult Brief - Large 48 -58 inch**Manufacturer:** Medline **Manufacturer Code:** FITEXTRALG**Features:** FitRight® Extra brief promotes a discreet, comfortable, garment like fit to enhance patient dignity, mobility and independence. 4D-Core® with odor protection helps wick fluid away to promote dryness and help maintain skin integrity. Skin-safe closures provide secure, safe and repeated refastenability. Ultra-soft cloth-like backsheet provides a discreet, garment-like, natural feel. Breathable side panels increase airflow and maximize patient comfort.**Usage:** The full absorbent brief designed for persons restricted by incontinence.**Precautions:** Moisture alert wetness indicator changes color when the garment has been soiled. Latex Free.**Size:** Large 48 -58 inches Blue**Standard Box Quantity:** 1

FITEXTRA-MD

Pads/Briefs

Adult Brief - Medium 32 - 42 inch**Manufacturer:** Medline **Manufacturer Code:** FITEXTRAMD**Features:** FitRight® Extra brief promotes a discreet, comfortable, garment like fit to enhance patient dignity, mobility and independence. 4D-Core® with odor protection helps wick fluid away to promote dryness and help maintain skin integrity. Skin-safe closures provide secure, safe and repeated refastenability. Ultra-soft cloth-like backsheet provides a discreet, garment-like, natural feel. Breathable side panels increase airflow and maximize patient comfort.**Usage:** The full absorbent brief designed for persons restricted by incontinence.**Precautions:** Moisture alert wetness indicator changes color when the garment has been soiled. Latex Free.**Size:** Medium 32-42 inches White**Standard Box Quantity:** 1

FITEXTRA-XLG

Pads/Briefs

Adult Brief - Xlarge 57 - 66 inch

Manufacturer: Medline **Manufacturer Code:** FITEXTRAXLG

Features: FitRight® Extra brief promotes a discreet, comfortable, garment like fit to enhance patient dignity, mobility and independence. 4D-Core® with odor protection helps wick fluid away to promote dryness and help maintain skin integrity. Skin-safe closures provide secure, safe and repeated refastenability. Ultra-soft cloth-like backsheet provides a discreet, garment-like, natural feel. Breathable side panels increase airflow and maximize patient comfort.

Usage: The full absorbent brief designed for persons restricted by incontinence.

Precautions: Moisture alert wetness indicator changes color when the garment has been soiled. Latex Free.

Size: Xlarge 57-66 inches Beige

Standard Box Quantity: 1



ML281224

Pads/Briefs

Incontinence Pad 44 x 61cm

Manufacturer: Medline **Manufacturer Code:** MSC281224

Features: Protection Plus Disposable Underpads with non-woven facing are soft against skin. Absorbent material helps contain fluid, and the polypropylene backing helps keep underpad in place.

Usage: Absorbent material for fluid containment.

Precautions: Perineal skin is particularly susceptible to damage commonly associated with incontinence.

Size: 44x61cm

Standard Box Quantity: 1



ML320201

Pads/Briefs

Mesh Panty Medium

Manufacturer: Medline **Manufacturer Code:** MBP3201

Features: Protection Plus® Short Term Use Mesh Incontinence Pants are disposable. Medium fits waist size 28"-40"

Usage: Ideal for maternity use.

Precautions: Latex-free.

Size: Medium

Standard Box Quantity: 1



SAN64211

Pads/Briefs

Mesh Pants Small/Medium

Manufacturer: SCA **Manufacturer Code:** 64211

Features: These washable reusable pants ensure that pads stay securely in place. The soft, elastic edges and latex-free material guarantee comfort and can be washed several times without losing their shape. Large fits waists 20" - 37"

Usage: Ideal for maternity use. Suitable for both men and women.

Precautions: Latex-free.

Size: Small/Medium 20-37in

Standard Box Quantity: Each



SAN64222

Pads/Briefs

Mesh Panty Large/Xlarge

Manufacturer: Molnlyke **Manufacturer Code:** 64222

Features: These washable reusable pants ensure that pads stay securely in place. The soft, elastic edges and latex-free material guarantee comfort and can be washed several times without losing their shape. Suitable for both men and women. Large fits waists 28" - 52"

Usage: Ideal for maternity use.

Precautions: Latex-free.

Size: Large

Standard Box Quantity: 1



COL306

Pediatric Catheters

6FR 10 inch Pediatric Latex Free Catheter

Manufacturer: Coloplast **Manufacturer Code:** 504410

Features: Straight tip intermittent catheter. Straight tip fire polished eyelets. Siliconized surface for smooth insertion.

Usage: For pediatric use.

Precautions: Latex free.

Size: 6FR 10 inches

Standard Box Quantity: 1



8FR 3cc Silicone Pediatric Catheter**Manufacturer:** Coloplast **Manufacturer Code:** AA6108**Features:** 8Fr Foley Catheter with 3cc Balloon, Coude Tip. The Folsil 100% silicone indwelling catheter features a significantly larger lumen than latex versions of the same FR size. Features a radiopaque stripe in a transparent shaft for easy detection on an x-ray.**Usage:** Helps to reduce obstruction and inhibit encrustation. The larger diameter eyelets allow for better drainage while an atraumatic tip adds to ease of insertion.**Precautions:** 100% Silicone**Size:** 8FR 3cc 12 inches**Standard Box Quantity:** 1**10FR 3cc Silicone Pediatric Catheter****Manufacturer:** Covidien **Manufacturer Code:** 8887603101**Features:** 100% silicone preferred for medical devices, Less urethral irritation and stricture versus latex material, Less encrustation compared to latex catheters. Smooth, separately molded, reinforced tip resists buckling.**Usage:** For bladder irrigation and drainage.**Precautions:** Do not over inflate balloon past recommended volume. Use sterile water or normal saline. Latex free.**Size:** 10FR 3cc**Standard Box Quantity:** 1**Catheter Tip Syringe 60 cc****Manufacturer:** BD **Manufacturer Code:** BD309620**Features:** 60 mL BD™ syringe catheter tip 2 oz.**Usage:** General Use Syringes.**Precautions:** Standard hygienic practices are needed to use this product safely. Contaminated syringes must be disposed of according to the nature of the contaminate.**Size:** 60cc**Standard Box Quantity:** 1

BARVUPD1012

Securements

Securement Statlock for Catheter 10-12FR

Manufacturer: Bard **Manufacturer Code:** VUPD1012

Features: Advanced technologies now combine with the StatLock® Stabilization Device design sophistication to deliver skin-safe, suture-free, releasable stabilization of even the most lubricious catheters. The StatLock® Universal Plus Stabilization Device comes in four sizes, to fit virtually all percutaneous drainage catheters.

Usage: For suture-free stabilization of biliary and nephrostomy percutaneous drainage, and peritoneal lavage catheters.

Precautions: More effective than traditional forms of tape and suture stabilization to achieve positive outcomes.

Size: 10-12FR

Standard Box Quantity: 1



BARVUPD1214

Securements

Securement Statlock for Catheter 12-14FR

Manufacturer: Bard **Manufacturer Code:** VUPD1214

Features: Advanced technologies now combine with the StatLock® Stabilization Device design sophistication to deliver skin-safe, suture-free, releasable stabilization of even the most lubricious catheters. The StatLock® Universal Plus Stabilization Device comes in four sizes, to fit virtually all percutaneous drainage catheters.

Usage: For suture-free stabilization of biliary and nephrostomy percutaneous drainage, and peritoneal lavage catheters.

Precautions: More effective than traditional forms of tape and suture stabilization to achieve positive outcomes.

Size: 12-14FR

Standard Box Quantity: 1



BARVUPD1416

Securements

Securement Statlock for Catheter 14-16FR

Manufacturer: Bard **Manufacturer Code:** VUPD1416

Features: Advanced technologies now combine with the StatLock® Stabilization Device design sophistication to deliver skin-safe, suture-free, releasable stabilization of even the most lubricious catheters. The StatLock® Universal Plus Stabilization Device comes in four sizes, to fit virtually all percutaneous drainage catheters.

Usage: For suture-free stabilization of biliary and nephrostomy percutaneous drainage, and peritoneal lavage catheters.

Precautions: More effective than traditional forms of tape and suture stabilization to achieve positive outcomes.

Size: 14-16FR

Standard Box Quantity: 1



BARVUPD68

Securements

Securement Statlock for Catheter 6-8FR

Manufacturer: Bard **Manufacturer Code:** VUPD68

Features: Advanced technologies now combine with the StatLock® Stabilization Device design sophistication to deliver skin-safe, suture-free, releasable stabilization of even the most lubricious catheters. The StatLock® Universal Plus Stabilization Device comes in four sizes, to fit virtually all percutaneous drainage catheters.

Usage: For suture-free stabilization of biliary and nephrostomy percutaneous drainage, and peritoneal lavage catheters.

Precautions: More effective than traditional forms of tape and suture stabilization to achieve positive outcomes.

Size: 6-8FR

Standard Box Quantity: 1



CON37449

Securements

Flexitrack Anchoring Device

Manufacturer: ConvaTec **Manufacturer Code:** 37449

Features: FLEXI-TRAK® Anchoring Device is a skin-friendly, modern acrylic adhesive that securely holds tubes, drains and catheters in place on the skin.

Usage: ConvaTec FLEXI TRAK Large Anchoring Device For anchoring tubes and catheters to skin, such as urinary catheters, gastrostomy tubes, chest tubes. Adhesive.

Precautions: This product does not contain latex.

Size: 38mm X 100mm

Standard Box Quantity: 1



BD306575

Solutions

Normal Saline 10cc Prefilled Syringe

Manufacturer: BD **Manufacturer Code:** BD306575

Features: 10 mL BD PosiFlush™ S.P. saline syringe

Usage: Pre-filled flush syringes specifically designed to reduce blood reflux that may be introduced by the plunger rod of standard syringes.

Precautions: Preservative free

Size: 10ml

Standard Box Quantity: 1



Sterile Water Bag For Irrigation (1000ml)**Manufacturer:** Baxter **Manufacturer Code:** JB7114**Features:** Baxter 1000ml sterile water for irrigation.**Usage:** For urethral irrigation given via an indwelling catheter.**Precautions:** Sterile. Not for injection.**Size:** 1000mL**Standard Box Quantity:** 1**Normal Saline, 3L bag****Manufacturer:** Baxter **Manufacturer Code:** JB7127**Features:** Normal Saline 0.9% NaCl.**Usage:** For Continuous Bladder Irrigation 3000 ml.**Precautions:****Size:** 3000ml**Standard Box Quantity:** 1**Swabs - Chlorhexidine (Packet)****Manufacturer:** Solu Med **Manufacturer Code:** 101.06**Features:** Individually sealed wipes for cleaning wounds when water is not readily available.**Usage:** For wound cleansing and pre-dressing preparation. Ensures greater adhesion for dressings.**Precautions:** Should redness or other signs of irritation appear, discontinue use.**Size:** Small 0.65mL**Standard Box Quantity:** Each

3M10207

Swabs/Wipes/Creams/Moisturizers

Chlorhexidine Swab Stick, Non-Alcohol**Manufacturer:** Solu Med **Manufacturer Code:** 102.07**Features:** Individually sealed swabsticks for cleaning wounds when water is not readily available.**Usage:** For wound cleansing and pre-dressing preparation. Ensures greater adhesion for dressings.**Precautions:** Should redness or other signs of irritation appear, discontinue use.**Size:** 1.6mL 4.75in**Standard Box Quantity:** Each

3M3344

Swabs/Wipes/Creams/Moisturizers

Cavilon No Sting Barrier Wipes**Manufacturer:** 3M **Manufacturer Code:** 3344E**Features:** A fast-drying, non-sticky, alcohol-free liquid barrier film that forms a breathable, transparent coating on the skin. Protects damaged or intact skin from body fluids, adhesive trauma, friction, and incontinence. The film is hypoallergenic, non-cytotoxic and will not sting even when applied to damaged or denuded skin.**Usage:** A hypoallergenic, non-cytotoxic liquid film which will not sting even when applied to damaged or denuded skin.**Precautions:** If Cavilon No Sting Barrier Film is applied to an area with skin folds or other skin-to-skin contact, ensure that these areas of skin are separated to allow coating to dry thoroughly before allowing skin to return to normal position.**Size:** 1 ml**Standard Box Quantity:** 1

3M3346

Swabs/Wipes/Creams/Moisturizers

Cavilon No Sting Barrier Spray 28ml btl**Manufacturer:** 3M **Manufacturer Code:** 3346E**Features:** A fast-drying, non-sticky, alcohol-free liquid barrier film that forms a breathable, transparent coating on the skin. Protects damaged or intact skin from body fluids, adhesive trauma, friction, and incontinence. The film is hypoallergenic, non-cytotoxic and will not sting even when applied to damaged or denuded skin.**Usage:** A hypoallergenic, non-cytotoxic liquid film which will not sting even when applied to damaged or denuded skin.**Precautions:** If Cavilon No Sting Barrier Film is applied to an area with skin folds or other skin-to-skin contact, ensure that these areas of skin are separated to allow coating to dry thoroughly before allowing skin to return to normal position.**Size:** 28ml**Standard Box Quantity:** 1


AMG118305

Swabs/Wipes/Creams/Moisturizers

Lubricant-Water Soluble 140 g Tube**Manufacturer:** AMG **Manufacturer Code:** 118-305**Features:** Sterile. Water Soluble. Non Greasy. Unscented.**Usage:** Ideal for both medical procedures and personal use.**Precautions:** Product is sterile until seal is opened, so as to ensure product safety for medical use.**Size:** 140 g**Standard Box Quantity:** 1

CHCS42050

Swabs/Wipes/Creams/Moisturizers

Betadine Swabsticks**Manufacturer:** PDI **Manufacturer Code:** S42050**Features:** Povidone iodine swabstick**Usage:****Precautions:****Size:** 4 inch**Standard Box Quantity:** Each

COL2041

Swabs/Wipes/Creams/Moisturizers

Coloplast Prep Wipe**Manufacturer:** Coloplast **Manufacturer Code:** 2041**Features:** Clear film forms a breathable barrier, protecting against skin irritants such as adhesives, urine and faeces.**Usage:** Used for cleansing and protecting the skin. Single application wipes.**Precautions:** Not to be used with Extra Extended Wear Skin Barriers.**Size:****Standard Box Quantity:** 54

Sween Cream 190ml**Manufacturer:** Coloplast **Manufacturer Code:** 7068

Features: The non-occlusive moisturizers and natural vitamins A&D, immediately help restore dry, fragile skin that's at risk for cracking, especially around contact points like heels and elbows. Non-occlusive formula absorbs quickly into the skin and will not interfere with the adhesion of tape or ostomy appliances.

Usage: Use Sween Cream® for: Cracked and flaky dry skin; Red skin due to minor burns; Radiation/chemotherapy skin reactions; Peri-wound skin; Perineal Erythema; Surgical incisions - reduces itching and pulling discomfort; Relieves itching under tape; Callouses; Any reddened skin.

Precautions: For external use only. If irritation occurs, consult a physician.

Size: 190ml**Standard Box Quantity:** 1**Critic Aid Cream Clear 71gm****Manufacturer:** Coloplast **Manufacturer Code:** 7573

Features: A clear petrolatum based moisture barrier ointment containing dimethicone. Absorbs skin moisture allowing it to adhere to denuded skin.

Usage: For light to heavy skin irritation due to incontinence.

Precautions: Do not use this medication if you are allergic to zinc, dimethicone, lanolin, cod liver oil, petroleum jelly, parabens, mineral oil, or wax. Zinc oxide topical will not treat a bacterial or fungal infection. Avoid getting this medication in your mouth or eyes.

Size: 71g**Standard Box Quantity:** 1**Convatec AllKare Adhesive Remover Wipe****Manufacturer:** ConvaTec **Manufacturer Code:** 37443

Features: Wipes contain an emollient-rich formula of Iso-alcohol-free solution that is non-irritating and non-stinging.

Usage: Wipes remove the film left from tape, wafer and skin adhesive that can reduce wearing time.

Precautions: AllKare wipes do not adversely affect Stomahesive or Durahesive skin barrier adhesion or wear time.

Size:**Standard Box Quantity:** 100

Convatec AllKare Protective Wipe**Manufacturer:** ConvaTec **Manufacturer Code:** 37444**Features:** Wipe is a non-water soluble protective barrier which resists breakdown and "breathes" for greater comfort.**Usage:** Wipes leave a thin, comfortable film, which protects against skin stripping, tension blisters, folliculitis, maceration and adhesive residue.**Precautions:** AllKare wipes do not adversely affect Stomahesive or Durahesive skin barrier adhesion or wear time.**Size:****Standard Box Quantity:** 100**Alcohol Swabs****Manufacturer:** Covidien **Manufacturer Code:** 63001**Features:** Used to wipe needles, thermometers, syringes and smaller body parts like finger tips, nails etc.**Usage:** Used to wipe needles, thermometers, syringes and smaller body parts like finger tips, nails etc.**Precautions:** Should redness or other signs of irritation appear, discontinue use.**Size:** Medium**Standard Box Quantity:** 1**Skin Prep Wipe Packet Pkg****Manufacturer:** Covidien **Manufacturer Code:** 6560**Features:** Forms protective barrier between skin and adhesive products. Helps reduce the chance of excoriation. Suitable for use in both hospital and long-term care. One wipe per packet helps reduce waste by using only what is needed.**Usage:** Helps protect skin from feces, urine, chemical irritants and gastrointestinal fluids.**Precautions:** Non-sterile.**Size:** Large**Standard Box Quantity:** Each

No-Sting Remover Wipes 50/Box**Manufacturer:** Hollister **Manufacturer Code:** 7760**Features:** Gentle solvent removes adhesives and cements; makes appliance removal fast and easy. Effective for removal of most adhesives and safe to use on adult skin. Easy to use in a convenient wipe package.**Usage:** Product has features that would be most appropriate for post-operative care.**Precautions:** Recommended for Home Use.**Size:****Standard Box Quantity:** 50**Skin Prep Wipe****Manufacturer:** Medline **Manufacturer Code:** MSC1500**Features:** Sureprep is a fast drying skin protectant that forms a protective film to help reduce friction during removal of tapes and films. Vapor permeable and delivers protection from friction and incontinence.**Usage:** The transparent barrier may be used on periwound, peristomal or areas that come in contact with bodily fluids. Frequency of application should be at each dressing change or 24-72 hours depending on frequency of incontinence.**Precautions:** Contains alcohol. Do not apply directly to open wounds or infected skin. Should redness or other signs of irritation appear, discontinue use.**Size:****Standard Box Quantity:** 50**Povidone-iodone Swab Sticks****Manufacturer:** Dynarex **Manufacturer Code:** S42325**Features:** Aqueous solution of 10% povidone-iodine. It is a fast-acting, broad-spectrum antiseptic that kills gram-positive and gram-negative bacteria (including antibiotic resistant organisms), as well as most fungi/yeasts, viruses and protozoa. Use swab to coat area with iodine solution.**Usage:** It is indicated for degerming skin, wounds and mucous membranes, and for preoperative skin preparation of patients. Helps reduce bacteria that potentially can cause skin infection.**Precautions:** May cause burning, redness, or a rash. Notify your doctor if you develop any of these symptoms. If you notice other effects not listed above, contact your doctor or pharmacist.**Size:****Standard Box Quantity:** PKG/3

SM1322N

Swabs/Wipes/Creams/Moisturizers

Individual Lubricants for Catheters

Manufacturer: Muko **Manufacturer Code:** SM1322N

Features: Bacteriostatic, non-conductive water soluble. Eases insertion of urological catheters, rectal thermometers, enemas, tampons. Lubricates condoms and provides personal lubrication.

Usage: Eases insertion of urological catheters, rectal thermometers, enemas, tampons. Lubricates condoms and provides personal lubrication.

Precautions: Sterile.

Size: 3.5g

Standard Box Quantity: 1



SN59431279

Swabs/Wipes/Creams/Moisturizers

Secura Cream (Zinc) 78gm(tube)

Manufacturer: Smith & Nephew **Manufacturer Code:** 59431279

Features: Restores moisture to very dry skin, Contains petrolatum and Vitamin E, Conditions and soothes dry skin.

Usage: Moisturizes and protects dry skin.

Precautions: For external use only. If irritation occurs, consult a physician.

Size: 78g

Standard Box Quantity: 1



3M1527/1

Tapes

Transpore Non-Allergic 2.5cm roll

Manufacturer: 3M **Manufacturer Code:** 1527-1

Features: A latex-free, hypoallergenic, transparent and perforated plastic tape that offers strong adhesion. Transpore tape is easy to tear into very thin strips, and is easy to use even with gloves.

Usage: Securing tubing and devices (e.g., catheters, I.V. tubing). Anchoring bulky dressings.

Precautions: Apply tape without tension. Assure that skin is clean and dry before applying tape. Protect skin with 3M™ Cavilon™ No Sting Barrier Film allowing it to dry before applying tape. If possible, avoid using tackifiers.

Size: 2.5cm

Standard Box Quantity: 1



3M15272

Tapes

Transpore Non-Allergic 5cm roll

Manufacturer: 3M **Manufacturer Code:** 1527-2

Features: A latex-free, hypoallergenic, transparent and perforated plastic tape that offers strong adhesion. Transpore tape is easy to tear into very thin strips, and is easy to use even with gloves.

Usage: Securing tubing and devices (e.g., catheters, I.V. tubing). Anchoring bulky dressings.

Precautions: Apply tape without tension. Assure that skin is clean and dry before applying tape. Protect skin with 3M™ Cavilon™ No Sting Barrier Film allowing it to dry before applying tape. If possible, avoid using tackifiers.

Size: 5cm

Standard Box Quantity: 1



3M1530/1

Tapes

Micropore (Paper) 2.5cm roll

Manufacturer: 3M **Manufacturer Code:** 1530-1

Features: A latex-free, hypoallergenic paper tape that is gentle to the skin yet adheres well and leaves minimal adhesive residue upon removal. General purpose, breathable surgical tape.

Usage: Securing small to medium dressings. Securing lightweight tubing. Securing ostomy pouches. Taping fragile, at-risk skin when repeated taping is needed.

Precautions: Apply tape without tension to clean, dry skin. Protect at-risk skin with 3M™ Cavilon™ No Sting Barrier Film allowing it to dry before applying tape. In general, do not encircle a limb completely with tape. Tape should extend at least one-half inch beyond the dressing. Tape should not be pulled or stretched when applied.

Size: 2.5cm

Standard Box Quantity: 1



3M15302

Tapes

Micropore (Paper) 5cm roll

Manufacturer: 3M **Manufacturer Code:** 1530-2

Features: A latex-free, hypoallergenic paper tape that is gentle to the skin yet adheres well and leaves minimal adhesive residue upon removal. General purpose, breathable surgical tape.

Usage: Securing small to medium dressings. Securing lightweight tubing. Securing ostomy pouches. Taping fragile, at-risk skin when repeated taping is needed.

Precautions: Apply tape without tension to clean, dry skin. Protect at-risk skin with 3M™ Cavilon™ No Sting Barrier Film allowing it to dry before applying tape. In general, do not encircle a limb completely with tape. Tape should extend at least one-half inch beyond the dressing. Tape should not be pulled or stretched when applied.

Size: 5cm

Standard Box Quantity: 1



Leukoplast Sleek LF 5cm x 3M waterproof pink tape, latex free**Manufacturer:** Smith & Nephew **Manufacturer Code:** 7235907**Features:** Plastic Waterproof Adhesive Tape. Often referred to as sleek or pink tape. Leukoplast® Sleek is easy to unroll, easy to apply and can stay on for hours and days. This tape resists sweat, grease, oil and water.**Usage:** Indicated for any situation where a reliable and conformable waterproof tape is required (for holding tubes or tracheotomy).**Precautions:** In general, do not encircle a limb completely with tape. Tape should extend at least one-half inch beyond the dressing. Tape should not be pulled or stretched when applied.**Size:** 5cm x 3m**Standard Box Quantity:** Each**Leukoplast Sleek 5cm x 3M waterproof pink tape****Manufacturer:** BSN **Manufacturer Code:** BSN7237906**Features:** Plastic Waterproof Adhesive Tape. Often referred to as sleek or pink tape. Leukoplast® Sleek is easy to unroll, easy to apply and can stay on for hours and days. This tape resists sweat, grease, oil and water.**Usage:** Indicated for any situation where a reliable and conformable waterproof tape is required (for holding tubes or tracheotomy).**Precautions:** In general, do not encircle a limb completely with tape. Tape should extend at least one-half inch beyond the dressing. Tape should not be pulled or stretched when applied.**Size:** 5cm x 3m**Standard Box Quantity:** 1**Pontoon Bedpan****Manufacturer:** Medline **Manufacturer Code:** DYND88205**Features:** Designed to be autoclavable at 250 degrees F, 15 lbs. pressure for 15-20 minutes. In autoclavable blue for easy identification.**Usage:** Bed pan.**Precautions:****Size:** 11 x 14 x 4 inches**Standard Box Quantity:** 1

ML80234

Urinals/Bedpans

Urinal 1L capacity

Manufacturer: Medline **Manufacturer Code:** DYND80234

Features: Contoured handle attaches easily to health care beds. Attached lid prevents spilling and reduces odors. Graduated to 32 oz/1000 mL

Usage: Urinal for non-ambulatory patients.

Precautions:

Size: 1L

Standard Box Quantity: 1



ML8521

Urinals/Bedpans

Slipper/Fracture Bedpan

Manufacturer: Medline **Manufacturer Code:** DYNC8521

Features: This smaller bedpan features a tapered front so that it slides easily under the body. Also used as a female urinal. Loop handle makes maneuvering easy

Usage: Non-stick surface.

Precautions: Not autoclavable.

Size: 9 x 12 x 3.25 inches

Standard Box Quantity: 1

